

A close-up photograph of a white, modern patient support rail. The rail is made of smooth, curved plastic or metal. On the left side, there is a prominent circular joint mechanism with a light blue or grey center. The rail extends horizontally across the frame. The background is a soft, out-of-focus light grey. A solid cyan-colored triangular shape is in the bottom right corner.

Ritter®

MEDICAL

## Patient Support Rails



# Better support throughout the entire exam

01



Visual Indicator Marks



02



03



## 01 ROTATION/LOCKING POINTS

Patient Support Rails can be locked into three positions—front horizontal or two elevated positions—or moved to a fully back horizontal position. Visual indicator marks help guide the operator to the different locking positions.

## 02 ON-CHAIR POSITIONING

When rotated to the back horizontal position, the support rails can be used as safety rails, designed to improve patient security and provide a gripping point during movement into supine, prone, left or right lateral positions.


## 03 CONTINUOUS GRIPPING SURFACE

Patient Support Rails provide patients with a 1 ¼ inch diameter continuous gripping surface for entering, exiting or repositioning on the exam chair. They can be used as an affixed gripping point for patients to independently raise themselves off the chair.


# SPECIFICATIONS

### PATIENT SUPPORT RAILS MODEL

9A600001 – Patient Support Rails (field installed)



FRONT VIEW



SIDE VIEW

# Better arm support when you need it most



01



02

03



Articulating Arm Pad

## 01 BETTER BP

For a more accurate blood pressure measurement, the patient's arm should be supported at heart level. If the upper arm is below heart level or unsupported, the readings will be too high. If the upper arm is above heart level, the readings will be too low.<sup>1</sup>

## 02 BETTER POSITIONING

The Patient Support Rails+ accessory functions the same as the Patient Support Rails with the added benefit of helping providers achieve proper positioning of the patient's arm for blood pressure capture and blood draw applications.


## 03 BETTER ARM SUPPORT

Easily support the patient's arm at heart level with the Patient Support Rails+ accessory. The unique articulating arm pad design allows for more precise positioning for blood pressure capture or other similar applications.


## SPECIFICATIONS

### PATIENT SUPPORT RAILS+ MODEL

9A600002 – Patient Support Rails+ (field installed)



FRONT VIEW



SIDE VIEW

<sup>1</sup> <https://www.aafp.org/afp/2005/1001/p1391.html>

Midmark is an ISO 13485 and ISO 9001 Certified Company.\*  
For more information, contact your Midmark dealer  
or call: 1.800.MIDMARK  
Outside the USA call: 1.937.526.3662 or visit our  
website: [midmark.com](http://midmark.com)

© 2019 Midmark Corporation, Miamisburg, Ohio USA  
Products subject to improvement changes without notice.  
Litho in USA.

007-10146-00 Rev. A1 (6/19)

\*Certain products not included.  
See the complete list at: [midmark.com/ISO](http://midmark.com/ISO)

# Ritter®