

CO₂ Sensor Best Practices


ALWAYS

"Zero" the CO₂ sensor after a new adapter or filter line is installed

NEVER

"Zero" a sensor without an adapter or filter line installed
"Zero" a sensor that is in use or attached to a patient

DISPOSE FREQUENTLY



Choose the Right Adapter or Filter Line

PEDIATRIC

≤4.0
ET Tubes

ADULT

≥4.5
ET Tubes

"Zero" the CO₂ Sensor*

- 1 Plug the CO₂ sensor into the Midmark Multiparameter Monitor.
- 2 Install the airway adapter or filter line.
- 3 Select the CO₂ waveform or CO₂ data in the parameter box to enter the CO₂ Setup Menu.
- 4 Once the module is warmed, press "CO₂ Zeroing" or "Start Zero Calibration".
- 5 Upon successful completion, the status bar will display the message "CO₂ Zero Successful".

*These steps are specific to zeroing Respirationics CO₂ sensors.



Still Have Questions?

Call 1.844.856.1232, press 2 to speak with a clinical support specialist who can guide you through specific questions.