Product Bulletin IQvitals® Zone™ Firmware Update Procedure

Applies to Models:
1-200-0310 (Kit 4-000-0620)
1-200-0320 (Kit 4-000-0622)
1-200-0330 (Kit 4-000-0624)
1-200-0340 (Kit 4-000-0720)
1-200-0350 (Kit 4-000-0722)
1-200-0360 (Kit 4-000-0724)

Items needed:
2-200-0235 Firmware Update Kit

Important: If you are using a PC running Midmark IQmanager® software or EMR system with the IQvitals Zone, make sure the PC is updated to IQvitals plug-in version 11.0.4 or higher before updating the firmware.

- Closely follow the update instruction below.
- Do not interrupt the process until the update is complete.

Language of origin:
English

Purpose: This bulletin provides the instruction to update the firmware of IQvitals Zone, which adds the following features and enhancements:

- *SPRINT/Custom blood pressure protocols*
- *Improves battery management*
- *Remote firmware update capability*
- *Event Log recording*
- *One-click to Step Deflation prompt*
- *Averaging, Interval, Continuous, SPRINT and Custom protocol NIBP modes available in Step-Deflation method only*


Step 1: Locate Firmware USB port

- Ensure device is connected to wall power.
- Ensure the SpO₂ sensor is connected to the device.
- Locate firmware USB port on back of device.

Step 2: Update with USB 1


- Ensure device is turned off.
- Insert USB 1 into firmware USB port.
- Turn device on.
- Firmware will update automatically.
- When "Program...DONE" message is displayed for DISPLAY PROCESSOR APPLICATION, USB 1 update is complete.
- Turn off device and remove USB 1.

USB 1 Firmware update complete

Step 3: Update with USB 2

- Ensure device is turned off.
- Insert USB 2 into firmware USB port.
- Turn device on.
- Firmware will update automatically.
- When complete, **Do Not Remove USB 2.**
- Turn off device with USB 2 still plugged in.


USB 2 Onscreen instructions

Product Bulletin

IQvitals® Zone™ Firmware Update Procedure

Step 4: Keep USB 2 in device

- Turn device on.
- Firmware update continues automatically.
- When “Program...DONE” messages are displayed for MAIN PROCESSOR APPLICATION, DISPLAY PROCESSOR APPLICATION, DISPLAY SPI FLASH GRAPHICS and NIBP MAIN PROCESSOR APPLICATION, the firmware update is complete.
- Turn device off.
- Remove USB 2.


USB 2 Firmware update complete

Step 5: Verify Firmware Version

- Turn device on.
- Select Settings.
- Select “More” tab.
- Enter password: 1590.
- Select FACTORY DEFAULT SETTINGS.
- Select Yes.
- Verify.

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