



Exam Assistant® Drawer System
To use the drawer front as a shelf...
 Pull drawer front down to release magnetic latch.

If drawer front folds past horizontal...
 Pull drawer front upward until it snaps back into the horizontal position.

Equipment Alert
 The maximum weight the shelf will support is 10 lbs (4.5 kg).

Optional Features:

The table may be equipped with a drawer heater and / or upholstery heater



Optional Features
To turn heater(s) ON / OFF...
 Press the desired heater switch. (Switch will illuminate when heater is ON)

Note: The upholstery heater will...

- reach operating temp. in 15 - 20 minutes.
- automatically shut off after 8 hours.

Drawer Heater Switch

Upholstery Heater Switch

MA8112p

Maintenance

Periodically inspect the following:

- Power cord must be free of visible damage.
- All fasteners must be present and secure.
- All mechanical functions should operate

Cleaning

Upholstery

Wash your upholstery weekly with a mild liquid soap and water mixture, rinse with clear water and dry completely to remove disinfection cleaner build-up.

Disinfect your upholstery with a solution of standard bleach and water mixed 1 in 10 (10%) or chlorine based cleaners (see current CDC Guideline for Disinfection & Sterilization in Healthcare Facilities). Follow this with a clear water rinse and thorough drying of material.

To minimize damage caused by disinfectant cleaner residue build-up, do not allow disinfectants to pool on the upholstery surface. Once the approved contact time has been obtained, remove and dry excess liquid remaining on the surface.

Detailed care and maintenance instructions are included with your product. This information is also available on www.midmark.com in the Technical library under the User Information tab for your product.

Painted Metal / Plastic Surfaces

Clean the painted metal and plastic surfaces weekly using a clean soft cloth, and mild cleaner. Periodic applications of common furniture wax will ease cleaning, and maintain the finished luster of the surfaces. Use of chlorine based disinfectant can cause damage to the hand control.

- Français : Référez-vous à www.midmark.com
- Deutsch: Siehe www.midmark.com
- Español: Refiera a www.midmark.com



Midmark Corporation
 60 Vista Drive
 Versailles, OH 45380-0286
 Phone: 937-526-3662
 Fax: 937-526-5542
www.midmark.com



Quick Reference Guide
625 (-001, -003, -005, -006)
Barrier-Free® Exam Table

Model Number:

Serial Number:

Dealer Information:

Purchase Date:


Upholstery Color:

003-2021-00 Rev D 7/10/15)

This guide outlines basic operation of the table. The enclosed User Guide provides additional instructions on operation and care, as well as regulatory compliance information. If service is required, contact your authorized Midmark dealer.



MA8113p





To operate the Base / Back functions...
 Press & hold the Up or Down arrow button for the desired function.

*Note: The table automatically stops at its maximum and minimum height.
 In the event of a malfunction, press the Stop button to terminate movement.*

To activate the Home -or- Quick Exam® function...
 Press & release the Home -or- Quick Exam® button.

Note: In the event of a malfunction, press the Stop button to terminate movement.

Equipment Alert
 If the table's weight limit is exceeded, you will hear a "beep" and the function will be disabled. If the problem occurs with a normal load, contact your service technician.



Quick Exam® Height
 Programmable
 (set by the user)

'Home' Height
 $20 \pm .5 \text{ in. (} 50.8 \pm 1.3 \text{ cm)}$



To program the Quick Exam® function...

- Press the Base Up button to move the table to the desired position.
- Press & hold the Quick Exam® & Base Up buttons simultaneously for three seconds.

Note: You will hear a brief chirp when the new Quick Exam® position is saved.