



CareExchange Essentials

Coiled Power Cord Plug Holder

Coiled power cord plug holder is conveniently placed within reach.



Locking Casters (2 front casters)

(2) Two locking 100mm casters. To lock casters, press lock down. To unlock, press lock up.



Codelock (if applicable)

To access medical drawer initially, enter access code #2244 on keypad and flip handle up or down.



CareExchange Essentials

Height Adjustment Hand Lever

To adjust the height, simply press paddle and adjust to preferred height.



Sit-to-Stand Height Adjustment

Flat Screen Monitor Arm (if applicable)

Flat screen monitor arm can be adjusted up or down with the easy-to-use cam lock.



Monitor Height Adjustment



LED Rotation



LED Tilt

Swivel-out Mouse

Swivel-out mousing gives the user ample area (8"W x 9"D)






CareExchange Essentials

Charge Before Deployment

LiFe (up to 2 hrs.) / SLA (up to 5 hrs.)

Once your Power System is fully set up with all electronic devices, it should be plugged into a medical grade outlet. **The Power System should fully charge before deploying.** (Note: All Powered Systems are shipped with a 30% charged LiFePO4 Battery and a 100% charged SLA. If storing in warehouse before deploying, charge fully and keep plugged in).

External User Interface



- Power Button (last button on right) - To turn power on, press and hold for 3 seconds. Indicator will beep and power on. To turn off power, press and hold for 2 second. Status Indicator will click and power off.
Note: Monitors, Laptops and other devices will need to be manually turned on.
- On/Off Indicator (first of five LED lights) - When lit indicates power system is on and technology component is receiving power from Power System.
- Battery Charge Level Indicator (Second through fifth LED Lights): Battery charge level indicator shows the battery state-of-charge and has a two-stage warning system that will alarm when battery charge is low. 1st stage level indicator (second LED light will flash Yellow) - 10% charge remaining. 2nd stage level indicator (second LED light will flash Red) - 5% charge remaining. This applies to both Lithium Iron Phosphate (LiFePO4) and Sealed Lead Acid (SLA) batteries.
- Mute Button (second to last button on right) - If the battery level falls below set alarm level, the charge an alarm will sound. This button will mute the sound.

CareExchange Essentials

External User Interface - LED Charge Level Overview

LED Position 1 (Power On/Off)	
Solid Green - System is on and all components will receive power.	
No Light - System is turned off. If cart is plugged in, unit can still be charged.	
LED Position 2-5 (Charge Level)	
Flashing Red - Less than 5% battery life remaining. Plug in unit immediately.	
Flashing Yellow - Less than 10% battery life remaining. Plug in unit as soon as possible.	
Solid Green - Battery is at least 90-100% charged. Component should still be plugged in when possible to maximize battery life.	
Flashing Green - Each LED flashing green reflects the level of charge achieved during recharging.	
No Light - Power System is turned off or battery is not properly connected to the power system.	
LED Position 5 (Charging Status when plugged into wall outlet and turned on)	
Solid Green - Battery is full	