

HEINE® Li-ion L rechargeable battery


 Please read and follow these instructions for use and keep them for future reference.


Intended Use

The HEINE® rechargeable batteries are intended solely for HEINE® rechargeable handles.


 **Federal law restricts this device to sale by or on the order of a Physician or Practitioner!**

Warnings and Safety Information

 **CAUTION!** Indicates potential hazardous situations. Ignoring the corresponding instructions may lead to dangerous situations of mild to moderate extent. (Background color, yellow; Foreground color, black for color print, or black in light-gray background.)

 **NOTE!** Indicates valuable advice in terms of set up, operation, maintenance or repair, as applicable. Notes are important, but not related to hazardous situations.

General Warnings


 Do not put into fire, do not short-circuit, do not modify, do not sterilize, do not drop, do not expose to severe shocks or to direct sunlight.

Dispose of batteries in the event of visible or presumed damage (e.g. leaking liquid, mechanical damage).

The batteries can be damaged by liquid. They must not be rinsed in running water or submerged under water.

After hygienic reprocessing of the base caps, allow 15 min for the cleaning agent to evaporate before using the handle.

General Notes


 The warranty for the entire product is invalidated if non-genuine HEINE products or non-original parts are used and if repairs or modifications are made to the device by persons not authorized by HEINE. For more information, please visit www.midmark.com.

Please observe the correct polarity or the instruction/arrow on the rechargeable battery.

Check the rechargeable batteries regularly. In case of any abnormalities in the operating or recharging times, replace the rechargeable batteries.

If the rechargeable battery has not been recharged for a long period it may discharge automatically. If this occurs a longer charging time will be required when recharging.

Disposal

 Dispose at your local collection point.

Technical specification














Environmental conditions for operation	+10 °C to +35 °C 30 % to 75 % rel. humidity 700 hPa to 1060 hPa
Environmental conditions for storage	+5 °C to +45 °C 45 % to 80 % rel. humidity 500 hPa to 1060 hPa
Environmental conditions for transport	-20 °C to +50 °C 45 % to 80 % rel. humidity 500 hPa to 1060 hPa
Recommended charging temperature of the rechargeable batteries	+10 °C to +35 °C
Cell type	Lithium-ion technology
Rated voltage [V]	3.7
Capacity [mAh]	2300
Charging time	approx. 2 h

CE




Explanation of utilized symbols

The following symbols are used on the device or on the packaging:

	The CE mark indicates that the product complies with the European medical device directive 93/42/EEC.
	Catalogue- or order number
	Manufacturer
	Date of manufacture
	Dispose at your local collection point
	Temperature limits in °C for storage and transport
	Temperature limits in °F for storage and transport
	Humidity limitation for storage and transport
	Pressure limitation for storage and transport
	Fragile, handle with care!
	Keep dry!
	Follow instructions for use! (Background color, blue; Foreground color, white for color print, or white on black background.)
	Instructions for use

Manufacturer

 **HEINE Optotechnik GmbH & Co. KG**
Kientalstr. 7 · 82211 Herrsching · Germany
www.heine.com

Distributed by

Midmark Corporation
60 Vista Drive
Versailles, Ohio 45380
www.midmark.com