

641 Oral Procedure Chair



Note: Chair shown with optional arm rests

PERFORM PROCEDURES MORE COMFORTABLY AND EFFICIENTLY WITH VERSATILE POSITIONING AND INCREASED ACCESSIBILITY.

The 641 Oral Procedure Chair is designed to help you better position patients while maintaining an ergonomically correct posture throughout procedures. Clear chair access and a large range of motion facilitate optimum positioning that can help patients remain calm and relaxed. Hand and foot controls with the auto-return option allow simple, efficient patient entry and exit. The adjustable head rest and variety of patient support accessories not only increase comfort for your patients, but also assist in ideal head and neck positioning.



01



02




03

- 01** The Elevated Flat position is designed to improve access and visibility to the patient's head, neck and torso.
- 02** The Semi-Fowler position is designed to improve patient comfort and provide better ergonomics for the dental clinician.
- 03** Choose from multiple headrest styles. An adjustable headrest can accommodate taller patients with an expanded range of patient positioning.

SPECIFICATIONS AND ACCESSORIES

Patient Weight Capacity	450 lb (204 kg)
Minimum Height	19" (48.3 cm)
Maximum Height	40" (101.6 cm)
Seat Tilt Angle Range	0° to 30°
Back Support Angle Range	0° to 85°
Foot Support Angle Range	0° to 90°
Optional Accessories	<ul style="list-style-type: none"> Chair Arms Articulating Armboard Caster Base Casters Wireless, Linear-Style Control 3-Foot Long Power Cord Foot Section Cover Double-Arm Instrument Tray Single-Arm Instrument Tray IV Pole Base Rail Security Straps
Headrests	<ul style="list-style-type: none"> Flat Headrest U-Shaped Headrest Oval-Shaped Headrest Rectangular-Shaped Headrest



FEATURED ACCESSORIES



Articulating Armboard



Adjustable Chair Arms



Security Straps



Double Arm Instrument Tray




Casters




IV Pole


HEADREST OPTIONS




Flat Headrest



Oval-Shaped Headrest



U-Shaped Headrest



Rectangular-Shaped Headrest