

VERSION MATRIX

Epic integration with Midmark relies on the correct version of Midmark software being installed on your servers and workstations. Refer to the chart below for software compatibility. Contact us to discuss purchasing software.

EPIC® SOLUTION	MIDMARK SOFTWARE
Hyperspace® (Before the February 2022 release)	Midmark Controls (IQiE) v8.6.1 – v10.1
Hyperspace (February 2022 – November 2022)	Midmark Controls (IQiE) v8.6.1 – v10.1 Midmark Connect for Epic
Hyperdrive	Midmark Connect for Epic

Epic and Hyperspace are registered trademarks of Epic Systems Corporation.

MKT-00512 | © 2022 Midmark Corporation, Miamisburg, Ohio USA

