

Monitoring and Critical Care

Midmark Cardell® Multiparameter Monitor Comparison Chart

Features	9405	9503	9500	MAX-1	MAX-12HD
Display:	5.7" Black & White	7.0" HD Color	7.0" HD Color	10.4" Color	12.1" HD Color
Printer:	Optional (external)	Standard (built-in)	Standard (built-in)	Standard (built-in)	Standard (built-in)
Waveform display:	3-Channel (user-defined)	4-Channel (static)	4-Channel (static)	4-Channel (user-defined)	6 to 8-Channel (user-defined)
ECG Vector display:	2-Channel	Single-Channel	Single-Channel	Single-Channel	4-Channel
ECG recall:	None	None	None	None	12 Minutes
Standard Cardell® BP, Nellcor SpO₂, ECG, Temp, Resp:	Yes	Yes	Yes	Yes	Yes
CO₂:	Yes	No	Optional Upgrade	Yes	Yes
Invasive Blood Pressure:	No	No	No	Yes	Optional
Anesthetic Agent Monitoring:	No	No	No	No	Optional
Graphical Trends:	No	Yes (24 hours)	Yes (24 hours)	Yes (24 hours)	Yes (120 hours)
Tabular Trends:	480	600	600	200	600
Parameter Alarms:	User-defined High/Low (2 settings)	User-defined High/Low	User-defined High/Low	User-defined High/Low	User-defined High/Low with 4 defaults (Cat, Dog, Horse, Other)
Download case record to PC:	Yes (Free Vital Record™ Software)	TBD	No	Yes (Free Vital Record™ Software)	TBD
Interface with Practice Management Software (Avimark, Cornerstone, et al):	Yes (Free Vital Record™ .pdf report)	TBD	No	Yes (Free Vital Record™ .pdf report)	TBD
Video Output:	No	No	No	No	Yes
Weight:	4.4 lbs	8.4 lbs	8.4 lbs	12.8 lbs	16 lbs
Battery Life:	4-6 hours	1.5 hours	1.5 hours	3.0 hours	3.0 hours

For product inquiries, please contact Midmark Veterinary Customer Service at 1-800-MIDMARK, extension 88951.
 For more information or a demonstration, contact Midmark, 60 Vista Drive, P.O. Box 286, Versailles, Ohio 45380-0286.
 Call: 1-800-MIDMARK Fax: 1-877-725-6495 midmark.com

