

---

# *Evaluation of an oscillometric blood pressure monitor for use in anesthetized cats*

---

*Journal of the American Veterinary Medical Association, Vol 221, No. 5,  
September 1, 2002, pp 646-650*

**Conclusions & Clinical Relevance:** “The oscillometric monitor was accurate for DAP and MAP throughout the entire pressure range and met the standards of the Association for the Advancement of Medical Instrumentation...Except for a minor underestimation of SAP during normo- and hypertension, the oscillometric monitor yielded reliable and easily obtainable blood pressure measurements in anesthetized cats.”