



INTEGRATED ACCESSORIES



Sopro Intraoral Cameras

- Sopro Touch allows for easy image capture without blurring
- Lightweight and easy to clean
- Analog and digital (USB 2.0) compatible

LIFE Camera (shown above)

- 3 modes to optimize your dental practice
- Diagnosis aid mode – shows the characterization of lesions at different clinical stages without losing consistency
- Treatment aid mode – images will show a differentiation of pre-operative healthy tissue versus diseased tissue
- “DAYLIGHT” mode – allows the dentist to use SOPRO LIFE as a white light intra-oral camera to produce high-quality images

617 Camera

- Wide 105° angle of view (easy to see distal areas)
- Curved design follows natural shape of the mouth
- High fidelity images
- LEDs utilize condensers to focus light with minimal shadowing
- Fixed-focus and large depth of field (portrait down to a single tooth)



Bien Air Optima MX2 LED

- The #1 manufacturer of electric micromotors in the world
- Brushless, autoclavable motor
- Maintenance-free
- Perfectly silent and very gentle to operate
- Motor speeds of 100 to 40,000 RPM for a broad range of procedures
- Programmable high/low RPM and low RPM torque limits with auto reverse
- Quick-connect coupling rotates over 400°
- Ideal for endodontics
- Accepts E type attachments
- Three year warranty



SATELEC® Mini LED Curing Light

- Second generation LED curing light combining power, efficiency and speed
- Only 6 to 12 seconds to polymerize 3mm of any composite
- Wide light spectrum (420-480 nm)
- Three modes of cure (fast, pulse and ramp)
- Comes standard with light shield and autoclavable 7.5 mm light guide
- Silent operation (no fans)
- Stylish handpiece is compact, lightweight and well balanced
- Shock resistance and easy to clean



SATELEC® SP Newtron® Ultrasonic

- 5th generation Piezo ultrasonic instrument
- Compact, lightweight and autoclavable
- Powerful, yet gentle for applications beyond scaling such as perio, endo and surgery
- Patented “Cruise Control” for automatic regulation of frequency and power
- Color coding system to ensure proper tip and power setting
- Most extensive range of tips in the market, over 80
- 3 tips included (#1, 10Z, H3)
- LED version is also available (shown above)



Dentsply Cavitron® SPS G98A

- The industry's most popular ultrasonic scaler
- Magnetostrictive technology
- SPS technology automatically adjusts system power to maintain tip stroke, compensating for load and conditions
- Water flow adjustment conveniently located on the handpiece
- 30 KHz insert included

Midmark is an ISO 13485 and ISO 9001 Certified Company. For more information or a demonstration, contact your Midmark dealer or call: 1-800-MIDMARK Fax: 1-877-725-6495 Outside the U.S.A. call: 1-937-526-3662 Fax: 1-937-526-8214 or visit our website at midmark.com

© 2009 Midmark Corporation
60 Vista Drive P.O. Box 286 Versailles, Ohio 45380-0286
Products subject to improvement changes without notice
Litho in U.S.A. 007-0804-00 Rev. C1 (2/13)