


SEEING CLEARLY


FOR A COMPLEX DIABETES DIAGNOSIS

Errors Impact Better Outcomes


With an industry-wide push toward value-based care reimbursement models, better outcomes are more important than ever. But, there may be errors and inconsistencies impacting your diagnoses.



of clinicians encounter diagnostic errors at least monthly.¹



of clinicians believe diagnostic errors are preventable.¹



of clinicians feel an inadequate physical examination is a common reason for diagnostic errors.¹

WHAT DOES THAT LOOK LIKE IN PRACTICE?

Think about it through the lens of your patients dealing with complex conditions such as diabetes.



A Better Way

With Midmark HEINE® diagnostic instruments you can see clearly, diagnose with confidence and provide your patients with the care they deserve.

LED HQ
LED NOW IN HEINE QUALITY.

This exclusive technology creates a light that shows a light that shows a light as they truly are.

Dustproof Design:
This design protects the optical components, ensures maintenance-free operation and makes the instrument shock-resistant.

Superior Optics:
Aspherical Optical System (AOS) redirects corneal reflections away from the examiner's path of sight, providing a crisp clear view of the fundus. It eliminates the need for filters and other systems that reduce resolution.

Better view. Better care. Better outcomes.

Sources:

¹ https://www.quantiamd.com/q-qcp/QuantiaMD_PreventingDiagnosticErrors_Whitepaper_1.pdf

² www.cdc.gov/diabetes/pdfs/data/statistics/national-diabetes-statistics-report.pdf

³ https://www.visionweb.com/content/consumers/dev_consumerarticles.jsp?RID=16

⁴ <https://nei.nih.gov/health/diabetic/retinopathy>

⁵ <https://jamanetwork.com/journals/jama/fullarticle/186384>



© 2020 Midmark Corporation, Miamisburg, Ohio USA
MKT00160

HEINE is a registered trademark of HEINE Optotechnik GmbH + Co. KG

midmark
Designing better care.

