



Midmark Corp.  
 675 Heathrow Drive  
 Lincolnshire, Illinois 60069-4206  
 1-847-415-9800  
 Fax 847-415-9801  
 midmark.com

## Progeny® Imaging Software

System Requirements:	Recommended:	Minimum:
<b>Windows Operating System</b>	<ul style="list-style-type: none"> <li>Windows 10 Pro or Enterprise (32 &amp; 64)</li> <li>Windows 2008 R2 Server &amp; above</li> </ul>	<ul style="list-style-type: none"> <li>Windows 8 Pro or Enterprise Editions (32 &amp; 64 bit)</li> <li>Windows 7® Ultimate, Enterprise or Professional (32 &amp; 64 bit)</li> </ul>
<b>Apple Support</b>	<ul style="list-style-type: none"> <li>Intel-based Apple hardware running native Windows Platform using Boot Camp</li> </ul>	
<b>Processor</b>	Intel i3 (or greater)	Intel® (or x86 compatible) Single Core 2 GHz or better
<b>Memory</b>	8 GB RAM (or greater)	4 GB RAM
<b>Storage</b>	250GB Hard Drive (or greater)	80GB Hard Drive
<b>Video</b>	32 bit, 1920 X 1080 Resolution Capable	32 bit, 1024 X 768 Resolution Capable
<b>Display</b>	1920 X 1080, 32 true bit color	1024 X 768, 32 true bit color
<b>USB Ports</b>	High Speed USB 3.0	High Speed USB 2.0

A back-up storage device is recommended (DVD, External Hard Drive, other)

Note: Windows 10 Home, Windows 8, Windows 8 RT, Windows 7® Home Premium, Windows Vista® Home (Basic, Premium, Ultimate), Microsoft XP® (Home or Pro) and Tablet PC Editions are NOT recommended. Windows Small Business Server is NOT supported. Any Windows Server acting as a Domain Controller is NOT supported.

For additional information, contact Midmark Imaging Technical Support at 847-415-9800 option 3.

Effective May 19, 2017

ML-00009 Rev. J1