



Midmark Corp.  
675 Heathrow Drive  
Lincolnshire, Illinois 60069-4206  
1-847-415-9800  
Fax 847-415-9801  
midmark.com

## Progeny® Imaging Software

System Requirements:	Recommended:	Minimum:
<b>Windows Operating System</b>	<ul style="list-style-type: none"><li>• Windows 10 Pro or Enterprise (32 &amp; 64)</li><li>• Windows 2008 R2 Server &amp; above</li></ul>	<ul style="list-style-type: none"><li>• Windows 8 Pro or Enterprise Editions (32 &amp; 64 bit)</li><li>• Windows 7® Ultimate, Enterprise or Professional (32 &amp; 64 bit)</li></ul>
<b>Apple Support</b>	<ul style="list-style-type: none"><li>• Intel-based Apple hardware running native Windows Platform using Boot Camp</li></ul>	
<b>Processor</b>	Intel i3 (or greater)	Intel® (or x86 compatible) Single Core 2 GHz or better
<b>Memory</b>	8 GB RAM (or greater)	4 GB RAM
<b>Storage</b>	250GB Hard Drive (or greater)	80GB Hard Drive
<b>Video</b>	32 bit, 1920 X 1080 Resolution Capable	32 bit, 1024 X 768 Resolution Capable
<b>Display</b>	1920 X 1080, 32 true bit color	1024 X 768, 32 true bit color
<b>USB Ports</b>	High Speed USB 3.0	High Speed USB 2.0

A back-up storage device is recommended (DVD, External Hard Drive, other)

Note: Windows 10 Home, Windows 8, Windows 8 RT, Windows 7® Home Premium, Windows Vista® Home (Basic, Premium, Ultimate), Microsoft XP® (Home or Pro) and Tablet PC Editions are NOT recommended. Windows Small Business Server is NOT supported. Any Windows Server acting as a Domain Controller is NOT supported.

For additional information, contact Midmark Imaging Technical Support at 847-415-9800 option 3.

Note:

If you decide to run additional resources or services, the Progeny Imaging Software performance may be impacted and could result in not meeting your expectations. Additionally, as technology changes it may present an increasing demand for processing power or other resources; as such our system requirements are subject to change over time.

Effective December 13, 2017

ML-00009 Rev. K1