



Seating

Midmark Veterinary
Products and Services



Let our seating products provide the comfort and support you deserve.

Midmark veterinary seating products were designed for unsurpassed durability, even under the most rigorous clinical conditions. With many choices, colors and accessories, we are certain to have a stool to meet your specific needs.

Classic Series Stools feature a polished aluminum/black accented 5-leg, star-base with heavy duty hooded casters, up to a 350 lb load capacity and a variety of backrest and lift options.



Midmark 272 Air Lift Stool



Midmark 277 Air Lift Stool



Midmark Dental Stools are not only engineered to deliver premium styling and exceptional comfort, they are also ergonomically refined to address musculoskeletal pain common when performing procedures. Promoting a more balanced posture, the seating system offers a contoured seat with raised front to improve weight distribution and overall comfort. Comfort is also enhanced with waterfall front edges that reduce pressure on the back of the thighs. Multiple seat sizes and seat heights, paired with unsurpassed adjustability, allow you to better personalize your seating solution.



Ergonomically Refined

Midmark stools are designed to allow you to maintain a proper working posture. The adjustable height allows users to keep feet flat on the floor. And lumbar support helps maintain the natural C-curve of the spine.

Doctor and Technician Stools

- Contoured seat with raised front to improve weight distribution and overall comfort
- Waterfall front edges reduce pressure on the back of the thighs
- One lever adjustment for locking seat and backrest tilt
- Ratcheted backrest for intuitive height adjustment without knobs
- Personalize your stool by selecting from multiple stool sizes, cylinder heights 15.25" - 27.5" and upholstery levels and colors
- Integrated footrest that automatically adjusts with seat height



Classic Series Stool Specifications

Seat height:
274 - 278
18" - 24.5" (45.7 - 62.3 cm)
279
19" - 25.5" (48.3 - 64.8 cm)
Seat diameter:
16.25" (41.3 cm)

Seat cushion thickness:
3.25" (8.3 cm)
Back cushion dimensions:
275, 277
14" W x 7.5" H (35.6 x 19.1 cm)
Caster base diameter:
21.5" (54.6 cm)
Weight capacity:
274, 276
350 lb (158.8 kg)
275, 277 - 279
225 lb (115.7 kg)

Value Series Stool Specifications

Seat height:
270 - 273
18" - 24.5" (45.7 - 62.2 cm)
Seat diameter:
16.25" (41.3 cm)
Seat cushion thickness:
3.25" (8.3 cm)
Back cushion dimensions:
271 - 273
14" W x 7.5" H (35.6 x 19.1 cm)
Caster base diameter:
23" (58.4 cm)
Weight capacity:
270, 272
300 lb (136.1 kg)
271, 273
225 lb (102 kg)

Classic and Value Series Options

Locking casters
Soft rubber casters
Glides
6" Height extension
Foot ring

680 Side Chair Specifications

Seat height:
17" (43.2 cm)
Seat cushion dimensions:
18" W x 19.5" D
(45.7 x 49.5 cm)
Seat cushion thickness:
2" (5.1 cm)
Back cushion dimensions:
18" W x 13.5" H (45.7 x 34.3 cm)
Inside arm width: 17.5" (44.5 cm)
on model with arms
Outside arm width: 23.5"
(59.7 cm)
on model with arms
Weight capacity:
400 lb (181.4 kg)

Doctor and Technician Stool Specifications








Seat size:
17.5" x 16.5" (44.5 x 41.9 cm)
standard
16.25" x 16" (41.3 x 40.6 cm)
small
Cylinder height:
16.5" - 22" (41.9 x 55.9 cm)
standard
15.25" - 18.3" (38.7 x 46.5 cm)
short
Short technician's stool cylinder height is not available with a footrest.

Available Colors








We have tried our very best to accurately represent product colors in this brochure. However, due to the nature of printed media, the color of actual products may slightly differ from how they appear on your piece.

To request color swatches, please contact the Midmark Veterinary Customer Experience Department at 1-800-MIDMARK.

Midmark Premium Solid Upholstery Selections



-  Blueberry
-  Cashmere Blue
-  Silver Sage
-  Terra Cotta
-  Fossil Grey
-  Tea Green
-  Black

Midmark Ultra-Premium Upholstery Selections*









-  Ultra-Pacific Blue
-  Ultra-Sky
-  Ultra-Spruce
-  Ultra-Canyon
-  Ultra-Arctic
-  Ultra-Olivine
-  Ultra-Raven Wing

* There will be an additional charge for Midmark Ultra-Premium.

Midmark Premium Pattern Upholstery Selections

-  Belagio
-  Capri
-  Siena
-  Treviso
-  Firenze
-  Tuscany

Ritter Upholstery Selections

-  Navy
-  Dusty Blue
-  Moss
-  Clay
-  Pebble Grey
-  Shadow
-  Perfect Plum
-  Black

Better Patient Care. Better Business.™

   midmarkanimalhealth.com

Manufactured by Midmark Corporation, Versailles, OH.

© 2011 Midmark Corporation • Products subject to improvement changes without notice • Litho in U.S.A. • 007-0778-00 Rev. B1 (1/16)
Midmark is an ISO 13485 Certified Company.



Because we care.