



DENTAL
Midmark Choice Package



The operator you've been imagining is within your reach.

Choose solutions designed to help you provide quality, safe and comfortable treatment while enabling more positive and meaningful interactions between you, your patients and your team. With the Midmark® Choice Package, you can make your dreams a reality without compromising quality for affordability.

CREATE YOUR PACKAGE

The Midmark Choice Package includes:* UltraComfort® or UltraTrim® chair, Asepsis 21® delivery unit, Asepsis 21 assistant's unit, LED lighting and stools

PAGE

4 CHAIRS AND SEATING



Midmark UltraComfort and UltraTrim dental chairs provide quality and comfort at an exceptional price. Midmark dental stools are also available as part of the Choice Package.

PAGE

6 DELIVERY SYSTEMS AND ACCESSORIES



The Asepsis 21 Delivery System offers superior reliability and peace of mind with removable covers that make the product easy to clean, promoting asepsis at the point of care. Procenter delivery systems are available as well, subject to standard pricing.

PAGE

8 LIGHTING AND MONITOR MOUNTS



The Midmark LED operator light offers true, reduced-shadow, white light that helps to properly identify shades, see details and make diagnoses with confidence.

PAGE

11 ORDERING GUIDE



Use this guide of the top-selling items in our operator product lines to find the solutions that are right for you.

*See Ordering Guide for details.

Chairs and Seating



Midmark dental chairs are designed to offer unrestricted clinician access to the patient, premium patient support and unparalleled comfort. The Ultra series features our patented integrated-armrest design for superior patient support and unencumbered entry and exit.

OPERATORY EXAMPLE #1

Chair: UltraComfort dental chair
Unit: Asepsis 21® rear delivery (doctor and assistant)
Light: Left/right mounted LED
Cabinetry: Artizan® Expressions base and overhead treatment station, side cabinet, dispenser and upper storage
Stools: Operator's and assistant's stools

COLORS

Cabinet: Silver Riftwood
Countertops: Hazelnut Corian® Solid Surface
Upholstery: Midnight Ultraleather®
Wall: Argos SW7065
Resin: Araña Coco
Floor: Mannington™ Lakeview Adura® Max (Rapid)

Designed to Be Better

Choose unrestricted clinical access, premium patient support, comfort and so much more.

01 The Midmark UltraComfort® dental chair offers a supportive, comfortable chairback with a cast-aluminum structure and integrated armrests that cradle the patient's elbows and inner shoulders for a secure, relaxed experience.



Both the UltraComfort and UltraTrim chairs come standard with programmable chair controls and a hydraulic drive system for smooth, quiet, precise operation. Available in your choice of limited-seams standard, plush Ultraleather® and a limited-seams version of Ultraleather® upholstery.

02 The Midmark UltraTrim® dental chair offers exceptional clinician access to the patient with a chairback that is narrower and 40% thinner than the UltraComfort® chairback. The patented chairback design features a cast-aluminum structure with integrated armrests that support unencumbered patient entry and exit.

03 The Midmark Operator's Stool and Assistant's Stool are engineered to deliver premium styling and exceptional comfort with ergonomic refinement that addresses the musculoskeletal strains common in dentistry.

The seat of the Operator's Stool can be tilted forward to help maintain the natural curve of the lower back, while the contoured backrest comfortably supports the lower back.

The Assistant's Stool offers an optional ratcheted backrest, a ratchet-release body support with height adjustment and an integrated footrest that automatically adjusts with the seat height.



Delivery Systems and Accessories

Midmark delivery systems are designed from the inside out to complement the way you work. The Asepsis 21® Delivery System provides quality and reliability in a sleek, compact, value-priced package. Or upgrade to our Procenter delivery system for additional features such as integrated controls and switches ergonomically designed for continuous use.

At the heart of each Midmark delivery system is our patented Kink Valve® technology, providing an independent control system for each handpiece and exceptional reliability—tested to over 1 million cycles. Midmark delivery systems are available in multiple configurations with accessories and options to meet your unique needs.



OPERATORY EXAMPLE #2

Chair: UltraTrim®
Unit: Asepsis 21 left/right delivery unit, chair-mounted rear assistant's unit
Light: Ceiling-mounted LED
Cabinetry: Synthesis® treatment station, side cabinet and dispenser



COLORS

Cabinet: Timber
Countertops: Lava Rock Corian® Solid Surface
Upholstery: Cyan Ultraleather® Limited Seams
Wall: Argos SW7065
Floor: BerryAlloc™ DreamClick Pro (River Oak Greige)




01 The Asepsis 21 Rear Delivery Unit offers the quality and reliability of the Asepsis 21 Delivery Unit with the flexibility of a rear delivery configuration.

02 The Asepsis 21 Delivery Unit provides air and water controls protected from contamination beneath an easily removable top.



03 The Procenter Assistant's Unit brings the thoughtful design, attention to detail and quality construction of the Procenter Delivery Unit to a configuration ideal for the dental assistant.

04 The Procenter Delivery Unit provides integrated controls and switches designed to deliver precision, flexibility and reliability.



ASEPSIS 21 DELIVERY UNIT: Available in cart, console-mount, or left/right delivery configurations, the Asepsis 21 is designed for flexibility. Continental-style delivery is also available.

PROCENTER ASSISTANT'S UNIT Available in chair, cabinet and wall-mounted options, the Procenter is designed to complement the way you work. Procenter delivery systems are available as part of the Choice Package, subject to standard pricing.

Lighting and Monitor Mounts

The advanced design of the Midmark LED light delivers the true, reduced-shadow white light you need to properly match shades, identify details and make diagnoses with confidence. The automatic on/off function and upfront touchpad control system make the light exceptionally easy to use.

Midmark lighting is available in various mounting configurations to meet the specific needs of your practice. Mounting options include ceiling-mounted, track-mounted, wall- or cabinet-mounted, left/right-mounted and track-mounted with monitor.



01 The LED Dental Operatory Light maximizes energy efficiency with only 22 watts of power consumption while producing a bright white light at 27,000, 21,000 or 15,000 Lux. The operatory light delivers true, reduced-shadow light that allows you to properly identify shades, see details and evaluate tissues while maintaining a soft, feathered edge to reduce eye strain associated with sharp differences in contrast. With no need for a noisy cooling fan, you can count on quiet operation.

The ergonomic design offers complete control with pull, tilt and twist options, as well as easy access to the touchpad. Conveniently select from three light intensities in operate mode, or easily switch to composite mode for filtered light that prevents premature curing.

02 The Track Light Monitor from Midmark combines excellent dental lighting with highly flexible monitor positioning, enhancing the patient experience and aiding patient acceptance of procedures.



OPERATORY EXAMPLE #3

Chair: UltraTrim®
Unit: Asepsis 21® left/right delivery unit, chair-mounted rear assistant's unit
Light: LED Track Light and Monitor
X-ray: Preva mounted in central station
Cabinetry: Artizan® Expressions treatment station, central station, and side cabinet, dispenser and upper storage

COLORS

Cabinet: Black Riftwood
Countertop: Canyon Corian® Solid Surface
Upholstery: Charcoal Ultraleather®
Wall: Argos SW7065
Floor: Mannington™ Restoration Wide Plank - Woodland Maple (Fawn)

What You Get with Midmark



EFFICIENCY. LONGEVITY. ATTAINABILITY.

When you choose Midmark, you get so much more than a product. You get a US-based partner with decades of hands-on manufacturing expertise, inspired design innovation and landmark thought leadership in workflow efficiency. Midmark continues its long tradition of improving patient care by infusing function-forward design, relentless quality and thoughtful technology into every product. And with an extensive warranty and dedicated technical support, we're right by your side the moment you need us.

WARRANTY AND SERVICE

We stand by the quality and reliability of our products. That is why we offer a standard five-year parts-only warranty for operatory products purchased through an authorized Midmark dealer. In addition to warranty protection, Midmark has a large team of expert customer experience and technical support representatives that will assist you with any product questions regarding installation and service.

RESOURCES



OPERATORY CABINETRY

Ask for our Operatory Cabinetry brochure to add beautiful, durable cabinetry to your operatory room design.

midmark.com/cabinetry



DENTAL DESIGN TOOL

Our easy-to-use, life-like digital tool helps you find your individual style and visualize what your exam room and operatory will look like.

apps.midmark.com/dental-design-tool



WEBSITE

For more information, visit midmark.com/choicepackage

Ask about our complete line of Midmark Choice Package items.

MIDMARK CHOICE PACKAGE

Ask for your custom quote with package savings today.

ORDERING GUIDE

Use this guide to find the products that are right for your practice. To receive the Choice Package discount, each package must include an Ultra-series chair and Asepsis 21 delivery unit at a minimum. The Asepsis 21 assistant's unit, LED light and stools are optional components for the discount. If you add any of these optional items to your package, you can then select the Procenter delivery unit instead of the Asepsis 21 delivery unit and still receive the Choice Package discount on all package items, paying standard price for the Procenter unit. This offer may not be combined with any other offers, promotions or contract pricing. For a complete list of Midmark solutions, call 1.800.MIDMARK or contact your Midmark sales representative.

CHAIRS***

153592-001	UltraComfort® Chair, Chair Only
153592-002	UltraComfort Chair, Console Mount
153592-003	UltraComfort Chair, Left/Right Mount
153758-001	UltraTrim® Chair, Chair Only
153758-002	UltraTrim Chair, Console Mount
153758-003	UltraTrim Chair, Left/Right Mount

DELIVERY (DOCTOR'S AND ASSISTANT'S) SYSTEMS

153600-010	Asepsis 21 Left/Right Unit
153600-011	Asepsis 21 Console Unit
153604LR	Left/Right Rear Hub Mounted Assistant's Unit
153604CE	Console Assistant's Unit
153604LRCA	Concept Left/Right Cuspidor and Assistant's Combination Unit
153603C	Cuspidor, Console Mounted
153603LR	Cuspidor, Concept Left/Right Mounted
153600-XXX	Asepsis 21 Continental*
KSRD-CM	Asepsis 21 Rear-Delivery Unit
KSRD-CMA	Asepsis 21 Rear-Delivery Unit with Assistant's Unit
153740-00X	Asepsis 21 Delivery Carts*
153600-003	Left/Right Procenter Instrument Delivery System**
153600-002	Procenter Console Mounted Delivery System Internal Umbilical**
153600-001	Procenter Console Mounted Delivery System External Umbilical (6')**
153678XX	Procenter Side Delivery Systems*/**
153659-XXX	Procenter Rear Delivery Systems*/**

MONITOR ACCESSORIES

153698	Tall Curved Arm Mount for Ultra-Series Chairs with Left/Right Units
153698-001	Tall Curved Arm Mount for Ultra-Series Chairs without Left/Right Units
153763	Short Curved Arm Mount for Ultra-Series Chairs with Left/Right Units
153763-001	Short Curved Arm Mount for Ultra-Series Chairs without Left/Right Units
153699	Tall Curved Arm for Procenter Console
153762	Short Curved Arm for Procenter Console

LED LIGHTS

153915-002	LED Left/Right Operatory Light
153901-004	LED Console Operatory Light
CEML-L	LED Ceiling-Mount Light
153964-XXX	LED Track Light and Monitor*
TRML-L	LED Track-Mount Operatory Light
WAML-L	LED Wall-Mount Operatory Light
CBML-L	LED Cabinet-Mount Light
153904-002	LED Universal-Mount Light

STOOLS***

153811	Operator's Stool
153820	Adjustable Armrests (Optional for Operator's Stool)
153812	Assistant's Stool
153822	Adjustable Lumbar Support (Optional for Assistant's Stool)

*Request specific item number based on configuration options

**Procenter Systems not eligible for Choice Discount but can be sold as part of the package

***Upholstery is eligible for Choice Package discount



Designing better care.®

midmark.com/dental

Ultraleather is a registered trademark owned by Ultrafabrics, Inc.
Corian is a DuPont™ Safety and Construction, Inc. registered trademark for its surfacing materials.
Mannington is a trademark owned by Mannington Mills, Inc.
Adura is a registered trademark owned by Mannington Mills, Inc.
Berry Alloc is a registered trademark owned by BerryAlloc, naamloze vennootschap.



CARB 93129.2 Phase 2 Compliant and TSCA Title VI Compliant

Visit midmark.com/dental/products for the most current colors and finishes available for each product. The color examples shown here and online are the best representation of the original material. Actual color may vary slightly. We strongly recommend that you contact Midmark Customer Experience at 1.800.MIDMARK to request a sample before placing your order.

Midmark is an ISO 13485 and ISO 9001 Certified Company. Certain products are not included. See the complete list at: midmark.com/ISO.

For more information, contact your Midmark dealer or call: 1.800.MIDMARK
Outside the USA call: 1.937.526.3662 or visit our website: midmark.com

© 2023 Midmark Corporation, Versailles, Ohio USA
Products subject to improvement changes without notice. Litho in USA.

007-10054-00 Rev D1 (4/23)

