


VETERINARY

Mobile Transport Unit



Ready to Roll,
**Wherever Care
Takes You**



Built for flexibility, a Mobile Transport Unit moves easily between spaces to support patient care wherever it's needed. With durable construction, insulated panels and smooth mobility, it's ready for daily use—whether for transport, temporary housing or emergency care.

- Durable stainless steel design with quiet, sound-dampening features
- Universal door swing, locking casters and full-width handle for easy positioning
- Optional Oxygen Therapy Door supports critical care needs
- Available in 36" or 48" widths

MOBILE TRANSPORT UNIT SPECIFICATIONS

MATERIALS		
Frame	1"-square 20-gauge Stainless Steel Tube Construction	
Base Framing	11-gauge Series 304 Stainless Steel	
Shell	22-gauge Type 304 Stainless Steel	
Sound-dampening Pads	76.2 x 203.2 x 1.5 mm (approx. 1/8")	
Hinges	Stainless Steel Reinforced Polymer	
Latch	Stainless Steel with Sound-dampening Design	
Casters	Reinforced Nylon with Polyurethane Locking Wheels	
	Mobile Transport Unit	XL Mobile Transport Unit
Part Number	902-0001-50	902-0000-57
Unit Dimensions	36" W x 30" H x 28¾" D	48" W x 30" H x 28¾" D
Unit Weights	119 lb	250 lb
Door Style	Standard Wire Door (swing is reversible during installation)	

CONSTRUCTION DETAILS

- One-piece stainless steel walls with smooth welds
- Interior corners are radiused to support easy cleaning and disinfection
- Insulated end panels help reduce noise and maintain internal temperature
- Full-width bar handle for two-sided access and control
- Four 5" medical-grade swivel casters, all locking

ACCESSORIES

- Card Holder
- Oxygen Therapy Door
- Quiet-Time Cover
- Record Holder
- Squeeze Restraint