

Barrier-Free® Examination Chair Model Matrix

COMPARISON CHART

	15.5" Seat Height (Uncompressed)	17" Seat Height (Uncompressed)	Powered Back	Clean Assist Roller System	Active Sensing Technology®	Lockout Hand + Foot Controls	Heated Upholstery	Integrated Digital Scale	650 lb Weight Capacity	500 lb Weight Capacity	Standard Rail Accessory Receivers	Tamper-Resistant Electrical Receptacles	Drawer Heater + Pelvic Tilt
626-001	•		•	•	•	•			•		•	•	•
626-002	•		•	•	•	•	•		•		•	•	•
626-003	•		•	•	•	•		•	•		•	•	•
626-004	•		•	•	•	•	•	•	•		•	•	•
626-005	•		•	•	•	•		•	•		•	•	•
626-006	•		•	•	•	•	•	•	•		•	•	•
229-002		•	•		•					•	•	•	•
229-003		•	•	•	•					•	•	•	•
228-001		•								•	•		
228-002		•			•					•	•	•	•
228-003		•		•	•					•	•	•	•

For product inquiries, please contact the Midmark Medical Customer Experience Department at 1.844.856.1233, option 1.
For more information, call 1.800.MIDMARK or visit: [midmark.com](https://www.midmark.com)