

Stainless Steel Cages

- All 20 Gauge 304 Stainless Steel Construction—23% thicker than industry average
- Heavy-duty reversible door
- One-piece stainless steel hinge plate—Multiple attachment points for stability—Prevents sagging
- 3/8" diameter frame
- 1/4" diameter crossmembers
- One-piece stainless steel corner/trim strips—Enhances strength and stability—Easier to clean
- 3/16" diameter vertical rods
- Full 1" radius corners
- Heavy-gauge latch—Padlock or door clip—Extremely quiet—Easily cleaned
- Formed Threading—Will not cross thread or strip out—Prevents self-Locking stainless steel screws from loosening with everyday use
- Easier to clean, less labor cost—Larger radius corners—Fluid retaining lip—Allows for a cleaner kennel environment
- Strong door latches and assemblies—Secure, yet easy to open and close




Single Doors

- 160-1824-00 Single Door 18"w x 24"h x 28"d
- 160-2418-00 Single Door 24"w x 18"h x 28"d
- 160-2424-00 Single Door 24"w x 24"h x 28"d
- 160-2430-00 Single Door 24"w x 30"h x 28"d
- 160-3024-00 Single Door 30"w x 24"h x 28"d
- 160-3030-00 Single Door 30"w x 30"h x 28"d
- 160-3630-00 Single Door 36"w x 30"h x 28"d

Double Doors


- 160-4830-00 Double Door 48"w x 30"h x 28"d
- 160-4830-01 Double Door w/ Spring-loaded Latch 48"w x 30"h x 28"d
- 160-4836-00 Double Door w/ Gravity Latch 48"w x 36"h x 28"d



Configure Your Cage Bank Online

Select an assembly width and browse through the options to find the configuration that fits your needs.

Popular Configurations:



BACK TO CONTENTS

