

At Midmark, we believe that who's using our technology is every bit as important as how and why it's being used. We are here to support you through these difficult times. With the influx of patients into healthcare facilities, we want to share with you the following exam room equipment best practices that could be helpful in managing current challenges.





Clean Assist Roller System

Available on select Midmark® and Ritter® Barrier-Free Examination Chairs, the Clean Assist Roller System is a built-in, retractable roller base that allows you to safely move the exam chair for cleaning and disinfecting.

[Watch Video](#)



Patient Support Rails

The Patient Support Rails exam chair accessory provides patients with a 1 ¼" diameter continuous gripping surface for entering, exiting or repositioning on the exam chair. When rotated to the back horizontal position, the support rails can be used as safety rails.

[Watch Video](#)

[Download Brochure](#)



Barrier-Free Exam Chair Low Height

Midmark and Ritter Barrier-Free Examination Chairs adjust in height 15.5 – 37 inches, improving accessibility to the exam chair and reducing the likelihood of patient and staff injury. Caregivers reported dramatically less perceived physical exertion (95% – 98% less) when using a height-adjustable exam chair versus a fixed-height table.

[Patient Safety Infographic](#)

[Staff Safety Infographic](#)



18" Low Flat Position



Midmark and Ritter Barrier-Free Examination Chairs and the Midmark 630 can be moved into a low, flat position. The low seat height and wide, flat transfer surface can enable patients to transfer onto and off of the exam or procedure chair with minimal assistance to help ensure patient and staff safety. Support rail and side panel accessories are available to help ensure the safety of patients in the supine position while resting on the exam or procedure chair.

[Patient Positioning Infographic](#)



Midmark 630 HUMANFORM[®] Procedure Chair

(Can be utilized in the exam space)

The security side panels accessory for the Midmark[®] 630 are designed to provide extra side-to-side security for the patient when in the vertical position. When in the flat position, side panels can serve as an extra work surface if desired.

[Watch Video](#)

Learn more about Midmark Medical products, [here](#).