



MEDICAL

Sphygmomanometers



Working toward better care.



As healthcare continues to shift, providers are rethinking the importance of using reliable instrument technologies for the early detection and observation of diseases. That's why we partnered with HEINE®, a world leader in high-quality primary diagnostic instruments, to bring you best-in-class, German-engineered instruments.



GAMMA™ Shockproof
Sphygmomanometers
are designed to provide
reliable, accurate measurement
values with intuitive operation.



GAMMA™ G5

The core of the GAMMA G5 is the quality and durability of the housing and spoon made of strong, impact-resistant, two-component thermoplastic. The integration of an intricate, shock-absorbing system protects the manometer.



GAMMA™ XXL LF

The classic GAMMA XXL LF is made without natural rubber latex, making it the ideal sphygmomanometer in most clinical settings. The oversized, adjustable scale is easy to read, even from a distance.



CONFIGURATIONS

Sphygmomanometers

GAMMA™ XXL LF Wall Kit

M-000-09-323-166
Includes GAMMA XXL LF sphygmomanometer, universal mounting plate, adult (29 – 41 cm) and adult large (35 – 47 cm) twin-tube cuffs



GAMMA™ G5 Hand-Held Kit

M-000-09-230-166
Includes GAMMA G5 sphygmomanometer, adult (29 – 41 cm) and adult large (35 – 47 cm) single-tube cuffs



Kit Options

Wall Set Kits for GAMMA XXL LF

002-11176-01
002-11177-01
002-11178-01
(Device combination of wallboard, GAMMA XXL LF, tip dispenser and Heine EN200 of their choice)




A-095-12-930-166
Includes GAMMA XXL LF + EN200 Wall Transformer with BETA® 400 LED Otoscope and BETA® 200 LED Ophthalmoscope + Tip Dispenser



Accessories

Wallboard

X-095-12-007-166
Wallboard for device combination EN 200 + GAMMA XXL LF
Sphygmomanometer + Tip Dispenser – no instruments



Instrument Panel

I2X.XXX
"34 x 20" Accommodates instrument configurations including: GAMMA XXL LF, EN200 and Tip Dispenser



Cuffs for GAMMA XXL LF

Name	Size (cm)	Circumference (cm)	Part Number
Infant	6 x 28	8 - 13	M-000-09-701-166
Child	10 x 34	13 – 20	M-000-09-703-166
Adult Small	14 x 53	20 – 29	M-000-09-705-166
Adult	14 x 58	29 – 41	M-000-09-629-166
Adult Large	17 x 70	35 – 47	M-000-09-631-166
Thigh	20 x 86	> 41	M-000-09-625-166

Cuffs for GAMMA G5

Name	Size (cm)	Circumference (cm)	Part Number
Infant	6 x 28	8 - 13	M-000-09-700-166
Child	10 x 34	13 – 20	M-000-09-702-166
Adult Small	14 x 53	20 – 29	M-000-09-704-166
Adult	14 x 58	29 – 41	M-000-09-619-166
Adult Large	17 x 70	35 – 47	M-000-09-706-166
Thigh	20 x 86	> 41	M-000-09-615-166

Specifications

Environmental conditions for operation:
50°F to 104°F
(+10°C to +40°C)
15% to 85% rel. humidity
700 hPa to 1060 hPa

Environmental conditions for storage:
41°F to 113°F
(+5°C to +45°C)
45% to 80% rel. humidity
500 hPa to 1060 hPa

Environmental conditions for transport:
-4°F to 122°F
(-20°C to +50°C)
45% to 80% rel. humidity
500 hPa to 1060 hPa

Measuring range:
0-300 mmHg

Accuracy of the pressure indicator:
± 3 mmHg



Designing better care.®

HEINE is an ISO 13485 Certified Company.

For more information, contact your authorized
Midmark distributor or call: 1.800.MIDMARK

Outside the USA call: 1.937.526.3662 or
visit our website: midmark.com

© 2023 Midmark Corporation, Versailles, Ohio USA
Products subject to improvement changes without notice.

HEINE and BETA are registered trademarks of
HEINE Optotechnik GmbH + Co. KG.

007-10392-00 Rev. C1 (6/25)

