

# Key Takeaways:

## Midmark HEINE LED Otoscopies + Ophthalmoscopes

Great outcomes start with great care powered by great equipment. With Midmark HEINE diagnostic instruments, clinicians can see clearly, diagnose with confidence and provide patients with the care they deserve. Engineered for proper visualization, this equipment provides superior optics and magnification, high-quality illumination and color rendering, and a dustproof design.



- 1 Superior Optics + Magnification**

The BETA® 200 LED Ophthalmoscope is the only ophthalmoscope to use glare-free technology called the Aspherical Optical System (AOS). AOS redirects corneal reflections away from the examiner's path of sight, providing a crisp, clear view of the entire fundus.

The BETA® 400 LED Otoscope provides 4.2x magnification—the maximum possible magnification needed to obtain a clear view while maintaining focus within the ear canal.
- 2 LED<sup>HO</sup> Technology**

HEINE was the first company to replicate halogen color-rendering index (CRI) properties and improve them using LED technology. The CRI indicates how colors are reproduced in comparison with daylight. Midmark HEINE diagnostic instruments have exceptionally high CRI values, ensuring that our customers can see colors as they truly are, contributing to a more accurate diagnosis and better outcomes.
- 3 Dustproof Design**

Sealed equipment heads in the BETA 200 LED Ophthalmoscopes protect the optical components, ensuring maintenance-free operation.

**Better view. Better care. Better outcomes.**  
To learn more, visit: [midmark.com/seeingisbelieving](https://midmark.com/seeingisbelieving)

HEINE and BETA are registered trademarks of HEINE Optotechnik GmbH + Co. KG.

MKT-00162

© 2020 Midmark Corporation, Miamisburg, Ohio USA

