

# Dental/Surgery Flat-Top Table

Designed  
for patient +  
care team  
comfort + safety



The Midmark Dental/Surgery Flat-Top table offers unencumbered access to the patient during procedures such as dentistry. Enhanced ergonomic and hygienic features promote an optimal care environment for patients and caregivers.

- Durable metal grate and built-in drainage basin help ensure fluids efficiently drain away from the patient.
- Raised marine edge prevents fluids from dripping onto the floor.
- Removable 3-quart fluid receptacle is easy to clean and disinfect.
- 50" and 60" table lengths accommodate a variety of procedure spaces.
- Adjustable table height and tilt support proper ergonomics in both electric and hydraulic models.
- Optional thermal heating helps maintain the patient's body temperature.
- Optional heavy-duty casters allow you to effortlessly move the table to where it is needed.



## CREATE YOUR OWN BEST-PRACTICE DENTAL SUITE

Combine the Midmark Dental/Surgery Flat-Top Table with our dental delivery system, digital dental radiography, anesthesia, multiparameter monitoring and Synthesis® Cabinetry.

# MIDMARK DENTAL/SURGERY FLAT-TOP TABLE SPECIFICATIONS



**Length:** 50" and 60" models

**Width:** 22" on both table models

**Height:** Adjustable from 29 ½" – 41 ½" on electric model and 29 ½" – 35" on hydraulic model


**Fluid Receptacle:** 96 oz at 0 tilt

**Optional:** Heated or Non-Heated Top

**Accessory:** Heavy-duty Caster Set

## ORDERING GUIDE

Item Number	Description:
100-2091-00	Dental/Surgery Flat-Top Table, 60", Electric, Heated
100-2092-00	Dental/Surgery Flat-Top Table, 60", Electric, No Heat
100-2093-00	Dental/Surgery Flat-Top Table, 50", Electric, Heated
100-2094-00	Dental/Surgery Flat-Top Table, 50", Electric, No Heat
100-2095-00	Dental/Surgery Flat-Top Table, 60", Hydraulic, Heated
100-2096-00	Dental/Surgery Flat-Top Table, 60", Hydraulic, No Heat
100-2097-00	Dental/Surgery Flat-Top Table, 50", Hydraulic, Heated
100-2098-00	Dental/Surgery Flat-Top Table, 50", Hydraulic, No Heat



For more information, contact your Midmark dealer or call: 1.800.MIDMARK  
Outside the USA call: 1.937.526.3662 or visit our website: midmark.com

© 2024 Midmark Corporation, Versailles, Ohio USA  
Products subject to improvement changes without notice. Litho in USA.

007-10464-00 Rev B1 (5/24)

