

Dental/Surgery Flat-Top Table

Designed for
comfort + safety



The Midmark Dental/Surgery Flat-Top Table is designed for unencumbered patient access, with ergonomic and hygienic features that support both patient care and caregiver comfort.

- Seamless top with rounded tabletop corners support patient and staff comfort.
- Larger drain basin and metal grate provide more usable workspace.
- Built-in drainage directs fluids into a removable 3-quart receptacle—no plumbing required.
- Raised marine edge helps prevent fluids from dripping onto the floor.
- 50" and 60" table lengths accommodate a variety of procedure spaces.
- Adjustable table height and tilt support proper ergonomics.
- Optional thermal heating helps maintain patient body temperature.
- Optional heavy-duty casters make the table easy to reposition.



MIDMARK DENTAL/SURGERY FLAT-TOP TABLE SPECIFICATIONS



Detail	Description
Material	18 gauge Stainless Steel Top
Length	50" and 60" models
Width	22" on both table models
Height	Adjustable from 28.5"– 40.5" on electric model and 29.5"– 41.5" on hydraulic model
Tilt Range	0-37°
Length	50" and 60" models
Removable Grate	9.5" L x 17.4" W
Optional	Heated or Non-Heated Top
Accessory	Heavy-duty Caster Set (adds 4" of Height)

ORDERING GUIDE

ELECTRIC COLUMN

100-2091-00 60" Surgical/ Dental Table, Heated Top, no Casters	100-2093-00 50" Surgical/ Dental Table, Heated Top, no Casters
100-2091-00-10 60" Surgical/ Dental Table, Heated Top, Casters	100-2093-00-10 50" Surgical/ Dental Table, Heated Top, Casters
100-2092-00 60" Surgical/ Dental Table, no Heat, no Casters	100-2094-00 50" Surgical/ Dental Table, no Heat, no Casters
100-2092-00-10 60" Surgical/ Dental Table, no Heat, Casters	100-2094-00-10 50" Surgical/ Dental Table, no Heat, Casters

HYDRAULIC COLUMN

100-2095-00 60" Surgical/ Dental Table, Heated Top, no Casters	100-2097-00 50" Surgical/ Dental Table, Heated Top, no Casters
100-2095-00-10 60" Surgical/ Dental Table, Heated Top, Casters	100-2097-00-10 50" Surgical/ Dental Table, Heated Top, Casters
100-2096-00 60" Surgical/ Dental Table, no Heat, no Casters	100-2098-00 50" Surgical/ Dental Table, no Heat, no Casters
100-2096-00-10 60" Surgical/ Dental Table, no Heat, Casters	100-2098-00-10 50" Surgical/ Dental Table, no Heat, Casters



CREATE YOUR OWN BEST-PRACTICE DENTAL SUITE

Combine the Midmark Dental/Surgery Flat-Top Table with our dental delivery system, digital dental radiography, anesthesia, multiparameter monitoring and Synthesis® Cabinetry.



Scan for additional
table solutions