



## OPTIMAL VIEWING. EASY DISINFECTING. ENHANCED PROTECTION.

The CANVYS® True Flat G Series 21.5" Display\* is a clinically certified (DICOM compliant) monitor made to withstand usage and cleaning procedures within a dental operator. It features a full HD display with a true flat design and a plane, non-reflecting glass surface for easy cleaning and disinfection. The CANVYS® True Flat G Series 21.5" Display uses anti-glare glass. With an IPX1 rating, it can resist water dripping vertically onto the unit, like condensation. The display's IP65-rated front is totally protected against dust and can resist low pressure jets of water from all directions. This display is LED backlit with low power consumption.

Use your monitor for patient entertainment, better chairside treatment planning and greater case acceptance.



Full HD display



Touch keys on front



True flat design with anti-glare protection glass and IP65-rated front, IPX1-rated whole unit



Plane non-reflecting glass surface for easy cleaning and disinfection



LED backlight with low power consumption



Medical certified (DICOM compliant)

## MONITOR SPECIFICATIONS

Overall Dimensions	Approx. 20.25" x 12.58"
Display Area	Approx. 18.78" x 10.62" (21.5" diagonal)
Pixels	1920 x 1080 pixels
Color	16.7 million
Backlight	LED
Interface	USB
Nominal Input Voltage	100–240 Vac
Nominal Input Frequency	50–60 Hz
Maximum Input Current	2 A
Nominal Output Voltage with Tolerance	24 VDC +/- 5%
Output Current	3.75 A

## INPUT TERMINALS

Power	DC power jack
VGA	D-SUB 15 pin connector
DVI	DVI 29 pin connector
HDMI	HDMI 19 pin connector
DP	DP 20 pin connector
USB	USB B type connector
Audio	3.5 mm 3 pole phone jack

\*Touchscreen functionality for the CANVYS True Flat G Series 21.5" display is not compatible with Midmark products.

CANVYS® is a registered trademark of Richardson Electronics, Ltd

For more information, contact your Midmark dealer or call 1.800.MIDMARK. Outside the USA, call 1.937.526.3662 or visit our website at [midmark.com](http://midmark.com)

© 2022 Midmark Corporation, Miamisburg, Ohio USA  
Products subject to improvement changes without notice. Litho in USA.  
007-10431-00 Rev A1 (8/22)

