



ENABLING HEALTHCARE FROM ANYWHERE

As healthcare becomes increasingly connected, organizations need flexible workstations that can support telehealth and other evolving care technologies. However, many healthcare organizations are ill-equipped to accommodate this shift. **We can help.** Midmark Workstations are made to order with customization that includes many standard features and accessories to assist with your unique telehealth program.

Easily incorporate technology such as cameras while maintaining an ergonomically correct working position whether seated or standing.

MONITOR ARM

- 180° of monitor rotation
- 9.5" of vertical monitor height adjustment
- 45° of monitor tilt adjustment

WORKSTATION HEIGHT

- 26" of work-surface vertical height adjustment accommodates the ergonomic height requirements of 95% of users in the clinical setting
- Adjust for a 5' person seated to a 6'4" person standing



MIDMARK® WORKSTATION ACCESSORIES

Design your workstations the way you need them with accessories such as tablet stations, single or dual monitors and AC Power Systems.

TELEHEALTH ACCESSORIES

TABLET DISPLAY HOLDER



LAPTOP DISPLAY HOLDER



DUAL MONITOR MOUNT



AC POWER SYSTEM



COILED POWER CORD



ADJUSTABLE UTILITY SHELF



VOIP SHELF



MOBILE PRINTER SHELF



Telehealth visits generally require cameras. To assist with this we offer tablet, laptop and monitor holders. Some monitors have cameras or one can easily be added. An AC Power System (includes power cord) or a power cord is recommended. Cameras, monitors, laptops and diagnostic devices not included.



For more information, visit:

[MIDMARK.COM/WORKSTATIONS](https://www.midmark.com/workstations)