






20" Extension Arm

Healthcare is changing. Are you changing with it?

The optional MIDMARK® Extension Arm is a perfect accessory that is installed directly to the workstation.

FEATURES & BENEFITS:

- 1. Intelligently Designed**
 - Provides an additional 20" of horizontal reach
 - 180° of lateral rotation which provides optimal privacy for clinician data entry
 - Internal pivot limiter option prevent over rotation
- 2. Installs in minutes**
 - Requires tightening of 2 bolts to secure
- 3. Built for Life**
 - Constructed of over 95% aluminium material, ensuring strength and durability



Model	Description	Rotation	Product Weight	Product Dimensions
9A588001	20" Extension Arm	270°	8 lbs.	22.00" x 2.00 x 5.80"