

Mobile Equipment

At Midmark, we believe that who's using our technology is every bit as important as how and why it's being used. We are here to support you through these difficult times. With the current influx of patients into healthcare facilities, we wanted to share with you some tips on how your Midmark equipment can be used to assist in the setup of pop-up triage and treatment stations, mobile field hospitals, quarantine tents and drive-up testing sites.





Midmark® Workstations

Midmark Workstations can be used in a variety of mobile and stationary applications. All workstations adjust in height from 23" to 49", accommodating the height requirements of the majority of users in the clinical setting. Locking casters can lock the wheel and swivel simultaneously for either mobile or stationary use.

Helpful Accessories:

- 01** AC Power System for up to 10 hours of power
- 02** Column + Rear Baskets for supply storage
- 03** Drawer for storage, available with a keyless locking system
- 04** Adjustable Utility Shelf for equipment or extra space
- 05** Vital Signs Accessory includes an equipment pole, BP cuff holder and vital signs equipment pole mount for mobile vital signs and IV fluids
- 06** CPU, Scale + Bar Code Holders for additional equipment
- 07** Mobile Printer Shelf for recordkeeping
- 08** Dual Monitor Mount, Laptop Display Holder, VESA Mounting Plate + Tablet Display to support mobile technology

[Download Brochure](#)

[Watch Video](#)



Midmark® Mobile Treatment Cabinets

Midmark Mobile Treatment Cabinets are a great way to get the supplies you need to the point of care, wherever that may be. Seamless panels and one-piece, molded drawer bodies with rounded corners aid in cleaning and disinfection. Keyed and keyless lock options are available for added security.

[Download Brochure](#)



Midmark® Digital Vital Signs Device

The mobile Midmark Digital Vital Signs Device can be used to capture vital signs in a variety of locations. Each device comes with a rechargeable battery or can be used with a powered workstation. Capture temperature, blood pressure, pulse and SpO₂ simultaneously. You can then import vital signs data directly to the EMR, eliminating transcription errors and reducing clicks.

[Download Brochure](#)

[Watch Video](#)

[Quick Reference Guide](#)



Midmark HEINE® Otosscopes + Ophthalmoscopes

Midmark HEINE desktop kits can be used in a variety of locations for diagnostic examinations. The rechargeable handles have an 8-hour continuous battery life and can be returned to the base to charge. Learn more at: midmark.com/seeingisbelieving

[Download Brochure](#)

[Watch Video](#)



Clean Assist Roller System

Available on select Midmark® and Ritter® Barrier-Free Examination Chairs, the Clean Assist Roller System is a built-in, retractable roller base that allows you to safely move the exam chair for cleaning and disinfecting.

[Watch Video](#)



Learn more about Midmark Medical products, [here](#).