



Recessed units not shown.

QUICKCLEAN™ ULTRASONIC CLEANER

POWERFUL CLEANING—CONSISTENT RESULTS

Providing powerful, effective cleaning and consistent results, the QuickClean™ Ultrasonic Cleaner is a piece of equipment your medical office can't afford to be without. Its cutting-edge technology can help you create a safer, more efficient work environment by decreasing worker exposure to contaminants and potential sharps injuries while reducing the time and effort needed for cleaning.

LEAP

POWERFUL TECHNOLOGY

The QuickClean™ unit cleans thoroughly and consistently throughout the entire bath using Frequency-Leap technology. This advanced software helps ensure all of your instruments are fully cleaned the first time, every time, no matter where they are placed in the basket.



SIMPLE TO OPERATE

The QuickClean™ is easy to use right out of the box, so you can get your staff up and running with minimal training time. Features like a digital LCD control panel that displays prompts, options for preset or customized settings, intelligent software that remembers the last cycle settings, and a fluid sensor that ensures the bath always has sufficient water levels help simplify your processes and increase productivity.



RANGE OF SIZES AND DESIGN OPTIONS

The QuickClean™ comes in three tabletop sizes (1.2, 3.3 and 6.6 gal.) and two recessed options (3.3 and 6.6 gal.), so you can choose the option that best fits your configuration and workflow. And the stainless steel finish and attractive design coordinates well with any décor.

150%
MORE EFFECTIVE
THAN THE SONICLEAN®

CUSTOMIZED
CYCLE TIMES FROM

1-99
MINUTES

FULL-SIZE
SAFETY BASKET
INCLUDED

3-YEAR
WARRANTY



Specifications, accessories and service part information on reverse.

SPECIFICATIONS:

QUICKCLEAN™ MODEL NUMBERS:	TABLETOP UNIT: QC1-01	TABLETOP UNIT: QC3-01 RECESSED UNIT: QC3R-01	TABLETOP UNIT: QC6-01 RECESSED UNIT: QC6R-01
Capacity:	1.2 Gal/4.5 L	3.3 Gal/12.5 L	6.6 Gal/25L
Overall Length:	13.5" (34.5 cm)	14.7" (37.5 cm)	21.6" (55 cm)
Overall Width:	10.2" (26 cm)	16.5" (42 cm)	16.5" (42 cm)
Overall Height:	10.4" (26.5 cm)	10.4" (26.5 cm)	14.3" (36.5 cm)
Tank Size:	11.8" L x 5.9" W x 5.9" H (30 x 15 x 15 cm)	13" L x 11.8" W x 5.9" H (33 x 30 x 15 cm)	19.9" L x 11.8" D x 7.8" H (50.5 x 30 x 20 cm)
Recessed Unit Cutout Size:	N/A	15.55" x 14.37" (39.5 x 36.5 cm)	22.48" x 14.37" (57.1 x 36.5 cm)
Table Top Unit Weight:	12.7 lb (5.76 kg)	18.6 lb (8.42 kg)	29.6 lb (13.44 kg)
Table Top Shipping Weight:	14.6 lb (6.64 kg)	20.9 lb (9.50 kg)	32.6 lb (14.78 kg)
Recessed Unit Weight:	N/A	22.3 lb (10.10 kg)	32.9 lb (14.92 kg)
Recessed Unit Shipping Weight:	N/A	24.6 lb (11.18 kg)	35.8 lb (16.26 kg)
Operating Frequency:	32-38 kHz	32-38 kHz	32-38 kHz
Temperature Up To:	140°F / 60°C	140°F / 60°C	140°F / 60°C
Stainless Steel Safety Basket:	Included	Included	Included
Electrical Requirments:	Input: 115VAC +/- 10%, 60 Hz, 110W, 1.0A Output: 100W	Input: 115V +/- 10%, 60 Hz, 210W, 1.8A Output: 200W	Input: 115V +/- 10%, 60 Hz, 410W, 3.6A Output: 400W
Standards:	UL 61010-1:2012; CAN/CSA C22.2 No. 61010-1:2012; EN 60601-1-2:2007; EN 61326:2013; FCC CFR 47 Part 15 part B	UL 61010-1:2012; CAN/CSA C22.2 No. 61010-1:2012; EN 60601-1-2:2007; EN 61326:2013; FCC CFR 47 Part 15 part B	UL 61010-1:2012; CAN/CSA C22.2 No. 61010-1:2012; EN 60601-1-2:2007; EN 61326:2013; FCC CFR 47 Part 15 part B
Warranty:	3 Years	3 Years	3 Years

ACCESSORIES & SERVICE PARTS:



1| BEAKER ACCESSORIES

- 9A612001 - QC1 – 2 Beaker Accessory Kit
- 9A613001 - QC3/QC3R – 4 Beaker Accessory Kit
- 9A614001 - QC6/QC6R – 6 Beaker Accessory Kit

2| SAFETY BASKETS

- 002-10007-00 – QC1 Safety Basket
- 002-10008-00 – QC3/QC3R Safety Basket
- 002-10009-00 – QC6/QC6R Safety Basket

3| CLEANING SOLUTIONS

- 9A296001 Midmark General Purpose Cleaner (32oz)
- 9A297001 Midmark Tartar and Stain Remover (32oz)
- 9A298001 Midmark Enzymatic Cleaner (32oz)

Midmark is an ISO 13485 and ISO 9001 Certified Company.
For more information or a demonstration, contact your Midmark
dealer or call: 1-800-MIDMARK Fax: 1-877-725-6495
Outside the U.S.A. call: 1-937-526-3662 Fax: 1-937-526-8214
or visit our website at midmark.com

© 2017 Midmark Corporation, Dayton, OH.
Products subject to improvement changes without notice
Litho in U.S.A. 007-10018-00 Rev. B1 (7/17)



Because we care.