

Ritter 203 & 95

S P E C I F I C A T I O N S



Ritter 203-012 Treatment Table

When it comes to caring for your patients, you need a table you can count on. The economical 203 is constructed of steel from its frame to its legs. And its flat top is designed to provide a stable, solid surface for patient manipulations or adjustments during examinations.

The metal-on-metal construction is also designed to give the 203 greater stability and structural integrity. Stirrups and lockable casters can be added as well as rounded drawers when you require additional storage options.

Have we mentioned how comfortable the 203 is? Extra padding and a built-in headrest pillow help keep patients comfortable. The Mist base is complemented by any of the 30 Ritter upholstery colors.

Features

Paper roll holder (hidden) with hold-down straps at the head and foot, designed to keep paper in place.

Ritter 203 Treatment Table Specifications

Patient Weight Capacity:
400 lb (182 kg)

Dimensions (top):
27" W x 74" L
(68.6 x 188 cm)

Height: 31" (78.7 cm)

Paper roll holder:
21" (53.3 cm)

Options

Infinitely adjustable (hideaway) stirrups from four lateral positions

Rounded drawers with color inserts on right or left side

Ritter 95 Examination Table Specifications

Patient Weight Capacity:
350 lb (159 kg)

Dimensions (top):
27" W x 72" L
(68.6 x 182.9 cm)

Height: 30" (76.2 cm)

Paper roll holder:
21" (53.3 cm)

Accessory

Vinyl upholstered shelf

Standard Upholstery Colors

- Mist-813
- Stone-814
- Restful Path-858
- Iron Ore-851
- Lunar Gray-845
- Obsidian-857
- Oasis-843
- Soothing Blue-855
- Dream-847
- Mineral-844
- Healing Waters-848
- Shaded Garden-853
- Citrus-849
- Sandy Retreat-850
- Curative Copper-846
- Cranberry-859
- Dark Linen-856
- Latte-860
- Robust Brown-852
- Deep Earth-854

UltraFree® Upholstery Colors

- UltraFree Stone-865
- UltraFree Restful Path-868
- UltraFree Obsidian-867
- UltraFree Harbor-861
- UltraFree Wheatgrass-862
- UltraFree Spice-863
- UltraFree Cranberry-869
- UltraFree Dark Linen-866
- UltraFree Latte-870
- UltraFree Branch-864

***UltraLeather® by Ultrafabrics**

CARB 93120.2 Phase 2 Compliant

Midmark is an ISO 13485 and ISO 9001 Certified Company.

For more information or a demonstration, contact your Midmark Dealer or Midmark Corporation, Dayton, Ohio.
Call: 1-800-MIDMARK
Fax: 1-800-365-8631
midmark.com



Ritter 95 Treatment Table

