



MEDICAL

Guide Book

Effective October 1, 2018

General Information

Ordering Notes.....	4
Freight and Shipping.....	4
Midmark and Ritter Medical Equipment Warranty.....	5

Midmark Chairs and Tables

Midmark 641 Procedure Chair.....	6-7
Midmark 630 HUMANFORM® Procedure Chair.....	8-10
Midmark 647 Podiatry Chair.....	11-12
Midmark 646 Podiatry Chair.....	13
Midmark 626 Barrier-Free® Examination Chair.....	14-16
Midmark 640 Pediatric Examination Table.....	17
New Upholstery Color Palette.....	18
Former Upholstery Color Palette.....	19

Ritter Chairs and Tables

Ritter 230 Procedure Chair.....	20-21
Ritter 244 Barrier-Free® Examination Table.....	22
Ritter 225 Barrier-Free® Examination Chair.....	23
Ritter 224 Barrier-Free® Examination Chair.....	24-25
Ritter 204 Manual Examination Table.....	26-27
Ritter 203 Manual Examination Table.....	28
Ritter 95 Manual Examination Table.....	28

Seating

Midmark 425 Physician Seating (Hand Operated).....	29
Midmark 427 Physician Seating (Foot Operated).....	29
Ritter 279 Physician Seating (Hand Operated).....	30
Ritter 278 Physician Seating (Foot Operated).....	30
Ritter 277 Physician Seating with Backrest.....	30
Ritter 276 Physician Seating.....	30
Ritter 275 Physician Seating with Backrest.....	31
Ritter 274 Physician Seating.....	31
Ritter 273 Physician Seating with Backrest.....	31
Ritter 272 Physician Seating.....	31
Ritter 271 Physician Seating with Backrest.....	32
Ritter 270 Physician Seating.....	32
Ritter 280 Patient Side Chair.....	33
Ritter 282 Patient Bariatric Side Chair.....	33
Midmark 680 Patient Side Chair.....	34
Ritter 281 Blood Draw Chair.....	34
Seating Upholstery Color Selections.....	35

Lighting

Ritter 255 LED Procedure Light.....	36
Ritter 253 LED Exam Light.....	37
Ritter 250 LED Exam Light.....	38

Diagnostic Products

Midmark Diagnostic Medical Products Warranty.....	39
Midmark Diagnostic Medical Products Ordering Information.....	40
Midmark HEINE Diagnostic Instruments Limited Warranty.....	41
Midmark Diagnostic Medical Products Support Information.....	42
Midmark Digital Vital Signs and Digital Vital Signs Monitor and Supplies.....	43-47
Midmark Digital Spirometer, Software and Supplies.....	48
Midmark Digital ECG, Software and Supplies.....	49
Midmark Digital Holter, Software and Supplies.....	50-51
Midmark HEINE Diagnostic Instrument Wall Kits and Accessories.....	52
Midmark HEINE Diagnostic Instrument Desktop Kits and Accessories.....	53
Midmark HEINE Diagnostic Instrument Replacement Parts and Supplies.....	54
Midmark HEINE Gamma® XXL LF Sphygmomanometer Kit and Accessories.....	55
Midmark HEINE Gamma® G5 Sphygmomanometer Kit and Accessories.....	56
Midmark Diagnostic Cart Accessories and Supplies.....	57

Midmark Care Exchange® Workstations

Midmark Workstation Warranty.....	58
Midmark 6211 Workstation.....	59
Midmark 6212 Workstation.....	60
Midmark 6213 Workstation.....	61
Midmark 6214 Workstation.....	62
Midmark 6215 Workstation.....	63
Midmark 6217 Workstation.....	64
Midmark 6218 Workstation.....	65
Midmark 6219 Workstation.....	66
Midmark 6281 Workstation.....	67
Midmark 6282 Workstation.....	67
Midmark 6283 Workstation.....	68
Midmark 6231 Workstation.....	69
Midmark 6201 Workstation.....	70
Midmark 6202 Workstation.....	71
Midmark 6203 Workstation.....	72
Midmark 6204 Workstation.....	73
Midmark 6205 Workstation.....	74
Midmark 6206 Workstation.....	75
Midmark 6207 Workstation.....	76
Midmark 6208 Workstation.....	77
Midmark 6261 Workstation.....	78
Midmark 6263 Workstation.....	79
Midmark 6271 Workstation.....	80
Midmark 6272 Workstation.....	81
Midmark 6273 Workstation.....	81
Midmark 6274 Workstation.....	82
Midmark 6275 Workstation.....	82
Midmark 6276 Workstation.....	83

Midmark Synthesis® Cabinetry Collection

Midmark Synthesis Cabinetry Warranty.....	84
Midmark Synthesis Cabinetry Color Selections.....	85
Midmark Synthesis Exam Room Packages E1-E7.....	86-92
Midmark Synthesis Mobile Treatment Cabinets M2; M4-M5.....	93-95
Midmark Synthesis Glove/ Towel/ Cup Dispenser D1-D4.....	96
Midmark Synthesis Dressing Nook N1-N2.....	97
Midmark Synthesis Waste Receptacle W1.....	98
Midmark Synthesis Instrument Panel I1-I2.....	99
Ritter 261 Waste Receptacle.....	100

Instrument Processing

Ritter M11 UltraClave® Automatic Sterilizer.....	101
Ritter M9 UltraClave® Automatic Sterilizer.....	101
Ritter M9D Manual Sterilizer.....	102
Ritter M3 UltraFast® Automatic Sterilizer.....	102
Sterilizer Accessories.....	103
VistaCool™ Direct-to-Drain System.....	104
QuickClean® Ultrasonic Cleaner QC6/QC6R.....	105
QuickClean® Ultrasonic Cleaner QC3/QC3R.....	105
QuickClean® Ultrasonic Cleaner QC1.....	106
Ultrasonic Parts and Accessories.....	106

Midmark Customer Experience

To place an order, email: Orders-Medical@midmark.com. For order status, requests and information, contact Medical Customer Experience at:

1-844-856-1233, Option 1

Email: allcustservice-medical@midmark.com

Midmark Technical Service

For service part identification, technical/operational questions, warranty and product issues/complaints, contact Technical Service at:

1-844-856-1230.

Sterilizer & Ultrasonic Products: sterilizers.ts@midmark.com

Diagnostic Products: techsupport@midmark.com

Procedure & Exam Chairs, Seating, Lighting, Diagnostic Instrument Wall/Desk Kits, Oscopes,

Ophthalmoscopes & Sphygmomanometers: examroomeq.ts@midmark.com

Medical Cabinetry, Carts & Workstations: casework.ts@midmark.com

Special Notes

1. **ORDERING**—Please provide the following information on all orders:

- Purchase order number
 - Bill to and ship to addresses
 - Product model and/or part number(s)
 - Description including upholstery color and/or panel/countertop color
 - Quantity needed
 - Pricing
 - Required shipping date
 - Special shipping requirements (see details below)
- MINIMUM ORDER — \$100.

DROP SHIP REQUIREMENTS

- Please specify if the customer has a loading dock or has special shipping requirements.
- Receiving contact name, email and phone number is required.

2. **SPECIAL DESIGN** — Orders for equipment of special design, construction, upholstery or special finish are not subject to cancellation after your acceptance of our acknowledgment.

3. **SPECIAL SHIPPING REQUIREMENTS** — These services are available to enhance delivery to your customers. Please call for pricing.

- **CALL BEFORE:** Carrier will call delivery contact to provide a drop-off timeframe once they have possession of the product from Midmark and are at the final terminal.
- **LIFTGATE:** Carrier uses liftgate/ramp to get product to ground level when facility does not have a dock. Facility is responsible for taking product once at ground level.
- **INSIDE DELIVERY:** Carrier unloads product and delivers inside the first set of doors only, provided entrance to facility can accommodate freight as packaged (including skids) and there are no steps.
- **MIDMARK DELIVERY SERVICE:** Carrier sets product in place per the delivery contact. Carrier is responsible for assembling product and removal of all packaging, trash and skids from the delivery site (custom casework does not include installation).

4. **FREIGHT DAMAGE** — Merchandise arriving with external damage should be so described on the bill of lading. In case of concealed damage, carrier must be notified within 15 days after receipt and claims for damage instituted.

5. **ORDER CHANGES** — Any changes made to an order within 48 hours of the scheduled ship date may affect the scheduled ship date, requiring Midmark to reschedule.

6. **ORDER CANCELLATION** — Any order, if cancelled within ten days before the scheduled shipping date, will be charged a 20% cancellation charge.

7. **RETURN OF MERCHANDISE** — Authorization must be obtained for the return of any item. If return is authorized, transportation charges must be prepaid. Risk of loss with respect to returned goods shall remain with purchaser until receipt by Midmark. Items returned for credit are subject to a charge of 25% or more to cover inspection, damage in transit, reconditioning and repacking. Items may be returned within 6 months (except service parts within 90 days) from the original ship date if in a good, unused condition. RMA # should be written on the outside of the packaging. Diagnostic products may be returned within 30 days if in a good, unused condition and in original sealed packaging.

Note: Certain items are non-returnable: Specials, Synthesis® Cabinetry, Midmark Workstations, Discontinued Items or Parts, Custom or Non-standard Colors, Service Upholstery and Service Parts if value is less than \$50. Diagnostic Consumables e.g. Electrodes, Prep Kits, Pneumotach Mouthpieces, Temperature Probe Covers, Nose Clips, Thermal Paper, All Disposable Items and Diagnostic Items value less than \$50 per.

8. **SPECIAL ORDER REQUIREMENTS** — Any order requiring any special deviation from standard products will carry the following disclaimer:

We hereby accept and acknowledge the above captioned order, provided however, because of the special nature of this order in that this equipment is being supplied upon special design specifications, special materials and/or field modifications, as per your request, the standard warranties of Midmark and Ritter on its products are not applicable to this order, and this acceptance and acknowledgment is conditioned upon your indemnification and promise to hold Midmark and Ritter, its officers, employees and agents, harmless in respect to any and all warranty claims, damages, both direct and consequential, expenses, costs and any liability arising from the use of or relating to this special equipment.

9. **IMPROVE PRODUCT QUALITY** — Midmark Corporation reserves the right to change design and construction without notice.

Midmark and Ritter Medical Equipment (Tables, Chairs, Physician Seating, Lights, Care Exchange Workstations, and Sterilizers) Limited Warranty

SCOPE OF WARRANTY Midmark Corporation ("Midmark") warrants to the original retail purchaser that it will, at Midmark's option, repair or replace components of the domestic and international medical products manufactured by Midmark (except for components not warranted under "Exclusions") that are defective in material or workmanship under normal use and service. The sole remedy under this limited warranty is the repair or replacement, at Midmark's option, of the applicable components. This limited warranty shall only apply to defects that are reported to Midmark within the applicable warranty period and which are determined to exist upon examination by Midmark. This warranty extends only to the original retail purchaser of a product and is not transferable or assignable. Replacement components or products may be used and/or refurbished components or products, provided they are of like quality and specifications as new components or products.

Midmark warrants to the original retail purchaser that during the applicable warranty period it will repair or replace software contained within the products manufactured by Midmark (except for those not warranted under "Exclusions") if: (1) the media on which the software is furnished exhibits defects in material or workmanship under normal use; or (2) the software does not substantially conform to its published specifications.

APPLICABLE WARRANTY PERIOD The applicable warranty period, measured from the date of invoice to the original retail purchaser of the product and shall be one (1) year for all warranted products and components.

OBTAINING WARRANTY SERVICE Warranty service must be obtained through either Midmark or an authorized dealer in the Midmark product line for which warranty service is requested. Midmark may be contacted for warranty service inquiries or issues via email at www.midmark.com; by phone at 1-800-MIDMARK or by mail to Midmark Corporation, 60 Vista Drive, Versailles, Ohio 45380.

It is the retail purchaser's obligation to arrange for delivery of a product to Midmark or one of its authorized dealers for warranty service, which delivery shall be at retail purchaser's expense. It is also the retail purchaser's obligation to comply with the warranty service instructions provided either by Midmark or its authorized dealer. The retail purchaser must provide Midmark with completed warranty registration information within thirty (30) days after purchase in order to obtain the benefits of this limited warranty.

EXCLUSIONS: This limited warranty does not cover and Midmark shall not be liable for the following:

- (1) Defects, damage or other conditions caused, in whole or in part, by misuse, abuse, negligence, alteration, accident, freight damage, negligent storage, tampering or failure to seek and obtain repair or replacement in a timely manner;
- (2) Products which are not installed, used, and properly cleaned and maintained as required or recommended in the Midmark "Installation" and/or "Installation/Operation Manual" for the applicable product, including the specified structural and operational environmental conditions and electrical requirements;
- (3) Products considered to be of a consumable or sterile nature;
- (4) Accessories or parts not manufactured by Midmark;
- (5) Charges by anyone for adjustments, repairs, replacement parts, installation or other work performed upon or in connection with such products which are not expressly authorized in writing in advance by Midmark;
- (6) Costs and expenses of routine maintenance and cleaning; and
- (7) Representations and warranties made by any person or entity other than Midmark.
- (8) Matching of color, grain or texture except to commercially acceptable standards;
- (9) Changes in color caused by natural or artificial light;
- (10) Custom manufactured products;
- (11) Alterations or modifications to the product by any person or entity other than Midmark; and
- (12) Products that would otherwise be covered under this limited warranty, but are acquired: (i) from a person or entity that is not Midmark or one of its authorized dealers; or (ii) from a Midmark dealer that is not authorized to sell the product at issue in the geographic territory where the purchaser is located, or is not authorized to sell the product at issue within the medical market.

SOFTWARE; WITH RESPECT TO SOFTWARE THAT IS A PRODUCT OR COMPONENT THEREOF, MIDMARK DOES NOT WARRANT THAT THE SOFTWARE: (1) IS ERROR FREE; (2) CAN BE USED WITHOUT PROBLEMS OR INTERRUPTIONS; OR (3) IS FREE FROM VULNERABILITY TO INTRUSION OR ATTACK BY VIRUSES OR OTHER METHODS.

EXCLUSIVE REMEDY; CONSEQUENTIAL DAMAGES DISCLAIMER MIDMARK'S ONLY OBLIGATION UNDER THIS LIMITED WARRANTY IS THE REPAIR OR REPLACEMENT OF DEFECTIVE PARTS. MIDMARK SHALL NOT BE LIABLE FOR AND HEREBY DISCLAIMS ANY DIRECT, SPECIAL, INDIRECT, INCIDENTAL, EXEMPLARY OR CONSEQUENTIAL DAMAGES OR DELAYS, INCLUDING, BUT NOT LIMITED TO, DAMAGES FOR LOSS OF PROFITS OR INCOME, LOSS OF USE, LOSS OF DATA, DOWNTIME, COVER AND EMPLOYEE OR INDEPENDENT CONTRACTOR WAGES, PAYMENTS AND BENEFITS. THIS DISCLAIMER SHALL SURVIVE ANY FAILURE OR ASSERTED FAILURE OF THE ESSENTIAL PURPOSE OF THIS LIMITED WARRANTY OR ITS REMEDIES SPECIFIED HEREIN.

WARRANTY DISCLAIMER THIS WARRANTY IS MIDMARK'S ONLY WARRANTY AND IS IN LIEU OF ALL OTHER WARRANTIES, EXPRESS OR IMPLIED. MIDMARK MAKES NO IMPLIED WARRANTIES OF ANY KIND INCLUDING ANY IMPLIED WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. THIS WARRANTY IS LIMITED TO THE REPAIR OR REPLACEMENT OF DEFECTIVE PARTS.

STATUTE OF LIMITATIONS No action may be brought against Midmark for breach of this limited warranty, an implied warranty, if any, or for any other claim arising out of or relating to the products, more than ninety (90) days following expiration of the limited warranty period.

NO AUTHORIZATION No person or firm is authorized to create or approve for Midmark any other obligation or liability in connection with the products.

*Warranty information included in product manuals



Midmark 641

Procedure Chair

- Power height, back, foot and tilt
- 450-pound patient weight capacity
- 19" to 40" seat height adjustment (non-rotation models)
- 22.5" to 43.5" seat height adjustment (rotation models)
- Perimeter contouring for maximum access to patient
- Flip-out foot extension
- Hand and foot control standard
- Paper roll holder standard
- Choice of Premium/Premium UltraFree* seamless or Ultra-Premium plush upholstery
- Meets UL 60601-1 and CAN/CSA C22.2 No. 601.1-M90 standards
- Table weight: 450 lb
- Shipping weight: 500 lb

Model

Description

Standard Models

641-002
002-0966-00-XXX² Procedure Chair Non-Programmable w/o receptacle, w/o top
Upholstery top - Premium

641-002
002-0968-00-XXX² Procedure Chair Non-Programmable w/o receptacle, w/o top
Upholstery top - Ultra-Premium

641-002
002-10147-XXX¹ Procedure Chair Non-Programmable w/o receptacle, w/o top
641 Premium UltraFree* 24" Upholstery

641-003
002-0966-00-XXX² Procedure Chair Programmable w/receptacle, w/o top
Upholstery top - Premium

641-003
002-0968-00-XXX² Procedure Chair Programmable w/receptacle, w/o top
Upholstery top - Ultra-Premium

641-003
002-10147-XXX¹ Procedure Chair Programmable w/o receptacle, w/o top
641 Premium UltraFree* 24" Upholstery

Rotation Models

641-004
002-0966-00-XXX² Procedure Chair Non-Programmable w/rotation, w/o top
Upholstery top - Premium

641-004
002-0968-00-XXX² Procedure Chair Non-Programmable w/rotation, w/o top
Upholstery top - Ultra-Premium

641-004
002-10147-XXX¹ Procedure Chair Non-Programmable w/rotation, w/o top
641 Premium UltraFree* 24" Upholstery

641-005
002-0966-00-XXX² Procedure Chair Programmable w/rotation, w/o top
Upholstery top - Premium

641-005
002-0968-00-XXX² Procedure Chair Programmable w/rotation, w/o top
Upholstery top - Ultra-Premium

641-005
002-10147-00-XXX¹ Procedure Chair Programmable w/rotation, w/o top
641 Premium UltraFree* 24" Upholstery

Accessories	Description
Arms/Armboards	
9A179002-XX ^{2,7}	Fixed armboard
9A396001-XXX ^{2,5,7}	Chair arms
9A398001-XXX ²	Chair arms, Ultra-Premium
9A81001-XXX ^{2,7}	Articulating armboard
9A82001-XXX ^{2,7}	Hand surgery armboard
Casters	
9A330001 ³	Caster base, swivel wheels
9A365001 ³	Casters
Controls	
9A408001 ⁴	Wireless control, non-programmable, linear-style, plant installed
9A408002 ⁴	Wireless control, programmable, linear-style, plant installed
9A408003 ⁴	Wireless control, non-programmable, linear-style, field installed
9A408004 ⁴	Wireless control, programmable, linear-style, field installed
Headrests	
9A378001-XXX ^{2,7}	Flat headrest, Premium upholstery
9A380001-XXX ²	Flat headrest, Ultra-Premium upholstery
9A382001-XXX ^{2,7}	U-shaped headrest, Premium upholstery
9A384001-XXX	U-shaped headrest, Ultra-Premium
9A386001-XXX ^{2,7}	Oval-shaped headrest, Premium upholstery
9A388001-XXX ²	Oval shaped headrest, Ultra-Premium
9A390001-XXX ^{2,7}	Rectangular-shaped headrest, Premium upholstery
9A392001-XXX ²	Rectangular-shaped headrest, Ultra-Premium upholstery
9A394001-XXX ^{2,7}	Magnetic headrest, Premium upholstery
9A395001-XXX	Magnetic headrest, Ultra-Premium upholstery
Instrument Trays/Work Shelves/IV Poles	
9A427001	Swing arm instrument tray, fixed top
9A427002	Swing arm instrument tray, extension top
9A428001	Swing arm PC work shelf
9A267003 ⁶	Instrument tray, double arm
9A77001 ⁶	IV pole
Miscellaneous	
9A152010	3' power cord
9A205003	Base rails, plant installed
9A205004	Base rails, field installed
9A406001	Security straps, plant installed
9A406002	Security straps, field installed
9A78001-XXX ^{2,6,7}	Vision block screen
9A416001	Foot section plastic cover
002-0966-00-999	Upholstery top - Premium Special Color
002-0968-00-999	Upholstery top - Ultra-Premium Special Color
002-0967-00-XXX ²	CAL 133 Upholstery top - Premium
002-0969-00-XXX ²	CAL 133 Upholstery top - Ultra-Premium
002-0967-00-999	CAL 133 Upholstery top - Premium Special Color
002-0969-00-999	CAL 133 Upholstery top - Ultra-Premium Special Color
002-10148-XXX	CAL 133 Upholstery Top - Premium UltraFree*

² - Former color palette available until 10/31/18. See page 19

³ - Not available on rotation models

⁴ - Accessory supplied with hand and foot control (batteries included)

⁵ - Interference with armboard accessories will likely occur

⁶ - Requires 9A20500X base rail

⁷ - Additional charge for Special, Premium UltraFree* or Ultra Premium Upholstery Color for Accessory

*ULTRALEATHER® by Ultrafabrics®



Midmark 630

HUMANFORM® Procedure Chair

- Power height, back, foot and tilt
- Motion Profiling® designed to provide more natural, human-like movement that works with the patient and provider for a clear view of the procedure site and improved comfort
- Active Sensing Technology® (AST) is designed to actively protect patient and provider by sensing and avoiding harmful impact during movement, such as objects being caught under the seat and footrest
- Barrier-Free design
- 650-pound patient weight capacity
- 18" to 40" seat height adjustment (non-rotational models)
- 21.5" to 43.5" seat height adjustment (rotational models)
- Hand and foot controls standard
- Available in 28" or 32" upholstery widths
- Meets UL 60601-1 and CAN/CSA C22.2 No.601.1-M90 standards
- Table weight: 450 lb
- Shipping weight: 525 lb

Model Non-Programmable

Description

630-010 002-10108-XXX ²	Procedure Chair Non-Programmable w/receptacle, w/o top 630 Premium 28" Upholstery
630-010 002-10114-XXX ²	Procedure Chair Non-Programmable w/receptacle, w/o top 630 Premium UltraFree* 28" Upholstery
630-011 002-10108-XXX ²	Procedure Chair Non-Programmable wireless, w/receptacle, w/o top 630 Premium 28" Upholstery
630-011 002-10114-XXX ²	Procedure Chair Non-Programmable wireless, w/receptacle, w/o top 630 Premium UltraFree* 28" Upholstery

Non-Programmable with Rotation

630-012 002-10108-XXX ²	Procedure Chair Non-Programmable w/rotation, w/o top 630 Premium 28" Upholstery
630-012 002-10114-XXX ²	Procedure Chair Non-Programmable w/rotation, w/o top 630 Premium UltraFree* 28" Upholstery
630-013 002-10108-XXX ²	Procedure Chair Non-Programmable wireless, w/rotation, w/o top 630 Premium 28" Upholstery
630-013 002-10114-XXX ²	Procedure Chair Non-Programmable wireless, w/rotation, w/o top 630 Premium UltraFree* 28" Upholstery

Model Programmable	Description
630-020 002-10108-XXX ²	Procedure Chair Programmable w/receptacle, w/o top 630 Premium 28" Upholstery
630-020 002-10114-XXX ²	Procedure Chair Programmable, w/receptacle, w/o top 630 Premium UltraFree* 28" Upholstery
630-021 002-10108-XXX ²	Procedure Chair Programmable wireless, w/receptacle, w/o top 630 Premium 28" Upholstery
630-021 002-10114-XXX ²	Procedure Chair Programmable wireless, w/receptacle, w/o top 630 Premium UltraFree* 28" Upholstery

Programmable with Rotation

630-022 002-10108-XXX ²	Procedure Chair Programmable, w/rotation, w/o top 630 Premium 28" Upholstery
630-022 002-10114-XXX ²	Procedure Chair Programmable, w/rotation, w/o top 630 Premium UltraFree* 28" Upholstery
630-023 002-10108-XXX ²	Procedure Chair Programmable wireless, w/rotation, w/o top 630 Premium 28" Upholstery
630-023 002-10114-XXX ²	Procedure Chair Programmable wireless, w/rotation, w/o top 630 Premium UltraFree* 28" Upholstery

Upholstery Options

002-10108-999 002-10114-999	Upholstery top - 28" wide Premium Special Color Upholstery top - 28" wide UltraFree* Special Color
002-10112-XXX ² 002-10112-999	Upholstery top - 32" wide Premium Upholstery top - 32" wide Premium Special Color
002-10116-XXX ² 002-10116-999	Upholstery top - 32" wide UltraFree* Upholstery top - 32" wide UltraFree* Special Color
002-10109-XXX ² 002-10109-999	CAL 133 Upholstery top - 28" wide Premium CAL 133 Upholstery top - 28" wide Premium Special Color
002-10115-XXX ² 002-10115-999	CAL 133 Upholstery top - 28" wide UltraFree* CAL 133 Upholstery top - 28" wide UltraFree* Special Color
002-10113-XXX ² 002-10113-999	CAL 133 Upholstery top - 32" wide Premium CAL 133 Upholstery top - 32" wide Premium Special Color
002-10117-XXX ² 002-10117-999	CAL 133 Upholstery top - 32" wide UltraFree* CAL 133 Upholstery top - 32" wide UltraFree* Special Color

*ULTRALEATHER® by Ultrafabrics®

Accessories

Arm Systems/Armboards

9A550001-XXX ^{2,7}	Chair arms, 28" wide Premium upholstery
9A550003-XXX	Chair arms, 28" wide Ultra-Premium upholstery
9A550005-XXX ^{2,7}	Chair arms, 32" wide Premium upholstery
9A550007-XXX ²	Chair arms, 32" wide Ultra-Premium upholstery
9A179002-XXX ^{2,6,7}	Fixed Armboard
9A81001-XXX ^{2,6,7}	Articulating Armboard
9A82001-XXX ^{2,6,7}	Hand surgery Armboard

Controls (additional)

9A552001	Wired hand control, non-programmable, field installed
9A553001	Wired foot control, non-programmable, field installed
9A552002	Wired hand control, programmable, field installed
9A553002	Wired foot control, programmable, field installed

Headrests

9A496001-XXX ^{2,7}	Rectangular-shaped headrest, Premium upholstery
9A497001-XXX ²	Rectangular-shaped headrest, Ultra-Premium upholstery
9A500001-XXX ^{2,7}	Oval-shaped headrest, Premium upholstery
9A501001-XXX ²	Oval-shaped headrest, Ultra-Premium upholstery
9A504001-XXX ^{2,7}	Round headrest, Premium upholstery
9A505001-XXX ²	Round headrest, Ultra-Premium upholstery
9A508001-XXX ^{2,7}	Flat headrest, Premium upholstery
9A509001-XXX ²	Flat headrest, Ultra-Premium upholstery
9A512001-XXX ^{2,7}	U-shaped headrest, Premium upholstery
9A513001-XXX ²	U-shaped headrest, Ultra-Premium upholstery
9A516001-XXX ^{2,7}	Magnetic headrest, Premium upholstery
9A517001-XXX ²	Magnetic headrest, Ultra-Premium upholstery

Instrument Trays/Work Shelves/IV Poles

9A363001	PC work surface
9A267003 ⁴	Instrument tray, double arm
9A77001 ⁴	IV pole

Knee Crutches

9A410001 ⁵	Standard knee crutches (crutches only)
9A410002 ⁵	Standard knee crutch system
9A411001 ⁵	Articulating knee crutches (crutches only)
9A411002 ⁵	Articulating knee crutch system
9A412001 ⁵	Knee crutch mounting brackets (brackets only)

Accessory Rails

9A521001	Seat Accessory Rails, Surgical Size, 28" wide upholstery, plant installed
9A522001	Seat Accessory Rails, Surgical Size, 32" wide upholstery, plant installed
9A567001	Seat Accessory Rails, Surgical Size, 28" wide upholstery, field installed
9A568001	Seat Accessory Rails, Surgical Size, 32" wide upholstery, field installed
9A557001	Accessory Rails, 28" wide upholstery, field installed
9A558001	Accessory Rails, Surgical Size, 28" wide upholstery, field installed
9A557002	Accessory Rails, 32" wide upholstery, field installed
9A558002	Accessory Rails, Surgical Size, 32" wide upholstery, field installed

Miscellaneous

9A554001 ³	Castors - set of 4
9A152010	3' power cord
9A157003 ⁴	Security side panels
9A559001 ⁶	Security straps, field installed
9A555001	Gynecology debris tray accessory
9A556001	Gynecology debris tray (tray only)
9A78001-XXX ^{2,7}	Anesthesia/Vision block screen
AD 1779 ³	Tamper-resistant duplex receptacle

² - Former color palette available until 10/31/18. See page 19

³ - Not available on rotation models

⁴ - Requires (1) 9A557001/9A557002/9A558001/9A558002 accessory rails

⁵ - Not compatible with 9A521001, 9A522001, 9A567001 or 9A568001

⁶ - Requires (1) 9A557001/9A557002/9A558001/9A558002 accessory rails and (1) 9A521001/9A522001 seat rails

⁷ - Additional charge for Special, UltraFree* or Ultra Premium Upholstery Color for Accessory

*ULTRALEATHER® by Ultrafabrics®



Midmark 647

Podiatry Chair

- Power height, tilt and back
- 450-pound patient weight capacity
- 19" to 40" seat height adjustment (non-rotation models)
- 22.5" to 43.5" seat height adjustment (rotation models)
- 9" foot section extension
- Sealed foot control standard
- Integrated touch controls standard
- Chair arms standard
- Choice of Premium seamless or Ultra-Premium plush upholstery
- Protective foot cover standard with upholstery
- Meets UL 60601-1 and CAN/CSA C22.2 No. 601.1-M90 standards
- Table weight: 433 lb (non-rotation models)
- Shipping weight: 515 lb (non-rotation models)
- Table weight: 539 lb (rotation models)
- Shipping weight: 621 lb (rotation models)

Model	Description
647-001 002-1119-XXX ²	Podiatry Chair Non-programmable w/receptacles, w/o top Upholstery top - Premium
647-001 002-1121-XXX ²	Podiatry Chair Non-programmable w/receptacles, w/o top Upholstery top - Ultra-Premium
647-001 002-10129-XXX	Podiatry Chair Non-programmable w/receptacles, w/o top 646/647 Premium UltraFree* 28" Upholstery
647-002 002-1119-XXX ²	Podiatry Chair Programmable w/receptacles, w/o top Upholstery top - Premium
647-002 002-1121-XXX ²	Podiatry Chair Programmable w/receptacles, w/o top Upholstery top - Ultra-Premium
647-002 002-10129-XXX	Podiatry Chair Programmable w/receptacles, w/o top 646/647 Premium UltraFree* 28" Upholstery
647-003 002-1119-XXX ²	Podiatry Chair Non-programmable w/rotation, w/o receptacles & top Upholstery top - Premium
647-003 002-1121-XXX ²	Podiatry Chair Non-programmable w/rotation, w/o receptacles & top Upholstery top - Ultra-Premium
647-003 002-10129-XXX	Podiatry Chair Non-programmable w/rotation, w/o receptacles & top 646/647 Premium UltraFree* 28" Upholstery
647-004 002-1119-XXX ²	Podiatry Chair Programmable w/rotation, w/o receptacles & top Upholstery top - Premium
647-004 002-1121-XXX ²	Podiatry Chair Programmable w/rotation, w/o receptacles & top Upholstery top - Ultra-Premium
647-004 002-10129-XXX	Podiatry Chair Programmable w/rotation, w/o receptacles & top 646/647 Premium UltraFree* 28" Upholstery

*ULTRALEATHER® by Ultrafabrics®

Accessories	Description
Instrument Trays/Work Shelves/IV Poles	
9A427001	Swing arm instrument tray, fixed top
9A427002	Swing arm instrument tray, extension top
9A428001	Swing arm PC work shelf
9A432001	Instrument tray, swing arm (tray only)
9A77001 ⁴	IV Pole
Casters	
9A330001 ³	Caster base, swivel wheels
9A365001 ³	Casters
Miscellaneous	
9A418001	Wireless foot control, non-programmable, linear-style, plant installed
9A418002	Wireless foot control, non-programmable, linear-style, field installed
9A418003	Wireless foot control, programmable, linear-style, plant installed
9A418004	Wireless foot control, programmable, linear-style, field installed
9A422001-XXX ^{2,5}	Vision block screen
9A420001	Plastic protective foot cover
9A423001	Paper roll holder
9A205003	Base rails, plant installed
9A205004	Base rails, field installed
9A152010	3' power cord
002-1119-999	Upholstery top - Premium Special Color
002-1121-999	Upholstery top - Ultra-Premium Special Color
002-1120-XXX ²	CAL 133 Upholstery top - Premium
002-1122-XXX ²	CAL 133 Upholstery top - Ultra-Premium
002-1120-999	CAL 133 Upholstery top - Premium Special Color
002-1122-999	CAL 133 Upholstery top - Ultra-Premium Special Color
002-10130-XXX	CAL 133 Upholstery top - UltraFree*

² - Former color palette available until 10/31/18. See page 19.

³ - Not available on rotation models

⁴ - Requires 9A20500X base rail

⁵ - Additional charge for Special, Ultra Premium or UltraFree* Accessory

*ULTRALEATHER® by Ultrafabrics®



Midmark 646 Podiatry Chair

- Power tilt and back
- 450-pound patient weight capacity
- 21" fixed seat height
- 9" foot section extension
- Sealed foot control standard
- Integrated touch controls standard
- Chair arms standard
- Choice of Premium Seamless or UltraFree* upholstery
- Protective foot cover standard with upholstery
- Meets UL 60601-1
- Table weight: 323 lb
- Shipping weight: 405 lb

***Not available for
sale in Canada**

Model	Description
646-001 002-1119-XXX ²	Podiatry Chair Non-programmable w/receptacles, w/o top Upholstery top - Premium
646-001 002-1121-XXX ²	Podiatry Chair Non-programmable w/receptacles, w/o top Upholstery top - Ultra-Premium
646-001 002-10129-XXX	Podiatry Chair Non-programmable w/ receptacles, w/o top 646/647 Premium UltraFree* 28in Upholstery

Accessories	Description
Instrument Trays/Work Shelves/IV Poles	
9A427001	Swing arm instrument tray, fixed top
9A427002	Swing arm instrument tray, extension top
9A428001	Swing arm PC work shelf
9A432001	Instrument tray, swing arm (tray only)
9A77001 ³	IV Pole
Casters	
9A330001	Caster base, swivel wheels
9A365001	Casters
Miscellaneous	
9A422001-XXX ^{2,4}	Vision block screen
9A420001	Plastic protective foot cover
9A205003	Base rails, plant installed
9A205004	Base rails, field installed
9A152010	3' power cord
002-1119-999	Upholstery top - Premium Special Color
002-1121-999	Upholstery top - Ultra-Premium Special Color
002-1120-XXX ²	CAL 133 Upholstery top - Premium
002-1122-XXX ²	CAL 133 Upholstery top - Ultra-Premium
002-1120-999	CAL 133 Upholstery top - Premium Special Color
002-1122-999	CAL 133 Upholstery top - Ultra-Premium Special Color

² - Former color palette available until 10/31/18. See page 19.

³ - Requires 9A20500X base rail

⁴ - Additional charge for Special, Ultra Premium or UltraFree* Upholstery Color for Accessory



Midmark 626

Barrier-Free® Examination Chair

- 650-pound weight capacity
- Low height of 15.5" (uncompressed) to 37" high
- Clean Assist Roller System
- Power base and back
- Hand and foot controls
- Digital Scale (-003 and -004 models)
- Heated upholstery options (-002 and -004 models)
- Quick Exam® and Return to Home
- Active Sensing Technology®
- Mist and Stone table base color lockout of hand and foot controls
- Table weight: 435 lb
- Shipping weight: 485 lb

Model	Description
Standard Models	
626-001	Barrier-Free Examination Chair w/Clean Assist Roller System, receptacles, drawer heater, pelvic tilt, wired controls, w/o top
002-2000-XXX ¹	626 Premium 28" Upholstery
626-001	Barrier-Free Examination Chair w/Clean Assist Roller System, receptacles, drawer heater, pelvic tilt, wired controls, w/o top
002-2002-XXX ¹	626 Premium 32" Upholstery
626-001	Barrier-Free Examination Chair w/Clean Assist Roller System, receptacles, drawer heater, pelvic tilt, wired controls, w/o top
002-2004-XXX ¹	626 Premium UltraFree* 28" Upholstery
626-002	Barrier-Free Examination Chair w/ Clean Assist Roller System, receptacles, drawer heater, pelvic tilt, wired controls, heated upholstery options, w/o top
002-2006-XXX ¹	626 Premium 28" Heated Upholstery
626-002	Barrier-Free Examination Chair w/ Clean Assist Roller System, receptacles, drawer heater, pelvic tilt, wired controls, heated upholstery options, w/o top
002-10201-XXX ¹	626 Premium 32" Heated Upholstery

*ULTRALEATHER® by Ultrafabrics®

Model	Description
Standard Models	
626-002	Barrier-Free® Examination Chair w/ Clean Assist Roller System, receptacles, drawer heater, pelvic tilt, wired controls, heated upholstery options, w/o top
002-2011-XXX ¹	626 Premium UltraFree Heated 28" Upholstery
626-003	Barrier-Free Examination Chair w/ Clean Assist Roller System, Digital Scale, receptacles, drawer heater, pelvic tilt, wireless controls, w/o top
002-2000-XXX ¹	626 Premium 28" Upholstery
626-003	Barrier-Free Examination Chair w/ Clean Assist Roller System, Digital Scale, receptacles, drawer heater, pelvic tilt, wireless controls, w/o top
002-2002-XXX ¹	626 Premium 32" Upholstery
626-003	Barrier-Free Examination Chair w/ Clean Assist Roller System, Digital Scale, receptacles, drawer heater, pelvic tilt, wireless controls, w/o top
002-2004-XXX ¹	626 Premium UltraFree* 28" Upholstery
626-004	Barrier-Free Examination Chair w/ Clean Assist Roller System, Digital Scale, heated upholstery options, receptacles, drawer heater, pelvic tilt, wireless controls, w/o top
002-2006-XXX ¹	626 Premium 28" Heated Upholstery
626-004	Barrier-Free Examination Chair w/ Clean Assist Roller System, Digital Scale, heated upholstery options, receptacles, drawer heater, pelvic tilt, wireless controls, w/o top
002-10201-XXX ¹	626 Premium 32" Heated Upholstery
626-004	Barrier-Free Examination Chair w/ Clean Assist Roller System, Digital Scale, heated upholstery options, receptacles, drawer heater, pelvic tilt, wireless controls, w/o top
002-2011-XXX ¹	626 Premium UltraFree* Heated 28" Upholstery

*ULTRALEATHER® by Ultrafabrics®

Arm Systems

Part #	Description
9A600002 ²	Patient Support Rails Plus, field installed
9A600001 ²	Patient Support Rails, field installed
9A582001	Patient Support Rails, receiver extension - pair; field install for 32" wide upholstery

Knee Crutches

9A643001	Standard knee crutch kit w/ acc. receiver, field install
9A645001	Articulating knee crutch kit w/ acc. receiver, field install
9A641001	Mounting brackets, hardware w/ acc. receiver, field install

Vitals Signs Accessories

9A478001	Midmark Digital Scale serial cable for Midmark Digital Vital Signs, 6 foot length, coiled
9A478002	Midmark Digital Scale serial cable for Midmark Digital Vital Signs, 15 foot length, straight
9A478003	Midmark Digital Scale serial cable for Midmark Digital Vital Signs, 30 foot length, straight
9A478004	Midmark Digital Scale serial cable for Midmark Digital Vital Signs, 50 foot length, straight
9A629001	Midmark Digital Vital Signs table mount with USB and serial cables
9A629002	Midmark Digital Vital Signs table mount with serial cable
9A629003	Midmark Digital Vital Signs table mount, non-IQhub®
9A479001	Welch Allyn Spot Vitals Lxi, 6 foot length coiled
9A479002	Welch Allyn Spot Vitals Lxi, 15 foot length straight
002-10064-00	626 chair connectivity kit

Miscellaneous

9A677001	Wireless hand & foot controls
9A610001	Stainless steel treatment pan
002-2015-XXX ¹	626 Premium UltraFree* 32" upholstery
002-10203-XXX ¹	626 Premium UltraFree* heated 32" upholstery
002-2003-XXX ¹	626 Premium 32" upholstery, CAL 133
002-2001-XXX ¹	626 Premium 28" upholstery, CAL133
002-2005-XXX ¹	626 Premium UltraFree* 28" upholstery, CAL133
002-2016-XXX ¹	626 Premium UltraFree* 32" upholstery, CAL133
002-2007-XXX ¹	626 Premium 28" heated upholstery, CAL133
002-10202-XXX ¹	626 Premium 32" heated upholstery, CAL133
002-2012-XXX ¹	626 Premium UltraFree* heated 28" upholstery, CAL133
002-10204-XXX ¹	626 Premium UltraFree* heated 32" upholstery, CAL133

¹See color palette on page 18

²9A582001 receiver extensions required for use with 32" wide upholstery

*ULTRALEATHER® by Ultrafabrics®



Midmark 640

Pediatric Examination Table


- Pebble grey base and panel color only
- An upholstered end panel and beveled sides designed to help keep infants from rolling off
- Flat sides with no sharp edges
- Accurate and easy-to-use, built-in scale
- 36-inch exam surface height
- Ample storage and a hidden paper roll holder
- Meets ANSI/AAMI ES 60601-1:2005 & CAN/CSA C22.2 No. 60601-1-08 standards
- 42" W x 51.5" H x 23"D
- Table weight: 170 lb
- Shipping weight: 210 lb

Model	Description
640-001-XXX ²	Pediatric Examination Table w/digital scale
640-001-999	Pediatric Examination Table w/digital scale - Special Color
Miscellaneous	
640-011-XXX ²	CAL 133 Pediatric Examination Table w/digital scale
640-012-XXX ²	CAL 133 Pediatric Examination Table with door and drawer lock, plant installed
640-011-999	CAL 133 Pediatric Examination Table w/digital scale - Special Color
640-012-999	CAL 133 Pediatric Examination Table w/door and drawer lock, plant installed - Special Color
640-002-XXX ²	Pediatric Examination Table w/door and drawer lock, plant installed
640-002-999	Pediatric Examination Table w/door and drawer lock, plant installed - Special Color

² - Former color palette available until 10/31/18. See page 19.

NEW! Color Palette

XXX - XXX - XXX



Example:
 630-010 *Chair*
 002-10108-843 *Top*
 630 HUMANFORM® *Procedure Chair - Oasis, Premium Upholstery Top*

Standard Upholstery Colors

Mist-813	Soothing Blue-855	Curative Copper-846
Stone-814	Dream-847	Cranberry-859
Restful Path-858	Mineral-844	Dark Linen-856
Iron Ore-851	Healing Waters-848	Latte-860
Lunar Gray-845	Shaded Garden-853	Robust Brown-852
Obsidian-857	Citrus-849	Deep Earth-854
Oasis-843	Sandy Retreat-850	

UltraFree* Upholstery Colors

UltraFree Stone-865	UltraFree Spice-863
UltraFree Restful Path-868	UltraFree Cranberry-869
UltraFree Obsidian-867	UltraFree Dark Linen-866
UltraFree Harbor-861	UltraFree Latte-870
UltraFree Wheatgrass-862	UltraFree Branch-864

*ULTRALEATHER® by Ultrafabrics®

Former Color Palette

Available on noted products until 10/31/18

(Not available on Midmark 626 and Ritter 225/224 Chairs)

Midmark Brand

Solid Colors

235 Blueberry
236 Silver Sage
237 Tea Green
238 Fossil Grey
239 Cashmere Blue
240 Terra Cotta
312 Black

Pattern Premium Colors

247 Belagio
248 Siena
249 Tuscany
250 Firenze
251 Capri
252 Treviso

Ultra-Premium Colors

228 Ultra-Raven Wing
241 Ultra-Pacific Blue
242 Ultra-Spruce
243 Ultra-Olivine
244 Ultra-Arctic
245 Ultra-Sky
246 Ultra-Canyon

¹See product price sheet
for Ultra-Premium color
availability and price.

Ritter Brand

216 Pebble Grey
229 Perfect Plum
230 Moss
231 Navy
232 Shadow
233 Dusty Blue
234 Clay
312 Black



Ritter 230

Procedure Chair

- Power height, back, foot and tilt
- 450-pound patient weight capacity
- 22.5" to 40" seat height adjustment (non-rotation models)
- 26" to 43.5" seat height adjustment (rotation models)
- Knee-chest configurable
- Hand control standard
- Meets UL60601-1 and CAN/CSA C22.2 No. 601.1 – M90 standards
- Table weight: 455 lb
- Shipping weight: 525 lb

Model	Description
230-001 002-0824-XXX ²	Procedure Chair w/o top 230 Seamless Upholstery top
230-001 002-10133-XXX	Procedure Chair w/o top 230 Seamless UltraFree* 28" Upholstery
230-002 002-0824-XXX ²	Procedure Chair w/receptacle, w/o top 230 Seamless Upholstery top
230-002 002-10133-XXX	Procedure Chair w/receptacle, w/o top 230 Seamless UltraFree* 28" Upholstery
230-003 002-0824-XXX ²	Procedure Chair w/rotation, w/o top 230 Seamless Upholstery top
230-003 002-10133-XXX	Procedure Chair w/rotation, w/o top 230 Seamless UltraFree* 28" Upholstery

*ULTRALEATHER® by Ultrafabrics®

Accessories	Description
Arm Systems/Armboards	
9A179002-XXX ^{2,6}	Fixed armboard
9A328001	Adjustable arm system, plant installed
9A328002	Adjustable arm system, field installed
9A350001	Adjustable arm system (receiver only), plant installed
9A350002	Adjustable arm system (receiver only), field installed
9A350003	Adjustable arm system (arms only)
9A81001-XXX ^{2,6}	Articulating armboard
9A82001-XXX ^{2,6}	Hand surgery armboard
Casters	
9A330001 ⁴	Caster base, swivel wheels
9A365001 ⁴	Casters
Controls	
9A371003	Foot control, linear style
Headrests/Footrests	
9A337001-XXX ^{2,6}	Larger 24" wide footrest
9A51001-XXX ^{2,6}	Facial pad
Instrument Trays/Work Shelves/IV Poles	
9A427001	Swing arm instrument tray, fixed top
9A427002	Swing arm instrument tray, extension top
9A428001	Swing arm PC work shelf
9A267003 ⁵	Instrument tray, double arm
9A77001 ⁵	IV pole
Knee Crutches	
9A410001	Standard knee crutches (crutches only)
9A410002	Standard knee crutch system
9A411001	Articulating knee crutches (crutches only)
9A411002	Articulating knee crutch system
9A412001	Knee crutch mounting brackets (brackets only)
Miscellaneous	
9A152010	3' Power Cord
9A157003	Security side panels
9A205003	Base rails, plant installed
9A205004	Base rails, field installed
9A214001	Security straps, plant installed
9A214002	Security straps, field installed
9A329003 ⁴	Urology drain system, plant installed
9A329004 ⁴	Urology drain system, field installed
9A78001-XXX ^{2,5,6}	Vision block screen
002-0824-999	Seamless Upholstery top Special Color
002-0825-XXX ²	CAL 133 Upholstery Top-Seamless
002-0825-999	CAL 133 Upholstery Top - Seamless Special Color
002-0951-XXX ²	CAL 133 Flat Upholstery set, 24" width
002-0951-999	CAL 133 Flat Upholstery set, 24" width Special Color
002-0950-XXX ²	Flat upholstery set, 24" width
002-0950-999	Flat Upholstery set, 24" width Special Color
002-10135-XXX	Flat UltraFree* Upholstery
002-10134-XXX	CAL 133 Seamless UltraFree Upholstery
002-10136-XXX	CAL 133 Flat UltraFree* Upholstery

² - Former color palette available until 10/31/18. See page 19.

³ - Not compatible with 9A337001

⁴ - Not available on rotation models

⁵ - Requires 9A20500X base rail

⁶ - Additional charge for Special, Ultra Premium or UltraFree* Accessory

*ULTRALEATHER® by Ultrafabrics®



Ritter 244

Barrier-Free® Examination Table

- Power height, back and tilt
- Reverse trendelenburg
- 850-pound patient weight capacity
- 18" to 34" seat height adjustment
- 32" wide seat
- Hand control standard
- Meets UL 60601-1 and CAN/CSA C22.2 No.601.1-M90 standards
- Table weight: 340 lb
- Shipping weight: 400 lb

Model

Description

244-001
002-0861-XXX²

Barrier-Free Examination Table w/o top
Seamless 32" Upholstery

244-001
002-10174-XXX

Barrier-Free Examination Table w/o top
Seamless UltraFree* 32" Upholstery

Accessories

9A627001	Adjustable arm system, field installed
9A346001	Casters
9A347001	Foot control
9A152010	3' Power Cord

Miscellaneous

002-0861-999	Seamless 32" Upholstery Top Special Color
002-1017-XXX ²	CAL 133 Seamless 32" Upholstery Top
002-1017-999	CAL 133 Seamless 32" Upholstery Top Special Color
002-10175-XXX	CAL 133 Seamless UltraFree* 32" Upholstery Top

² - Former color palette available until 10/31/18. See page 19.

*ULTRALEATHER® by Ultrafabrics®



Ritter 225

Barrier-Free® Examination Chair

- 500-pound weight capacity
- Low height of 18" to 37" high
- Power base and back
- Wired foot control
- Active Sensing Technology®
- Mist table base color table weight: 409 lb
- Shipping weight: 465 lb

Model	Description
Standard Models	
225-002	Barrier-Free Examination Chair w/ receptacles, drawer heater, pelvic tilt, w/o top
002-2009-XXX ¹	224/225 Seamless 28" Upholstery
225-002	Barrier-Free Examination Chair w/ receptacles, drawer heater, pelvic tilt, w/o top
002-2008-XXX ¹	224/225 Flat 28" Upholstery
225-002	Barrier-Free Examination Chair w/ receptacles, drawer heater, pelvic tilt, w/o top
002-2013-XXX ¹	224/225 Seamless UltraFree* 28" Upholstery
225-003	Barrier-Free Examination Chair w/Clean Assist Roller System, receptacles, drawer heater, pelvic tilt, w/o top
002-2009-XXX ¹	224/225 Seamless 28" Upholstery
225-003	Barrier-Free Examination Chair w/Clean Assist Roller System, receptacles, drawer heater, pelvic tilt, w/o top
002-2008-XXX ¹	224/225 Flat 28" Upholstery
225-003	Barrier-Free Examination Chair w/Clean Assist Roller System, receptacles, drawer heater, pelvic tilt, w/o top
002-2013-XXX ¹	224/225 Seamless UltraFree* 28" Upholstery
Arm Systems	
9A600002	Patient Support Rails Plus, Field Installed
9A600001	Patient Support Rails, Field Installed
9A615001	Assist Arms
Knee Crutches	
9A643001	Standard Knee Crutch Kit w/ Acc. Receiver, Field Install
9A645001	Articulating Knee Crutch Kit w/ Acc. Receiver, Field Install
9A641001	Mounting Brackets, Hardware w/ Acc. Receiver, Field Install
Miscellaneous	
9A610001	Stainless steel treatment pan
9A638001	Supply Assistant
9A681001	225 Hand Control Accessory
002-2010-XXX ¹	224/225 Seamless 28" Upholstery, CAL133
002-2014-XXX ¹	224/225 Seamless UltraFree* 28" Upholstery, CAL133
002-10149-XXX ¹	Premium Thick Footrest Pad, CAL 117/133
002-10214-XXX ¹	Premium UltraFree* Thick Footrest Pad, CAL 117/133

¹ - See color palette on page 18
*ULTRALEATHER® by Ultrafabrics®



Ritter 224

Barrier-Free® Examination Chair

- 500-pound weight capacity
- Low height of 18" to 37" high
- Power base
- Wired foot control
- Active Sensing Technology® (-002 and -003 models)
- Table weight: 409 lb
- Shipping weight: 465 lb

Model	Description
Standard Models	
224-001 002-2009-XXX ¹	Barrier-Free Examination Chair - excludes drawer and debris tray 224/225 Seamless 28" Upholstery
224-001 002-2008-XXX ¹	Barrier-Free Examination Chair - excludes drawer and debris tray 224/225 Flat 28" Upholstery
224-001 002-2013-XXX ¹	Barrier-Free Examination Chair - excludes drawer and debris tray 224/225 Seamless UltraFree* 28" Upholstery
224-002 002-2009-XXX ¹	Barrier-Free Examination Chair w/receptacles, drawer heater, pelvic tilt, w/o top 224/225 Seamless 28" Upholstery
224-002 002-2008-XXX ¹	Barrier-Free Examination Chair w/receptacles, drawer heater, pelvic tilt, w/o top 224/225 Flat 28" Upholstery
224-002 002-2013-XXX ¹	Barrier-Free Examination Chair w/receptacles, drawer heater, pelvic tilt, w/o top 224/225 Seamless UltraFree* 28" Upholstery
224-003 002-2009-XXX ¹	Barrier-Free Examination Chair w/Clean Assist Roller System, receptacles, drawer heater, pelvic tilt w/o top 224/225 Seamless 28" Upholstery
224-003 002-2008-XXX ¹	Barrier-Free Examination Chair w/Clean Assist Roller System, receptacles, drawer heater, pelvic tilt w/o top 224/225 Flat 28" Upholstery
224-003 002-2013-XXX ¹	Barrier-Free Examination Chair w/Clean Assist Roller System, receptacles, drawer heater, pelvic tilt w/o top 224/225 Seamless UltraFree* 28" Upholstery

Model	Description
Standard Models	
Arm Systems	
9A600002	Patient Support Rails Plus, field installed
9A600001	Patient Support Rails, field installed
9A615001	Assist Arms
Knee Crutches	
9A643001	Standard knee crutch kit w/ acc. receiver, field install
9A645001	Articulating knee crutch kit w/ acc. receiver, field install
9A641001	Mounting brackets, hardware w/ acc. receiver, field install
Miscellaneous	
9A610001	Stainless steel treatment pan
9A638001	Supply Assistant
002-2010-XXX ¹	224/225 Seamless 28" Upholstery, CAL133
002-2014-XXX ¹	224/225 Seamless UltraFree* 28" Upholstery, CAL133
002-10149-XXX ¹	Premium Thick Footrest Pad, CAL 117/133
002-10214-XXX ¹	Premium UltraFree* Thick Footrest Pad, CAL 117/133

¹ - See color palette on page 18



Ritter 204

Manual Examination Table

- 500-pound patient weight capacity
- Electrical receptacle standard
- Styled Seamless or Soft Touch stitched upholstery
- Pass-thru drawer design, no need to specify RH or LH
- Manual back with air-spring assist
- Meets UL 60601-1 and CAN/CSA C22.2 No.601.1-M90 standards
- Please specify upholstery/drawer panel color
- Mist table base color table weight: 250 lb
- Shipping weight: 300 lb

Model	Description
204-011 002-10137-XXX ²	Manual examination table w/o top 204 Soft Touch Uph, 28"
204-011 10138-XXX ²	Manual examination table w/o top 204 Upholstery top - Styled, Seamless
204-012 10137-XXX ²	Manual examination table w/drawer heater, pelvic tilt w/o top 204 Soft Touch Uph, 28"
204-012 10138-XXX ²	Manual examination table w/drawer heater, pelvic tilt w/o top 204 Upholstery top - Styled, Seamless

² - Former color palette available until 10/31/18. See page 19.

Accessories	Description
Arm/Handle Systems	
9A627001 ³	Adjustable arm system (arms only)
9A632001 ³	Assist handles (handles only), 28" wide upholstery only, field installed
9A632002	Accessory receiver only, 28" wide upholstery only, plant installed
9A632003	Accessory receiver only, 28" wide upholstery only, field installed
Knee Crutches	
9A410001	Standard knee crutch only (204/604/22X/62X/230/630) - field installed
9A410004	Standard knee crutch kit with accessory receiver (204/604/22X/62X) - field installed
9A411001	Articulating knee crutch only (204/604/22X/62X/230/630) - field installed
9A411004	Articulating knee crutch kit with accessory receiver (204/604/22X/62X) - field installed
9A412003	Mounting brackets, hardware with accessory receiver (204/604/22X/62X) - field installed
Miscellaneous	
9A104001	Urology drain pan w/hose
9A372001	Pass-through drawer stop (set of 3)
9A70000	Stainless steel treatment pan
9A629003	Midmark Digital Vital Signs table mount, Non-IQhub®
002-10143-XXX ²	204 Seamless UltraFree* Upholstery, 28"
002-10141-XXX	Flat Stitched Upholstery w/o pillow
002-10145-XXX	Flat UltraFree* Upholstery top w/o pillow
002-10193-00-XXX	Thick footrest pad kit (for Soft-Touch & Flat uph)
002-10194-00-XXX	Thick footrest pad kit (for Seamless/Ultra-Prem/UltraFree* uph)

² - Former color palette available until 10/31/18. See page 19.

³ - 9A632002 or 9A632003 accessory receiver required

*ULTRALEATHER® by Ultrafabrics®



Ritter 203

Manual Examination Table

- 400-pound patient weight capacity
- 27" W x 74" D x 31" H
- Please specify upholstery/panel color
- Panels match upholstery color
- Table weight: 135 lb
- Shipping weight: 170 lb
- 203-001 through 008 models discontinued 4/30/18

Model	Description
203-011-XXX ²	Manual examination table standard
203-012-XXX ²	Manual examination table w/drawers
203-013-XXX ²	Manual examination table w/stirrups*
203-014-XXX ²	Manual examination table w/drawers and stirrups*
Miscellaneous	
203-011-999	Manual examination table standard - special color
203-012-999	Manual examination table w/drawers - special color
203-013-999	Manual examination table w/stirrups* - special color
203-014-999	Manual examination table drawers and stirrups* - special color

* Stirrups factory installed only

² - Former color palette available until 10/31/18. See page 19.



Ritter 95

Manual Examination Table

- 350-pound patient weight capacity
- 27" W x 72" D x 30" H
- Please specify upholstery color
- Frame is a dark oak finish only
- Paper roll holder included
- Table weight: 80 lb
- Shipping weight: 100 lb

Model	Description
95-001-XXX ^{2,3}	Manual examination Table w/o shelf
95-001-999	Manual examination Table w/o shelf Special Color
Accessories	
9A364001-XXX ^{2,3}	Description Shelf (field installed only)
9A364001-999	Shelf Special Color
Options	
95-002-XXX	CAL 133 Manual Examination Table w/o shelf
95-002-999	CAL 133 Manual Examination Table w/o shelf - Special Color

² - Former color palette available until 10/31/18. See page 19.

³ - Additional charge for UltraFree* Upholstery/Accessory Color.



Midmark 425

Physician Seating (Hand Operated)

- Weight capacity: 225 lb
- Stool weight: 20 lb
- Shipping weight: 25 lb

Model

425-001
002-0273-XXX^{1,2,3}

Description

Basic stool w/o seat cushion
Seat cushion

Accessories

002-0274-XXX^{1,2,4}
002-0275-XXX^{1,2,4}
002-1546-00

Back kit
Arm kit
Lab height w/black foot ring (425 only)

Miscellaneous

002-0944-XXX^{1,2}
002-0945-XXX^{1,2}
002-0946-XXX^{1,2}

CAL 133 Seat Cushion-Premium Upholstery
CAL 133 Back Kit-Premium Upholstery
CAL 133 Arm Kit-Premium Upholstery

¹ - See color code chart on page 40

² - Additional charge for Special or Ultra Premium/UltraFree* Upholstery Color

³ - Additional charge for Special or Ultra-Premium/UltraFree* Upholstery Color for Accessory



Midmark 427

Physician Seating (Foot Operated)

- Weight capacity: 225 lb
- Stool weight: 20 lb
- Shipping weight: 25 lb

Model

427-001
002-0273-XXX^{2,3}

Description

Basic stool w/o seat cushion
Seat cushion

Accessories

002-0274-XXX^{2,4}
002-0275-XXX^{2,4}

Back kit
Arm kit

Miscellaneous

002-0944-XXX²
002-0945-XXX²
002-0946-XXX²

CAL 133 Seat Cushion-Premium Upholstery
CAL 133 Back Kit-Premium Upholstery
CAL 133 Arm Kit-Premium Upholstery

² - Former color palette available until 10/31/18. See page 19.

³ - Additional charge for Special or Ultra Premium/UltraFree* Upholstery Color

⁴ - Additional charge for Special or Ultra-Premium/UltraFree* Upholstery Color for Accessory

*ULTRALEATHER® by Ultrafabrics®



Model

279-001-XXX^{2,3}

Ritter 279

Physician Seating (Hand Operated)

- Chrome five leg base
- Weight capacity: 225 lb
- Stool weight: 30 lb
- Shipping weight: 35 lb

Description

Basic stool w/hand release - pneumatically adjustable



Model

278-001-XXX^{2,3}

Ritter 278

Physician Seating (Foot Operated)

- Chrome five leg base
- Weight capacity: 225 lb
- Stool weight: 30 lb
- Shipping weight: 35 lb

Description

Basic stool w/foot release - pneumatically adjustable



hand release

Model

277-001-XXX^{2,3}

277-201-XXX^{2,3}

Ritter 277

Physician Seating with Backrest

- Chrome five leg base
- Weight capacity: 225 lb
- Stool weight: 20 lb
- Shipping weight: 25 lb

Description

Basic stool w/back - pneumatically adjustable, hand release

Basic stool w/back - pneumatically adjustable, foot release



foot release



hand release

Model

276-001-XXX^{2,3}

276-201-XXX^{2,3}

Ritter 276

Physician Seating

- Chrome five leg base
- Weight capacity: 350 lb
- Stool weight: 15 lb
- Shipping weight: 20 lb

Description

Basic stool - pneumatically adjustable, hand release

Basic stool - pneumatically adjustable, foot release



foot release

² - Former color palette available until 10/31/18. See page 19.

³ - Additional charge for Special or Ultra Premium/UltraFree* Upholstery Color



Model
275-001-XXX^{2,3}

Ritter 275

Physician Seating with Backrest

- Chrome five leg base
- Weight capacity: 225 lb
- Stool weight: 20 lb
- Shipping weight: 25 lb

Description

Basic stool w/back - manually adjustable



Model
274-001-XXX^{2,3}

Ritter 274

Physician Seating

- Chrome five leg base
- Weight capacity: 350 lb
- Stool weight: 15 lb
- Shipping weight: 20 lb

Description

Basic stool - manually adjustable



Model
273-001-XXX^{2,3}

Ritter 273

Physician Seating with Backrest

- Black composite five leg base
- Weight capacity: 225 lb
- Stool weight: 20 lb
- Shipping weight: 25 lb

Description

Basic stool w/back - pneumatically adjustable



Model
272-001-XXX^{2,3}

Ritter 272

Physician Seating

- Black composite five leg base
- Weight capacity: 300 lb
- Stool weight: 15 lb
- Shipping weight: 20 lb

Description

Basic stool - pneumatically adjustable

² - Former color palette available until 10/31/18. See page 19.

³ - Additional charge for Special or Ultra Premium/UltraFree* Upholstery Color

*ULTRALEATHER® by Ultrafabrics®



Model
271-001-XXX^{2,3}

Ritter 271

Physician Seating with Backrest

- Black composite five leg base
- Weight capacity: 225 lb
- Stool weight: 20 lb
- Shipping weight: 25 lb

Description
Basic stool w/back - manually adjustable



Model
270-001-XXX^{2,3}

Ritter 270

Physician Seating

- Black composite five leg base
- Weight capacity: 300 lb
- Stool weight: 15 lb
- Shipping weight: 20 lb

Description
Basic Stool - manually adjustable

Model	Description
Miscellaneous	
27X-X11-XXX ²	CAL 133 seat upholstery
9A366001-XXX ^{2,4}	Back rest kit (weight capacity, 225 lb)
9A366002-XXX ²	CAL 133 back rest kit (weight capacity, 225 lb)
27X-002-XXX ^{2,3}	Set of soft rubber casters
27X-003-XXX ^{2,3}	Set of glides (not available on foot release models)
27X-004-XXX ^{2,3}	Set of locking casters
27X-005-XXX ^{2,3}	Auto locking casters
9A367001	Foot ring, field installed (not available on foot release models)
9A368001	6" lab height ext. kit (not available on foot release models)

² - Former color palette available until 10/31/18. See page 19.
³ - Additional charge for Special or Ultra Premium/UltraFree* Upholstery Color
⁴ - Additional charge for Special or Ultra-Premium/UltraFree* Upholstery Color for Accessory

*ULTRALEATHER® by Ultrafabrics®



Ritter 280

Patient Side Chair
(Available 11/1/2018)

- 400-pound patient weight capacity
- Wall saver, kickback rear leg design
- Leveling feet
- Chair weight: 25 lb
- Shipping weight: 35 lb

Model

Description

280-001-XXX	Basic chair (without arms)
280-001-XXX	Basic chair (without arms) UltraFree* Upholstery
280-002-XXX	Basic chair (with arms)
280-002-XXX	Basic chair (with arms) UltraFree* Upholstery

Miscellaneous

280-001-999	Basic chair (without arms) Special Color
280-002-999	Basic chair (with arms) Special Color
280-003-XXX	CAL 133 Basic Chair (without arms)
280-004-XXX	CAL 133 Basic Chair (with arms)
280-003-999	CAL 133 Basic Chair (without arms) Special Color
280-004-999	CAL 133 Basic Chair (with arms) Special Color

See color code chart on page 35



Ritter 282

Patient Bariatric Side Chair
(Available 11/1/2018)

- 650-pound patient weight capacity
- Wall saver, kickback rear leg design
- Leveling feet
- Chair weight: 37 lb
- Shipping weight: 47 lb

Model

Description

282-001-XXX	Basic chair
282-001-XXX	Basic chair UltraFree* Upholstery

Miscellaneous

282-001-999	Basic chair Special Color
282-002-XXX	CAL 133 Basic Chair
282-002-999	CAL 133 Basic Chair Special Color

See color code chart on page 35



Midmark 680

Patient Side Chair (Discontinued 12/31/2018)

- 400-pound patient weight capacity
- Wall saver, kickback rear leg design
- Leveling feet
- Chair weight: 20 lb
- Shipping weight: 30 lb

Model

680-001-XXX²

680-001-XXX²

680-002-XXX²

680-002-XXX²

Description

Basic chair (without arms)

Basic chair (without arms) UltraFree* Upholstery

Basic chair (with arms)

Basic chair (with arms) UltraFree* Upholstery

Miscellaneous

680-001-999

Basic chair (without arms) Special Color

680-002-999

Basic chair (with arms) Special Color

680-003-XXX²

CAL 133 Basic Chair (without arms)

680-004-XXX²

CAL 133 Basic Chair (with arms)

680-003-999

CAL 133 Basic Chair (without arms) Special Color

680-004-999

CAL 133 Basic Chair (with arms) Special Color

² - Former color palette available until 10/31/18. See page 19.



Ritter 281

Blood Draw Chair

- 400-pound patient weight capacity
- Optional storage drawer
- Chair weight: 60 lb
- Shipping weight: 75 lb



Model

281-011-XXX²

281-011-XXX²

281-012-XXX²

281-012-XXX²

Description

Basic chair (without storage drawer)

Basic chair (without storage drawer) UltraFree* Upholstery

Basic chair (with storage drawer)

Basic chair (with storage drawer) UltraFree* Upholstery

Miscellaneous

281-011-999

Basic chair (without storage drawer) Special Color

281-012-999

Basic chair (with storage drawer) Special Color

281-013-XXX²

CAL 133 Basic Chair (without storage drawer)

281-014-XXX²

CAL 133 Basic Chair (with storage drawer)

281-013-999


CAL 133 Basic Chair (without storage drawer) Special Color

281-014-999

CAL 133 Basic Chair (with storage drawer) Special Color

² - Former color palette available until 10/31/18. See page 19.

Seating Color-Specific Product Number Chart



Example: 276-001-814
276 Physician Seating, Stone Upholstery

Standard Upholstery Colors

Mist-813	Healing Waters-848
Stone-814	Shaded Garden-853
Restful Path-858	Citrus-849
Iron Ore-851	Sandy Retreat-850
Lunar Gray-845	Curative Copper-846
Obsidian-857	Cranberry-859
Oasis-843	Dark Linen-856
Soothing Blue-855	Latte-860
Dream-847	Robust Brown-852
Mineral-844	Deep Earth-854

UltraFree* Upholstery Colors

UltraFree Stone-865
UltraFree Restful Path-868
UltraFree Obsidian-867
UltraFree Harbor-861
UltraFree Wheatgrass-862
UltraFree Spice-863
UltraFree Cranberry-869
UltraFree Dark Linen-866
UltraFree Latte-870
UltraFree Branch-864

*ULTRALEATHER® by Ultrafabrics®



Mobile Procedure Light



Single Procedure Light



Dual Procedure Light



Wall Mount

Ritter 255

LED Procedure Light

- 4,400 K color temperature
- Color Rendering Index (CRI) greater than 95
- 7,500 fc at 1 meter (39 inches) focal length
- Effortless positioning
- IEC/ES 60601-1, CAN/CSA C22.2 60601-1, 60601-2-41
- Weight: 65 lb (118 lb for dual light)
- Shipping weight: 120 lb (190 lb for dual light)

Model	Description
255-001	Single mount, 8' ceiling
255-002	Single mount, 9' ceiling
255-003	Dual mount, 8' ceiling
255-004	Dual mount, 9' ceiling
255-005	Mobile Procedure Light
255-006	Wall Mount

Accessories

002-1865-00 ^{2, 3}	10' ceiling mount adapter kit
9A584001 ¹	Conversion kit w/cover (single to double light)

¹ - To convert an existing single light to a dual light, this conversion kit must be ordered in addition to a second 8' or 9' light
² - Converts an 8' or 9' light to a 10' ceiling height
³ - One required per light arm

The following model numbers are recommended based on the following ceiling heights (measured from the bottom of the ceiling surface to the floor):

Less than 8'	255-005	Mobile Procedure Light
8' to 9'	255-001, 255-003	8' Single or Dual Procedure Light
9' to 10'	255-002, 255-004	9' Single or Dual Procedure Light
Greater than 10'	255-001, 255-003, 255-005	8' Single or Dual Light with 10' ceiling mount adapter kit(s) or Mobile Procedure Light



Ritter 253

LED Exam Light

- LED technology for up to 50,000 hours of useful life
- Designed to consume 45-60% less energy than other lights
- 4,500 K color temperature
- 1,700 fc intensity (18,298 lux) at 16-inch focal length
- Sleek, ergonomic and easy-to-clean design
- Mist color weight of light unit: 6 lb (253-002 light kit including bracket)
- Shipping weight: 8 lb (253-002 light kit including bracket)
- UL6060-1 and CAN/CSA C22.2 No. 601.1-M90



Exam Table Kit

Mobile Caster Stand Kit

Model	Description
253-011	LED Exam Light - light only - excludes mounting hardware
253-012 ³	LED Exam Light - exam table kit - models 204/22X*/604/62X* - includes light, splines and table mounting hardware, field installed
253-013 ²	LED Exam Light - exam table kit - model 626 for 32" upholstery and Patient Support Rails - includes light, splines and table mounting hardware, field installed
253-015	LED Exam Light - procedure chair kit - models 230/641/646/647 - includes light, base rail and chair mounting hardware, field installed
253-016	LED Exam Light - mobile caster stand kit - includes light & caster stand
253-017	LED Exam Light - wall mount kit - includes light, wall mounting bracket and hardware, field installed

Miscellaneous

9A619001	LED exam table mounting hardware (only) - models 204/22X*/604/62X* - field installed
9A620001	LED exam table mounting hardware (only) - model 626 w/ 32" upholstery - field installed
9A622001	LED procedure chair mounting hardware (only) - models 230/641/646/647 - field installed
9A624001	LED 630 procedure table spline mounting hardware (only), 28" wide upholstery, field installed
9A625001	LED 630 procedure table spline mounting hardware (only), 32" wide upholstery, field installed
9A626001 ¹	LED 630 procedure table rail mount hardware (only), field installed
9A459001	Disposable Sterility Sleeve; 5 boxes of 25 per carton
9A628001	LED light mobile caster stand (only)
9A623001	LED wall mount bracket (only)

* Mounts to tables with a 2006 or newer manufacture date

¹ - Requires (1) accessory rails 9A557001 or 9A557002

² - Use 253-013 with Patient Support Rails™

³ - Slight interference with Patient Support Rails, use 253-013 to minimize



Ritter 250

LED Exam Light

- LED technology for up to 25,000 hours of useful life
- Designed to consume 45-60% less energy than other lights
- 4,000 K color temperature
- 1,800 fc intensity (19,375 lux)
- UL6060-1 and CAN/CSA C22.2 No. 601.1-M90
- Mist color exam light weight: 20 lb
- Shipping weight: 25 lb



Caster Base Accessory

Model	Description
250-003	LED exam light
9A621001	Caster base accessory

Midmark Diagnostic Medical Products (ECG, Holter, Spirometer, Stress, and Vital Signs) Limited Warranty

SCOPE OF WARRANTY

Midmark Corporation ("Midmark") warrants to the original retail purchaser that it will repair or replace devices manufactured by Midmark (except for those not warranted under "Exclusions") that are defective in material or workmanship under normal use and service. Midmark's obligation under this warranty is limited to the repair or replacement, at Midmark's option, of applicable devices. This limited warranty shall only apply to defects that are reported to Midmark within the applicable warranty period and which, upon examination by Midmark, prove to be defective. This warranty extends only to the original retail purchaser of a product, and is not transferable or assignable.

Midmark warrants to the original retail purchaser that it will repair or replace software contained within the products manufactured by Midmark (except for those not warranted under "Exclusions") if: (1) the media on which the software is furnished exhibits defects in material or workmanship under normal use; or (2) the software does not substantially conform to its published specifications. Midmark does not warrant that the software: (1) is error free; (2) can be used without problems or interruptions; or (3) is free from vulnerability to intrusion or attack by viruses or other methods. Midmark also agrees to provide the original retail purchaser with software updates, if any, that become generally available to the public within the warranty period stated, for the specific product owned by the purchaser.

APPLICABLE WARRANTY PERIOD

The applicable warranty period, measured from the date of delivery to the original retail purchaser, shall be as follows:

- (a) One (1) year for software and diagnostic devices (unless otherwise noted below and except items listed under "Exclusions"),
- (b) Two (2) years for vital signs devices,
- (c) One (1) year for reusable blood pressure cuffs and SpO₂ sensors,
- (d) Ninety (90) days for all other accessories (sold in a product kit or purchased separately),

OBTAINING WARRANTY SERVICE

Warranty service must be obtained through Midmark. Contact Midmark Technical Service for warranty service inquiries or issues by phone at 1-800-624-8950, option 2, or via email at techsupport@midmark.com.

It is the original retail purchaser's obligation to arrange for delivery of a product to Midmark for warranty service or exchange. It is also the original retail purchaser's obligation to comply with the warranty service instructions provided by Midmark.

EXCLUSIONS

This warranty does not cover and Midmark shall not be liable for the following:

- (1) any product that has been repaired, modified, or tampered with by any person or entity other than Midmark;
- (2) any product for which any label or tag attached to the product has been removed or tampered with;
- (3) defects, damage, or other conditions caused, in whole or in part, by misuse, abuse, negligence, alteration, accident, freight damage, tampering, or failure to seek and obtain repair or replacement in a timely manner;
- (4) products which are not installed, used, and properly cleaned and maintained as stated in the Midmark Installation, User Guide and/or Operation Manual for the applicable product;
- (5) products considered to be of a consumable nature;
- (6) accessories or parts covered by the respective original manufacturers' warranty, including but not limited to treadmills, computers, printers, scales, carts, mounting options, and carrying cases, after the initial 90 days;
- (7) charges by anyone for adjustments, repairs, replacement parts, installation, or other work performed upon or in connection with such products which are not expressly authorized in writing in advance by Midmark;
- (8) costs and expenses of routine maintenance and cleaning;
- (9) representations and warranties made by any person or entity other than Midmark; and
- (10) items provided through sales promotions.

EXCLUSIVE REMEDY; CONSEQUENTIAL DAMAGES DISCLAIMER

MIDMARK'S ONLY OBLIGATION UNDER THIS WARRANTY IS THE REPAIR OR REPLACEMENT OF DEFECTIVE PARTS. MIDMARK SHALL NOT BE LIABLE FOR AND HEREBY DISCLAIMS ANY DIRECT, SPECIAL, INDIRECT, INCIDENTAL, EXEMPLARY, OR CONSEQUENTIAL DAMAGES OR DELAYS, INCLUDING, BUT NOT LIMITED TO, DAMAGES FOR LOSS OF PROFITS OR INCOME, LOSS OF DATA, LOSS OF USE, DOWNTIME, COVER, AND EMPLOYEE OR INDEPENDENT CONTRACTOR WAGES, PAYMENTS, AND BENEFITS.

NO AUTHORIZATION

No person or firm is authorized to create or approve for Midmark any other obligation or liability in connection with the products.

WARRANTY DISCLAIMER

THIS WARRANTY IS MIDMARK'S ONLY WARRANTY AND IS IN LIEU OF ALL OTHER WARRANTIES, EXPRESS OR IMPLIED. MIDMARK MAKES NO IMPLIED WARRANTIES OF ANY KIND INCLUDING ANY IMPLIED WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. THIS WARRANTY IS LIMITED TO THE REPAIR OR REPLACEMENT OF DEFECTIVE PARTS.

THIS WARRANTY WILL BE VOID IF ANY COVERED PRODUCT IS REPAIRED, MODIFIED, OR TAMPERED WITH BY ANY PERSON OR ENTITY OTHER THAN MIDMARK, OR IF ANY LABEL OR TAG ATTACHED TO A COVERED PRODUCT IS REMOVED OR TEMPERED WITH.

STATUTE OF LIMITATIONS

No action may be brought against Midmark for breach of this limited warranty, an implied warranty, if any, or for any other claim arising out of or relating to the products, more than ninety (90) days following expiration of the limited warranty period. In the event multiple warranty periods exist with respect to a product, the ninety (90) day period provided for herein shall begin to run from expiration of the warranty period for the component to which the claim relates.

Midmark Digital Diagnostic Devices Part Number Ordering Information

Note: Midmark recommends that all Midmark devices drop ship directly to the customer

TO ORDER:

Midmark devices require either an electronic medical record (EMR) that has integrated Midmark software or Midmark IQmanager® Software to operate. Devices include, but are not limited to, the part numbers shown below:

4-000-0025
4-000-0026
4-000-0062
4-000-0110
4-000-0113
4-000-0116

If the customer *currently has* an EMR software package that has integrated Midmark Software, no additional software is required.

Place an order for the appropriate device(s) and please provide:

- EMR name and version
- Customer facility & contact name, address and phone
- Customer computer environment (Thin Client, Fat Client)

If the customer *does not have* an EMR, or using Allscripts TouchWorks™ EHR, Epic® EMR or GE Centricity®, please contact your local Midmark Sales Representative to order the appropriate Midmark software solution. Additional charges or license fees will apply.

For ordering questions regarding Midmark devices, please contact your local Midmark Sales Representative or call 800-624-8950, option 2.

Midmark HEINE Diagnostic Instruments Limited Warranty

SCOPE OF WARRANTY

Midmark Corporation ("Midmark") warrants to the original retail purchaser that it will, at Midmark's option, repair or replace diagnostic instruments manufactured by HEINE and sold by Midmark (except for those not warranted under "Exclusions") that are defective in material or workmanship under normal use and service. The sole remedy under this limited warranty is the repair or replacement, at Midmark's option, of the applicable components. This limited warranty shall only apply to defects that are reported to Midmark within the applicable warranty period and which are determined to exist upon examination by Midmark. This warranty extends only to the original retail purchaser of a product and is not transferable or assignable. Replacement components or products may be used and/or refurbished components or products, provided they are of like quality and specifications as new components or products.

APPLICABLE WARRANTY PERIOD

The applicable warranty period, measured from the date of invoice to the original retail purchaser of the product shall be as follows:

- (1) Ten (10) years for Diagnostic instruments (other than GAMMA G5 Sphygmomanometers) manufactured by HEINE and sold by Midmark;
- (2) Five (5) years for GAMMA G5 Sphygmomanometers manufactured by HEINE and sold by Midmark; and
- (3) Twelve (12) months on parts and labor.

OBTAINING WARRANTY SERVICE

Warranty service must be obtained through Midmark. Midmark may be contacted for warranty service inquiries or issues via email at www.midmark.com; by phone at 1-800-MIDMARK or by mail to Midmark Corporation, 60 Vista Drive, Versailles, Ohio 45380. It is the retail purchaser's obligation to arrange for delivery of a product to Midmark or one of its authorized dealers for warranty service, which delivery shall be at retail purchaser's expense. It is also the retail purchaser's obligation to comply with the warranty service instructions provided either by Midmark or its authorized dealer. The retail purchaser must provide Midmark with completed warranty registration information within thirty (30) days after purchase in order to obtain the benefits of this limited warranty.

EXCLUSIONS:

This limited warranty does not cover and Midmark shall not be liable for the following:

- (1) Defects, damage or other conditions caused, in whole or in part, by misuse, abuse, negligence, alteration, accident, freight damage, negligent storage, tampering or failure to seek and obtain repair or replacement in a timely manner;
- (2) Products which are not installed, used, and properly cleaned and maintained as required or recommended in the Midmark "Installation" and/or "Installation/Operation Manual" for the applicable product, including the specified structural and operational environmental conditions and electrical requirements;
- (3) Normal wear and tear
- (4) Fire or theft
- (5) Products considered to be single-use items or consumables (e.g. bulbs, cuffs, and batteries)
- (6) Rechargeable batteries and batteries
- (7) Use of non-original HEINE spare parts
- (8) Exposed and accessible control and supply cables with corresponding plug connection
- (9) Charges by anyone for adjustments, repairs, replacement parts, installation or other work performed upon or in connection with such products which are not expressly authorized in writing in advance by Midmark;
- (10) Costs and expenses of routine maintenance and cleaning;
- (11) Representations and warranties made by any person or entity other than Midmark.
- (12) Custom manufactured products;
- (13) Alterations or modifications to the product by any person or entity other than Midmark; and
- (14) Products that would otherwise be covered under this limited warranty, but are acquired: (i) from a person or entity that is not Midmark or one of its authorized dealers; or (ii) from a Midmark dealer that is not authorized to sell the product at issue in the geographic territory where the purchaser is located, or is not authorized to sell the product at issue within the medical market

EXCLUSIVE REMEDY; CONSEQUENTIAL DAMAGES DISCLAIMER

MIDMARK'S ONLY OBLIGATION UNDER THIS LIMITED WARRANTY IS THE REPAIR OR REPLACEMENT OF DEFECTIVE PARTS. MIDMARK SHALL NOT BE LIABLE FOR AND HEREBY DISCLAIMS ANY DIRECT, SPECIAL, INDIRECT, INCIDENTAL, EXEMPLARY OR CONSEQUENTIAL DAMAGES OR DELAYS, INCLUDING, BUT NOT LIMITED TO, DAMAGES FOR LOSS OF PROFITS OR INCOME, LOSS OF USE, LOSS OF DATA, DOWNTIME, COVER AND EMPLOYEE OR INDEPENDENT CONTRACTOR WAGES, PAYMENTS AND BENEFITS. THIS DISCLAIMER SHALL SURVIVE ANY FAILURE OR ASSERTED FAILURE OF THE ESSENTIAL PURPOSE OF THIS LIMITED WARRANTY OR ITS REMEDIES SPECIFIED HEREIN.

WARRANTY DISCLAIMER

THIS WARRANTY IS MIDMARK'S ONLY WARRANTY AND IS IN LIEU OF ALL OTHER WARRANTIES, EXPRESS OR IMPLIED. MIDMARK MAKES NO IMPLIED WARRANTIES OF ANY KIND INCLUDING ANY IMPLIED WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. THIS WARRANTY IS LIMITED TO THE REPAIR OR REPLACEMENT OF DEFECTIVE PARTS.

STATUTE OF LIMITATIONS

No action may be brought against Midmark for breach of this limited warranty, an implied warranty, if any, or for any other claim arising out of or relating to the products, more than ninety (90) days following expiration of the limited warranty period.

NO AUTHORIZATION

No person or firm is authorized to create or approve for Midmark any other obligation or liability in connection with the products.

*Warranty information included in product manuals

Midmark Medical Diagnostic Products Support Information

For technical assistance with a Midmark medical diagnostic product, please contact Technical Service at 844-856-1230.

Select: Option 2 – for Midmark Digital Vital Signs & Midmark Digital Vital Signs Monitor, Digital Spirometer, Digital ECG, Digital Holter, Software, Supplies and Midmark Diagnostic Cart. Available Monday-Friday, 8:00 am – 8:00 pm, EST.

Select: Option 3 – for Diagnostic Instrument Wall/Desktop Kits, Instrument Replacement Parts and Supplies (Otoscope/Ophthalmoscope) and Sphygmomanometer Kits. Available 8:00 am – 7:00 pm Monday-Thursday and 8:00 am – 5:00 pm Friday, EST.

For troubleshooting tips and product use information, refer to the operation manual. Additional help topics, self-guided user tools and documentation can be accessed at any time through our support section and knowledge base at midmark.com.

Return Material Authorization

To return any product for repair, an RMA (Return Material Authorization) number must be obtained from a Midmark Technical Service Specialist. The RMA number should be referenced on the package(s) containing the items to be returned.

Shipping

Before shipping any product to Midmark, be certain that an RMA number has been issued and that all guidelines regarding this authorization are followed. We highly recommend that you follow all guidelines for the shipment of medical products set forth by the shipping company you choose to use. If a question should arise regarding the appropriate method of shipment, please feel free to ask when calling for your RMA number. It is the responsibility of the customer when shipping a product to ensure that all packages and their contents get to Midmark safely.

Midmark will not assume responsibility for damage incurred by improper packaging, shipment or product use. Such actions will void all applicable warranties.

Midmark Digital Vital Signs Monitor

with Touchscreen, Serenity BP™ and Exergen TemporalScanner™

- includes Spot, Averaging, Continuous, Interval and SPRINT Protocol BP Modes



Shown with Midmark Digital
Vital Signs Monitor mobile
cart option*

- 1 - Midmark Digital Vital Signs Monitor
- 1 - Reusable Dual Tube BP Cuff – Adult (26-35cm)
- 1 - Reusable Dual Tube BP Cuff – Large Adult (32-42cm)
- 1 - Dual Tube BP Hose – 6.5'
- 1 - Exergen TemporalScanner Thermometer
- 1 - Exergen Disposable Caps (10 ea)
- 1 - Midmark Digital Vital Signs Wireless USB Dongle
- 1 - Midmark Digital Vital Signs USB Cable – 10'
- 1 - Midmark Digital Vital Signs Laptop Tray – 15"
- 1 - Front Storage Bin Kit
- 1 - Device to Temperature Module Cable
- 1 - Three-Arm Knobs (Set of 2)
- 1 - AC Power Supply
- 1 - Setup Guide
- 1 - Quick Reference Guide
- 1 - Operation Manual CD

*Mobile cart not included in pricing.

Model	Description
4-000-0724	Midmark Digital Vital Signs Monitor, with Serenity BP, Exergen TemporalScanner and Masimo® SpO ₂
4-000-0722	Midmark Digital Vital Signs Monitor, with Serenity BP, Exergen TemporalScanner and Nellcor™ SpO ₂
4-000-0720	Midmark Digital Vital Signs Monitor, with Serenity BP and Exergen TemporalScanner

Midmark Digital Non-Wireless Vital Signs Monitor

with Touchscreen, Serenity BP™ and Exergen TemporalScanner™
- includes Spot, Averaging and SPRINT Protocol BP Modes



Shown with Midmark Digital Vital Signs Monitor mobile cart option*

- 1 - Midmark Digital Vital Signs Monitor
- 1 - Reusable Dual Tube BP Cuff – Adult (26-35cm)
- 1 - Reusable Dual Tube BP Cuff – Large Adult (32-42cm)
- 1 - Dual Tube BP Hose – 6.5'
- 1 - Exergen TemporalScanner Thermometer
- 1 - Exergen Disposable Caps (10 ea)
- 1 - Midmark Digital Vital Signs Monitor USB Cable – 10'
- 1 - Midmark Digital Vital Signs Monitor Laptop Tray – 15"
- 1 - Front Storage Bin Kit
- 1 - Device to Temperature Module Cable
- 1 - Three-Arm Knobs (Set of 2)
- 1 - AC Power Supply
- 1 - Setup Guide
- 1 - Quick Reference Guide
- 1 - Operation Manual CD

*Mobile cart not included in pricing.

Model	Description
4-000-0624	Midmark Digital Vital Signs Monitor, non-Wireless, with Serenity BP, Exergen TemporalScanner and Masimo® SpO ₂
4-000-0622	Midmark Digital Vital Signs Monitor, non-Wireless, with Serenity BP, Exergen TemporalScanner and Nellcor™ SpO ₂
4-000-0620	Midmark Digital Vital Signs Monitor, non-Wireless, with Serenity BP and Exergen TemporalScanner

Midmark Digital Vital Signs Monitor Mounting Options

Model	Description
2-200-0072	Midmark Digital Vital Signs Monitor Mobile Cart
2-200-0082	Midmark Digital Vital Signs Monitor Wall Mount Articulating Arm

Midmark Digital Vital Signs Monitor Accessories and Supplies

Model	Description
1-100-1603	Fairbanks® TeleWeigh® Digital Scale
3-200-0030	Blood Pressure Cuff, Dual Tube, Reusable, Child (13-20cm)
3-200-0031	Blood Pressure Cuff, Dual Tube, Reusable, Small Adult (18-26cm)
3-200-0032	Blood Pressure Cuff, Dual Tube, Reusable, Adult (26-35cm)
3-200-0033	Blood Pressure Cuff, Dual Tube, Reusable, Adult Long (29-38cm)
3-200-0034	Blood Pressure Cuff, Dual Tube, Reusable, Large Adult (32-42cm)
3-200-0035	Blood Pressure Cuff, Dual Tube, Reusable, Large Adult Long (35-44cm)
3-200-0036	Blood Pressure Cuff, Dual Tube, Reusable, Thigh (42-50cm)
2-200-0020	Blood Pressure Cuff, Dual Tube, Disposable, Child (13-20cm), 5/pk
2-200-0021	Blood Pressure Cuff, Dual Tube, Disposable, Small Adult (18-26cm), 5/pk
2-200-0022	Blood Pressure Cuff, Dual Tube, Disposable, Adult (26-35cm), 5/pk
2-200-0023	Blood Pressure Cuff, Dual Tube, Disposable, Adult Long (29-38cm), 5/pk
2-200-0024	Blood Pressure Cuff, Dual Tube, Disposable, Large Adult (32-42cm), 5/pk
2-200-0025	Blood Pressure Cuff, Dual Tube, Disposable, Large Adult Long (35-44cm), 5/pk
2-200-0026	Blood Pressure Cuff, Dual Tube, Disposable, Thigh (42-50cm), 5/pk
2-200-0030	SpO ₂ Sensor, Reusable, Masimo®, Pediatric (M-LNCS DCI-P)
2-200-0032	SpO ₂ Sensor, Reusable, Masimo®, Adult (M-LNCS DC-I)
2-200-0031	SpO ₂ Patient Cable, Masimo (RC-4) - 4'
3-200-0007	SpO ₂ Sensor Kit, Reusable, Nellcor™, Pediatric (D-YS and D-YSPD)
2-200-0233	SpO ₂ Sensor, Reusable, Nellcor, Adult (DS100A)
2-200-0234	SpO ₂ Interface Cable, Nellcor (DOC4) - 4'
2-200-0115	Disposable Caps, Exergen Thermometer, 1000pc/box
2-200-0116	Disposable Caps, Exergen Thermometer, 100pc/box
2-200-0074	Storage Bin Kit, Dual, Cart Mounted
2-200-0084	Storage Bin Kit, Wall Mounted
3-200-0009	Tray, Midmark Digital Vital Signs Monitor - 12"
3-200-0012	USB Dongle, Midmark Digital Vital Signs Monitor Wireless

Connectivity Supplies for Midmark 626

Model	Description
9A478001	Midmark Digital Scale serial cable for Midmark Digital Vital Signs, 6' length, coiled
9A478002	Midmark Digital Scale serial cable for Midmark Digital Vital Signs, 15' length, straight
9A478003	Midmark Digital Scale serial cable for Midmark Digital Vital Signs, 30' length, straight
9A478004	Midmark Digital Scale serial cable for Midmark Digital Vital Signs, 50' length, straight
002-10064-00	626 Chair Connectivity Kit



Midmark Digital Vital Signs Device

with Touchscreen - Blood Pressure and Temperature

- 1 - Midmark Digital Vital Signs Device w/Blood Pressure and Temperature
- 1 - AC Power Supply
- 1 - Reusable Blood Pressure Cuff - Adult
- 1 - Reusable Blood Pressure Cuff - Large Adult
- 1 - Blood Pressure Hose - 6.5'
- 1 - Oral/Axillary Temperature Probe
- 1 - Temperature Probe Covers (box of 20)
- 1 - USB Cable - 10'
- 1 - Operation Manual CD
- 1 - Quick Reference Guide

Model

4-000-0500

Description

Midmark Digital Vital Signs with Blood Pressure and Alaris® Turbo Temp® Thermometer

Midmark Digital Vital Signs Device with Touchscreen - Blood Pressure, Temperature and SpO₂

- 1 - Midmark Digital Vital Signs Device w/Blood Pressure, Temperature and SpO₂
- 1 - AC Power Supply
- 1 - Reusable Blood Pressure Cuff - Adult
- 1 - Reusable Blood Pressure Cuff - Large Adult
- 1 - Blood Pressure Hose - 6.5'
- 1 - Reusable SpO₂ Adult Finger Sensor
- 1 - SpO₂ Extender Cable - 4'
- 1 - Oral/Axillary Temperature Probe
- 1 - Temperature Probe Covers (box of 20)
- 1 - USB Cable - 10'
- 1 - Operation Manual CD
- 1 - Quick Reference Guide

Model

4-000-0510

Description

Midmark Digital Vital Signs with Blood Pressure, Alaris Turbo Temp Thermometer and SpO₂

Midmark Digital Vital Signs Accessories and Supplies

Model	Description
1-100-1603	Fairbanks® TeleWeigh® Digital Scale
1-100-1605	Midmark Digital Vital Signs Thermal Printer
3-004-2004	Midmark Digital Vital Signs Mobile Cart Printer Bracket
3-004-2010	Midmark Digital Vital Signs Mobile Cart Scale Mount
3-009-0010	AC Power Supply
3-009-0012	Power Supply Adapter Kit
3-009-0014	Midmark Digital Vital Signs & Midmark Digital Vital Signs PC Replacement lithium ion battery
3-009-0016	USB Cable - 10'
3-009-0018	Serial Cable - 10'
3-009-0022	Blood Pressure Hose - 6.5'
3-009-0100	Blood Pressure Hose - 10'
3-009-0021	Reusable SpO ₂ Small Finger Sensor, Envitec®
3-009-0020	Reusable SpO ₂ Adult Finger Sensor, Envitec
3-009-0026	SpO ₂ Extender Cable, 4', Envitec®
3-009-0024	Alaris® Turbo Temp® Oral/Axillary Temperature Probe
3-009-0058	Turbo Temp Probe Covers, Carton of 10 boxes (20 covers per box)
3-009-0068	Reusable Blood Pressure Cuff, Infant (8-14 cm)
3-009-0070	Reusable Blood Pressure Cuff, Child (13-20 cm)
3-009-0062	Reusable Blood Pressure Cuff, Small Adult (18-26 cm)
3-009-0064	Reusable Blood Pressure Cuff, Adult, Midmark (26-35 cm)
3-009-0066	Reusable Blood Pressure Cuff, Large Adult, Midmark (32-42 cm)
3-009-0072	Reusable Blood Pressure Cuff, Adult Long (29-38 cm)
3-009-0074	Reusable Blood Pressure Cuff, Large Adult Long (35-44 cm)
3-009-0076	Reusable Blood Pressure Cuff, Thigh (42-50 cm)
3-004-2000	Midmark Digital Vital Signs Mobile Cart
3-009-0001	Midmark Digital Vital Signs Countertop Mount
3-009-0003	Midmark Digital Vital Signs Wall Mount with Basket
3-004-2008	Midmark Digital Vital Signs Equipment Pole Mount

Connectivity Supplies for Midmark 626

Model	Description
9A551002	Midmark Digital Vital Signs Table Mount, Non-IQhub
9A478001	Midmark Digital Scale serial cable for Midmark Digital Vital Signs, 6' length, coiled
9A478002	Midmark Digital Scale serial cable for Midmark Digital Vital Signs, 15' length, straight
9A478003	Midmark Digital Scale serial cable for Midmark Digital Vital Signs, 30' length, straight
9A478004	Midmark Digital Scale serial cable for Midmark Digital Vital Signs, 50' length, straight
9A629001	Midmark Digital Vital Signs table mount with USB and serial cables
9A629002	Midmark Digital Vital Signs table mount with serial cable
9A629003	Midmark Digital Vital Signs table mount, non-IQhub®
002-10064-00	626 Chair Connectivity Kit*

*Not recommended for connecting Midmark Digital Vital Signs
to PC running Windows® 10 64-bit OS



Midmark Digital Spirometer

- 1 - Midmark Digital Spirometer
- 10 - Disposable Spirometer Mouthpieces
- 1 - 10/pack Disposable Nose Clips
- 1 - Mouse Pad
- 1 - Operation Manual CD
- 1 - Quick Reference Guide
- 1 - Carrying Case

Model

4-000-0026
4-000-0025

Description

Midmark Digital Spirometer - USB version with syringe
Midmark Digital Spirometer - USB version

Software

Please contact your local Midmark Sales Representative to order the appropriate Midmark software solution, including the new IQmanager® version 10.0.0 for non-EMR users.

Midmark Digital Spirometer Accessories and Supplies

Model

1-100-0007
2-100-1205
2-100-1206
2-100-1207
2-244-0001

Description

3-Liter Spirometer Calibration Syringe
Disposable Spirometer Mouthpieces (box of 25)
Disposable Spirometer Mouthpieces (case - 4 boxes)
Box of 500 Mouthpieces, individually bagged, bulk
Disposable Nose Clips (box of 25)



Midmark Digital ECG with Lead Management (Interpretive)

- 1 - Midmark Digital ECG Acquisition Module
- 1 - 10-Lead Patient Cable, with Lead Management
- 1 - 10/pack Clear ECG Clips
- 1 - 100/pack Disposable EZ-Trode ECG Electrodes
- 1 - Mouse Pad
- 1 - Operation Manual CD
- 1 - Quick Reference Guide
- 1 - Carrying Case

Model

4-000-0062

Description

Midmark Digital ECG - with Lead Management

Software

Please contact your local Midmark Sales Representative to order the appropriate Midmark software solution, including the new IQmanager® version 10.0.0 for non-EMR users.

Digital ECG Supplies

Model

3-047-0001
3-047-0005
3-100-0203
2-100-0205
2-100-0206

Description

*Universal ECG Clips (3mm, 4mm & Snap) 10/pack
Clear ECG Clips (4mm & Snap) 10/pack
Patient Cable, with Lead Management, IQecg
Disposable EZ-Trode ECG Electrodes (box of 500)
Disposable EZ-Trode ECG Electrodes (case - 10 boxes)

*Available to the U.S. market only

Midmark Digital Holter



Holter Kit with Digital Recorder

- 1 - Midmark Digital Holter Recorder
 - Includes: lead wire set, secure digital card
- 1 - Reusable Digital Holter Recorder Pouch and Belt
- 1 - Secure Digital Card Reader
- 1 - Software Security Key, USB Version
- 4 - Holter Prep Kits
- 1 - Mouse Pad
- 1 - Operation Manual CD
- 1 - Quick Reference Guide

Model	Description
4-000-0110	Midmark Digital Holter with one (1) Recorder
4-000-0113	Midmark Digital Holter EX with one (1) Recorder
4-000-0116	Midmark Digital Holter EP with Pacemaker detection and one (1) Recorder

Midmark Digital Holter Recorder

- 1 - Midmark Digital Holter Recorder
- 1 - Lead wire set
- 1 - Secure Digital Card
- 1 - Operation Manual CD
- 4 - Holter Prep Kits
- 1 - Reusable Digital Holter Recorder Pouch and Belt

Model	Description
4-000-0119	Midmark Digital Holter Recorder
1-370-0020	Midmark Digital Holter Recorder Upgrade Kit (24, 48, 72 Hours) One Upgrade Kit Required per Recorder

Software

Please contact your local Midmark Sales Representative to order the appropriate Midmark software solution, including the new IQmanager® version 10.0.0 for non-EMR users.

Secure Memory Card

Model	Description
1-370-0010	Midmark Digital Holter Secure Digital Memory Card
1-370-0004	USB Secure Digital Card Reader

Digital Holter Supplies

Model	Description
2-100-0080	Prep Kit - Midmark Digital Holter, 5-lead (each) Each kit contains: 1 diary, 2 AA batteries, 1 razor 1 Abrading Pad, 2 alcohol wipes, 5 electrodes
2-100-0081	Prep Kit - Midmark Digital Holter, 5-lead (case) Each kit contains: 1 diary, 2 AA batteries, 1 razor 1 Abrading Pad, 2 alcohol wipes, 5 electrodes (case - 30 kits)
2-100-0090	Prep Kit - Midmark Digital Holter, 7-lead (each) Each kit contains: 1 diary, 2 AA batteries, 1 razor 1 Abrading Pad, 2 alcohol wipes, 7 electrodes
2-100-0091	Prep Kit - Midmark Digital Holter, 7-lead (case) Each kit contains: 1 diary, 2 AA batteries, 1 razor 1 Abrading Pad, 2 alcohol wipes, 7 electrodes (case - 30 kits)
90-30-0205	Lead Wire, 5-lead, Midmark Digital Holter (Blue)
90-30-0206	Lead Wire, 7-lead, Midmark Digital Holter (Black)
2-100-0060	Midmark Digital Holter Electrodes (box of 50)
2-100-0050	Disposable Digital Holter Recorder Pouches (pack of 25)
2-100-0051	Reusable Digital Holter Recorder Pouch with Belt



Midmark HEINE Diagnostic Instrument Wall Kits

- 1 - EN200® Wall Transformer
- 1 - Otoscope with polycarbonate, 11.5' (3.5 m) corded handle
- 1 - Ophthalmoscope with polycarbonate, 11.5' (3.5 m) corded handle

Model	Description
A-095-12-208-166	EN200 Wall Transformer with BETA® 400 LED Otoscope and BETA® 200 LED Ophthalmoscope
A-095-12-209-166	EN200 Wall Transformer with BETA 200 LED Otoscope and BETA 200 LED Ophthalmoscope
A-095-12-204-166	EN200 Wall Transformer with BETA® 400 XHL Otoscope and BETA® 200 XHL Ophthalmoscope
A-095-12-205-166	EN200 Wall Transformer with BETA 200 XHL Otoscope and BETA 200 XHL Ophthalmoscope
A-095-12-210-166	EN200 Wall Transformer with K180® XHL Otoscope and K180® XHL Ophthalmoscope

Midmark HEINE Diagnostic Instrument Wall Kit Accessories

Model	Description
X-095-12-006-166	Wallboard (for device combination of EN200 + Tip Dispenser - no instruments)
X-095-12-007-166	Wallboard (for device combination of EN200 + Tip Dispenser + GAMMA® XXL Sphygmomanometer - no instruments)
B-000-11-149-166	AllSpec® Tip Dispenser - white



Midmark HEINE Diagnostic Instrument Desktop Kits

- 1 – NT4 Desktop Charger
- 1 – Otoscope with BETA®4 NT, Lithium-ion rechargeable handle
- 1 – Ophthalmoscope with BETA®4 NT, Lithium-ion rechargeable handle

Model	Description
A-853-24-420-166	NT4 Desktop Charger with BETA® 400 LED Otoscope and BETA® 200 LED Ophthalmoscope
A-832-24-420-166	NT4 Desktop Charger with BETA 200 LED Otoscope and BETA 200 LED Ophthalmoscope
A-853-23-420-166	NT4 Desktop Charger with BETA® 400 XHL Otoscope and BETA® 200 XHL Ophthalmoscope
A-832-23-420-166	NT4 Desktop Charger with BETA 200 XHL Otoscope and BETA 200 XHL Ophthalmoscope
A-879-23-420-166	NT4 Desktop Charger with K180® XHL Otoscope and K180 XHL Ophthalmoscope

Midmark HEINE Diagnostic Instrument Desktop Kit Accessories

Model	Description
B-000-11-149-166	AllSpec® Tip Dispenser - white



Midmark HEINE Diagnostic Instrument Replacement Parts

Model	Description
B-008-11-400-166	BETA® 400 LED Otoscope (instrument head only)
B-002-11-400-166	BETA® 400 XHL Otoscope (instrument head only)
B-008-11-500-166	BETA® 200 LED Otoscope (instrument head only)
B-002-11-500-166	BETA® 200 XHL Otoscope (instrument head only)
B-002-11-550-166	K180® XHL Otoscope (instrument head only)
C-008-30-100-166	BETA 200 LED Ophthalmoscope (instrument head only)
C-002-30-100-166	BETA 200 XHL® Ophthalmoscope (instrument head only)
C-002-30-207-166	K180 XHL Ophthalmoscope (instrument head only)
X-002-99-494-166	NT4 Desktop Charger (without instruments)
X-007-99-396-166	BETA® 4 NT Rechargeable Handle (with battery and handle bottom cap)
X-002-99-394-166	BETA® 4 NT Handle Bottom Cap
X-095-12-250-166	EN200® Wall Transformer (without instruments)



Midmark HEINE Diagnostic Instrument Supplies

Model	Description
B-000-11-127-166	AllSpec® Disposable Tips 4 mm (box of 1,000)
B-000-11-137-166	AllSpec® Disposable Tips 4 mm (case of 10 boxes, 1,000/box)
B-000-11-128-166	AllSpec® Disposable Tips 2.5 mm (box of 1,000)
B-000-11-138-166	AllSpec® Disposable Tips 2.5 mm (case of 10 boxes, 1,000/box)
X-002-88-078-166	XHL® Xenon Halogen 3.5 V Bulb for BETA 400, BETA 200 and K180 Otoscopes (package of 6)
X-002-88-070-166	XHL Xenon Halogen 3.5 V Bulb for BETA 200 Ophthalmoscope (package of 6)
X-002-88-086-166	XHL Xenon Halogen 3.5 V Bulb for K180 Ophthalmoscope (package of 6)
X-007-99-383-166	Rechargeable Lithium-ion Battery



Midmark HEINE GAMMA® XXL LF Sphygmomanometer Kit

- 1 – GAMMA XXL LF Professional Sphygmomanometer
- 1 – Twin-tube adult cuff (29 – 41 cm circ.)
- 1 – Twin-tube adult large cuff (35 – 47 cm circ.)
- 1 - Universal mounting plate

Model	Description
M-000-09-323-166	GAMMA XXL LF Sphygmomanometer with universal mounting plate and two twin-tube cuffs (adult and adult large)

Midmark HEINE GAMMA XXL Accessories

Model	Description
M-000-09-701-166	Twin-tube LF blood pressure cuff, Infant (6 x 28 cm; 8 – 13 cm circumference)
M-000-09-703-166	Twin-tube LF blood pressure cuff, Child (10 x 34 cm; 13 - 20 cm circumference)
M-000-09-705-166	Twin-tube LF blood pressure cuff, Adult Small (14 x 53 cm; 20 – 29 cm circumference)
M-000-09-629-166	Twin-tube LF blood pressure cuff, Adult (14 x 58 cm; 29 – 41 cm circumference)
M-000-09-631-166	Twin-tube LF blood pressure cuff, Adult Large (17 x 70 cm; 35 - 47 cm circumference)
M-000-09-625-166	Twin-tube LF blood pressure cuff, Thigh (20 x 86 cm; >41 cm circumference)



Midmark HEINE GAMMA® G5 Sphygmomanometer Kit

- 1 – GAMMA G5 Standard Sphygmomanometer
- 1 – Single-tube adult cuff (29 – 41 cm circ.)
- 1 – Single-tube adult large cuff (35 – 47 cm circ.)

Model	Description
M-000-09-230-166	GAMMA G5 Sphygmomanometer with two single-tube cuffs (adult and adult large)

Midmark HEINE GAMMA G5 Accessories

Model	Description
M-000-09-700-166	Single-tube LF blood pressure cuff, Infant (6 x 28 cm; 8 – 13 cm circumference)
M-000-09-702-166	Single-tube LF blood pressure cuff, Child (10 x 34 cm; 13 - 20 cm circumference)
M-000-09-704-166	Single-tube LF blood pressure cuff, Adult Small (14 x 53 cm; 20 – 29 cm circumference)
M-000-09-619-166	Single-tube LF blood pressure cuff, Adult (14 x 58 cm; 29 – 41 cm circumference)
M-000-09-706-166	Single-tube LF blood pressure cuff, Adult Large (17 x 70 cm; 35 - 47 cm circumference)
M-000-09-615-166	Single-tube LF blood pressure cuff, Thigh (20 x 86 cm; >41 cm circumference)



Midmark Diagnostic Cart

(Some assembly required)

- 1 - Storage Basket
- 1 - Digital ECG Arm
- 1 - Digital Spirometer Holster
- 1 - Set of Locking Disks

Model	Description
3-004-1000	Midmark Diagnostic Cart with Arm

Midmark Diagnostic Cart Accessories and Supplies

Model	Description
3-004-2008	Midmark Digital Vital Signs Equipment Pole Mount
3-004-6028	Diagnostic Cart Basket (8" x 6" x 8")

Midmark Workstations: Limited Warranty

SCOPE OF WARRANTY

Midmark Corporation warrants to the original purchaser its Midmark Workstation products manufactured by Midmark to be free from defects in material and workmanship under normal interior use and service. Midmark's obligation under this warranty is limited to the repair or replacement, at Midmark's option, of the parts or the products the defects of which are reported to Midmark within the applicable warranty period and which upon examination by Midmark prove to be defective.

APPLICABLE WARRANTY PERIOD

The applicable warranty period, measured from the date of delivery to the original user, shall be as follows:

- A. Five years for all Midmark Workstation products and components (except for the shorter periods in B), including drawer fronts and casters
- B. One year for Power system components and coiled power cords

EXCLUSIONS

This warranty does not cover and Midmark shall not be liable for: (1) Failure to install or use according to the Midmark user and installation manuals (2) Normal wear and tear, repairs and replacements required because of misuse, abuse, negligence, alteration, accident, freight damage or tampering (including, without limitation, use of product in unsuitable environments or conditions); (3) matching of color, grain or texture except to commercially acceptable standards; (4) changes in color caused by natural or artificial light; (5) products which are

not installed, used or properly cleaned as required in Midmark's written Installation Manual; (6) specially manufactured products; (7) products considered to be of a consumable nature such as light bulbs and surge suppression products; (8) accessories or parts not manufactured by Midmark or use of non-authorized Midmark components; (9) charges by anyone (including Midmark's authorized dealers) for adjustments, repairs, replacement parts, installation or other work performed upon or in connection with such products which is not expressly authorized in writing in advance by Midmark. 10) Care Exchange products and accessories (pgs 69-83) follow warranty on page 5.

EXCLUSIVE REMEDY

Midmarks' only obligation under this warranty is the repair or replacement of defective parts. Midmark shall not be liable for any direct, special, indirect, incidental, exemplary or consequential damages or delay.

NO AUTHORIZATION

No person or firm is authorized to create for Midmark any other obligation or liability in connection with the products.

THIS WARRANTY IS MIDMARK'S ONLY WARRANTY AND IS IN LIEU OF ALL OTHER WARRANTIES, EXPRESSED OR IMPLIED. MIDMARK MAKES NO IMPLIED WARRANTIES OF ANY KIND INCLUDING ANY WARRANTIES OF MERCHANTABILITY OR FITNESS FOR ANY PARTICULAR PURPOSE. THIS WARRANTY IS LIMITED TO THE REPAIR OR REPLACEMENT OF DEFECTIVE PARTS.



Midmark 6211 Secure Laptop Workstation

- 26" of height adjustment from 23" to 49"
- Keyed laptop security system
- 23" D x 26" W worksurface
- 11" D x 20.5" W secondary worksurface
- Two total-locking casters in the front lock both the wheel and swivel simultaneously (two non-locking casters in the rear)
- 19.5" D x 20.5" W base dimension
- Weight capacity: 20 lb
- Unit weight: 58 lb
- Shipping weight: 103 lb

Model

6211-001-xxx

Description

Secure Laptop Workstation

¹ - See color code chart on page 85

Accessories

Plant Installed	Field Installed	
BCHL	9A659001	Bar Code Holder-LH
BCHR	9A659001	Bar Code Holder-RH
BPCH	9A658001	Blood Pressure Cuff Holder
CCWT	9A657001	Cart Counterweight
CPUB	9A663001	CPU Holder-Back
CPUF	9A663001	CPU Holder-Front
DBCM	9A653001	Double Basket
DWLW	9A647001	Drawer w/Lock-Work Surf Mnt
DWNW	9A646001	Drawer-Work Surf Mnt
EPWS	9A651001	Equipment Pole
IQEP	3-004-2008	Digital Vital Signs Equipment Pole Mount
MPSF	9A648001	Mobile Printer Shelf
PMSL	9A655001	Power Mount w/Steel Cable Mgmt-LH
PMSR	9A655001	Power Mount w/Steel Cable Mgmt-RH
RBWS	9A650001	Rear Basket
SBCB	9A652001	Single Basket-Back
SBCF	9A652001	Single Basket-Front
SCML	9A668001	Steel Brkt + MG Coiled Cord-LH
SCMR	9A668001	Steel Brkt + MG Coiled Cord-RH
SHEP	3-004-6030	Equipment Pole Mount
SHLH	9A654001	Scale Holder-LH
SHRH	9A654001	Scale Holder-RH
USCB	9A656001	Adjustable Utility Shelf
VIPL	9A649001	VOIP Shelf-LH Mount
VIPR	9A649002	VOIP Shelf-RH Mount



Midmark 6212

Flat Panel PC Workstation

- 26" of height adjustment from 23" to 49"
- 23" D x 26" W worksurface
- Two total-locking casters in the front lock both the wheel and swivel simultaneously (two non-locking casters in the rear)
- Monitor arm w/ 75 mm and 100 mm mounting patterns
- Monitor arm w/ 9.5" independent vertical adj. and up to 40° positive tilt
- Monitor arm capacity range 5 lb min - 20 lb max
- Monitor arm independent rotation of 180°
- 22.5" D x 20.5" W base dimension
- Weight capacity: 30 lb
- Unit weight: 55 lb
- Shipping weight: 100 lb

Model

6212-001-xxx

Description

Flat Panel PC Workstation

¹ - See color code chart on page 85

Accessories

Plant Installed	Field Installed	
AIOM	N/A	Monitor over 10 lb
BCHD	9A664001	Bar Code Holder Mount
BCHL	9A659001	Bar Code Holder-LH
BCHR	9A659001	Bar Code Holder-RH
BPCH	9A658001	Blood Pressure Cuff Holder
CCWT	9A657001	Cart Counterweight
CPUB	9A663001	CPU Holder-Back
CPUF	9A663001	CPU Holder-Front
DBCM	9A653001	Double Basket
DMMT	9A667001	Dual Monitor Mount
DWLW	9A647001	Drawer w/Lock-Work Surf Mnt
DWNW	9A646001	Drawer-Work Surf Mnt
EPWS	9A651001	Equipment Pole
IQEP	3-004-2008	Digital Vital Signs Equipment Pole Mount
LPDH	9A662001	Laptop Display Holder
MPSF	9A648001	Mobile Printer Shelf
PMSL	9A655001	Power Mount w/Steel Cable Mgmnt-LH
PMSR	9A655001	Power Mount w/Steel Cable Mgmnt-RH
RBWS	9A650001	Rear Basket
SBCB	9A652001	Single Basket-Back
SBCF	9A652001	Single Basket-Front
SCML	9A668001	Steel Brkt + MG Coiled Cord-LH
SCMR	9A668001	Steel Brkt + MG Coiled Cord-RH
SHEP	3-004-6030	Spiro Holder EP Mount
SHLH	9A654001	Scale Holder-LH
SHRH	9A654001	Scale Holder-RH
TBDH	9A661001	Tablet Display Holder
USCB	9A656001	Adjustable Utility Shelf
VESA	9A660001	VESA Mounting Plate
VIPL	9A649001	VOIP Shelf-LH Mount
VIPR	9A649002	VOIP Shelf-RH Mount



Midmark 6213

Flat Panel Secure PC Workstation

- 26" of height adjustment from 23" to 49" adj. and up to 40° positive tilt
- Keyed PC security system
- 23" D x 26" W worksurface
- 11" D x 20.5" W secondary worksurface
- Two total-locking casters in the front lock both the wheel and swivel simultaneously (two non-locking casters in the rear)
- Monitor arm w/ 75 mm and 100 mm mounting patterns
- Monitor arm w/ 9.5" independent vertical
- Monitor arm capacity range 5 lb min - 20 lb max
- Monitor arm independent rotation of 180°
- 22.5" D x 20.5" W base dimension
- Weight capacity: 30 lb
- Unit weight: 64 lb
- Shipping weight: 109 lb

Model

6213-001-xxx

¹ - See color code chart on page 85

Description

Flat Panel Secure PC Workstation

Accessories

Plant Installed	Field Installed	
A10M	N/A	Monitor over 10 lb
BCHD	9A664001	Bar Code Holder Mount
BCHL	9A659001	Bar Code Holder-LH
BCHR	9A659001	Bar Code Holder-RH
BPCH	9A658001	Blood Pressure Cuff Holder
CCWT	9A657001	Cart Counterweight
CPUB	9A663001	CPU Holder-Back
CPUF	9A663001	CPU Holder-Front
DBCM	9A653001	Double Basket
DMMT	9A667001	Dual Monitor Mount
DWLW	9A647001	Drawer w/Lock-Work Surf Mnt
DWNW	9A646001	Drawer-Work Surf Mnt
EPWS	9A651001	Equipment Pole
IQEP	3-004-2008	Digital Vital Signs Equipment Pole Mount
LPDH	9A662001	Laptop Display Holder
MPSF	9A648001	Mobile Printer Shelf
PCML	9A669001	Plastic Bin - MG Coiled PWR Cord-LH
PCMR	9A669001	Plastic Bin - MG Coiled PWR Cord-RH
PMBL	9A655002	Power Mount w/Bin Cable Mgmt-LH
PMBR	9A655002	Power Mount w/Bin Cable Mgmt-RH
RBWS	9A650001	Rear Basket
SBCB	9A652001	Single Basket-Back
SBCF	9A652001	Single Basket-Front
SHEP	3-004-6030	Spiro Holder EP Mount
SHLH	9A654001	Scale Holder-LH
SHRH	9A654001	Scale Holder-RH
TBDH	9A661001	Tablet Display Holder
USCB	9A656001	Adjustable Utility Shelf
VESA	9A660001	VESA Mounting Plate
VIPL	9A649001	VOIP Shelf-LH Mount
VIPR	9A649002	VOIP Shelf-RH Mount



Midmark 6214

Procedure Workstation

- 26" of height adjustment from 23" to 49"
- 21" D x 26" W worksurface
- Two total-locking casters in the front lock both the wheel and swivel simultaneously (two non-locking casters in the rear)
- Includes: rear basket, BP cuff holder, integrated spiro holsters and equipment pole
- 19.5" D x 20.5" W base dimension
- Weight capacity: 30 lb
- Unit weight: 56 lb
- Shipping weight: 101 lb

Model

6214-001-xxx

Description

Procedure Workstation

¹ - See color code chart on page 85

Accessories

Plant Installed	Field Installed	
CCWT	9A657001	Cart Counterweight
CPUB	9A663001	CPU Holder-Back
CPUF	9A663001	CPU Holder-Front
DBCM	9A653001	Double Basket
DWLW	9A647001	Drawer w/Lock-Work Surf Mnt
DWNW	9A646001	Drawer-Work Surf Mnt
IQEP	3-004-2008	Digital Vital Signs EP Mount
PMSL	9A655001	Power Mount w/Steel Cable Mgmnt-LH
PMSR	9A655001	Power Mount w/Steel Cable Mgmnt-RH
SBCB	9A652001	Single Basket-Back
SBCF	9A652001	Single Basket-Front
SHEP	3-004-6030	Spiro Holder Equipment Pole Mount
SHLH	9A654001	Scale Holder-LH
SHRH	9A654001	Scale Holder-RH
USCB	9A656001	Adjustable Utility Shelf



Midmark 6215 Compact Workstation

- 26" of height adjustment from 23" to 49"
- 19.5" D x 20.5" W base dimension
- 21" D x 26" W worksurface
- Two total-locking casters in the front lock both the wheel and swivel simultaneously (two non-locking casters in the rear)
- Weight capacity: 30 lb
- Unit weight: 48 lb
- Shipping weight: 93 lb

Model

6215-001-xxx

Description

Compact Workstation

¹ - See color code chart on page 85

Accessories

Plant Installed	Field Installed	
BCHL	9A659001	Bar Code Holder-LH
BCHR	9A659001	Bar Code Holder-RH
BPCH	9A658001	Blood Pressure Cuff Holder
CCWT	9A657001	Cart Counterweight
CPUB	9A663001	CPU Holder-Back
CPUF	9A663001	CPU Holder-Front
DBCM	9A653001	Double Basket
DWLW	9A647001	Drawer w/Lock-Work Surf Mnt
DWNW	9A646001	Drawer-Work Surf Mnt
EPWS	9A651001	Equipment Pole
IQEP	3-004-2008	Digital Vital Signs Equipment Pole Mount
MPSF	9A648001	Mobile Printer Shelf
PMSL	9A655001	Power Mount w/Steel Cable Mgmt-LH
PMSR	9A655001	Power Mount w/Steel Cable Mgmt-RH
RBWS	9A650001	Rear Basket
SBCB	9A652001	Single Basket-Back
SBCF	9A652001	Single Basket-Front
SCML	9A668001	Steel Brkt + MG Coiled Cord-LH
SCMR	9A668001	Steel Brkt + MG Coiled Cord-RH
SHEP	3-004-6030	Spiro Holder EP Mount
SHLH	9A654001	Scale Holder-LH
SHRH	9A654001	Scale Holder-RH
USCB	9A656001	Adjustable Utility Shelf
VIPL	9A649001	VOIP Shelf-LH Mount
VIPR	9A649002	VOIP Shelf-RH Mount



Midmark 6217

Secure Laptop Workstation + Mouse Tray

- 26" of height adjustment from 23" to 49"
- 19.5" D x 20.5" W base dimension
- Keyed laptop security system
- 23" D x 22" W worksurface w/ slideout mouse surface
- 11" D x 20.5" W secondary worksurface
- Two total-locking casters in the front lock both the wheel and swivel simultaneously (two non-locking casters in the rear)
- Weight capacity: 30 lb
- Unit weight: 60 lb
- Shipping weight: 105 lb

Model

6217-001-xxx

Description

Secure Laptop Workstation + Mouse Tray

¹ - See color code chart on page 85

Accessories

Plant Installed	Field Installed	
BCHL	9A659001	Bar Code Holder-LH
BCHR	9A659001	Bar Code Holder-RH
BPCH	9A658001	Blood Pressure Cuff Holder
CCWT	9A657001	Cart Counterweight
CPUB	9A663001	CPU Holder-Back
CPUF	9A663001	CPU Holder-Front
DBCM	9A653001	Double Basket
DWLM	9A647001	Drawer w/Lock-Mouse Trk Mnt
DWNM	9A646001	Drawer-Mouse Track Mnt
EPWS	9A651001	Equipment Pole
IQEP	3-004-2008	Digital Vital Signs Equipment Pole Mount
MPSF	9A648001	Mobile Printer Shelf
RBWS	9A650001	Rear Basket
SBCB	9A652001	Single Basket-Back
SBCF	9A652001	Single Basket-Front
SCML	9A668001	Steel Brkt + MG Coiled Cord-LH
SCMR	9A668001	Steel Brkt + MG Coiled Cord-RH
SHEP	3-004-6030	Spiro Holder Equipment Pole Mount
SHLH	9A654001	Scale Holder-LH
SHRH	9A654001	Scale Holder-RH
USCB	9A656001	Adjustable Utility Shelf
VIPL	9A649001	VOIP Shelf-LH Mount
VIPR	9A649002	VOIP Shelf-RH Mount



Midmark 6218

Flat Panel PC Workstation + Mouse Tray

- 26" of height adjustment from 23" to 49"
- Monitor arm capacity range 5 lb min - 20 lb max
- Monitor arm independent rotation of 180°
- 23" D x 22" W worksurface w/ slideout mouse surface
- 22.5" D x 20.5" W base dimension
- Two total-locking casters in the front lock both the wheel and swivel simultaneously (two non-locking casters in the rear)
- Monitor arm w/ 75 mm and 100 mm mounting patterns
- Monitor arm w/ 9.5" independent vertical adj. and up to 40° positive
- Weight capacity: 30 lb
- Unit weight: 55 lb
- Shipping weight: 100 lb

Model

6218-001-xxx

Description

Flat Panel PC Workstation + Mouse Tray

¹ - See color code chart on page 85

Accessories

Plant Installed	Field Installed	
AIOM	N/A	Monitor over 10 lb
BCHD	9A664001	Bar Code Holder Mount
BCHL	9A659001	Bar Code Holder-LH
BCHR	9A659001	Bar Code Holder-RH
BPCH	9A658001	Blood Pressure Cuff Holder
CCWT	9A657001	Cart Counterweight
CPUB	9A663001	CPU Holder-Back
CPUF	9A663001	CPU Holder-Front
DBCM	9A653001	Double Basket
DMMT	9A667001	Dual Monitor Mount
DWLM	9A647001	Drawer w/Lock-Mouse Trk Mnt
DWNM	9A646001	Drawer-Mouse Track Mnt
EPWS	9A651001	Equipment Pole
IQEP	3-004-2008	Digital Vital Signs Equipment Pole Mount
LPDH	9A662001	Laptop Display Holder
MPSF	9A648001	Mobile Printer Shelf
RBWS	9A650001	Rear Basket
SBCB	9A652001	Single Basket-Back
SBCF	9A652001	Single Basket-Front
SCML	9A668001	Steel Brkt + MG Coiled Cord-LH
SCMR	9A668001	Steel Brkt + MG Coiled Cord-RH
SHEP	3-004-6030	Spiro Holder Equipment Pole Mount
SHLH	9A654001	Scale Holder-LH
SHRH	9A654001	Scale Holder-RH
TBDH	9A661001	Tablet Display Holder
USCB	9A656001	Adjustable Utility Shelf
VESA	9A660001	VESA Mounting Plate
VIPL	9A649001	VOIP Shelf-LH Mount
VIPR	9A649002	VOIP Shelf-RH Mount



Midmark 6219

Flat Panel Secure PC Workstation + Mouse Tray

- 26" of height adjustment from 23" to 49"
- Monitor arm capacity range 5 lb min - 20 lb max
- Keyed PC security system
- Monitor arm independent rotation of 180°
- 23" D x 22" W worksurface w/ slideout mouse surface
- 22.5" D x 20.5" W base dimension
- 11" D x 20.5" W secondary worksurface
- Two total-locking casters in the front lock both the wheel and swivel simultaneously (two non-locking casters in the rear)
- Monitor arm w/ 75 mm and 100 mm mounting patterns
- Monitor arm w/ 9.5" independent vertical adj. & up to 40° positive tilt
- Weight capacity: 30 lb
- Unit weight: 64 lb
- Shipping weight: 109 lb

Model

6219-001-xxx

Description

Flat Panel Secure PC Workstation + Mouse Tray

¹ - See color code chart on page 85

Accessories

Plant Installed	Field Installed	
AIOM	N/A	Monitor over 10 lb
BCHD	9A664001	Bar Code Holder Mount
BCHL	9A659001	Bar Code Holder-LH
BCHR	9A659001	Bar Code Holder-RH
BPCH	9A658001	Blood Pressure Cuff Holder
CCWT	9A657001	Cart Counterweight
CPUB	9A663001	CPU Holder-Back
CPUF	9A663001	CPU Holder-Front
DBCM	9A653001	Double Basket
DMMT	9A667001	Dual Monitor Mount
DWLM	9A647001	Drawer w/Lock-Mouse Trk Mnt
DWNM	9A646001	Drawer-Mouse Track Mnt
EPWS	9A651001	Equipment Pole
IQEP	3-004-2008	Digital Vital Signs Equipment Pole Mount
LPDH	9A662001	Laptop Display Holder
MPSF	9A648001	Mobile Printer Shelf
PCML	9A669001	Plastic Bin - MG Coiled PWR Cord-LH
PCMR	9A669001	Plastic Bin - MG Coiled PWR Cord-RH
RBWS	9A650001	Rear Basket
SBCB	9A652001	Single Basket-Back
SBCF	9A652001	Single Basket-Front
SHEP	3-004-6030	Spiro Holder Equipment Pole Mount
SHLH	9A654001	Scale Holder-LH
SHRH	9A654001	Scale Holder-RH
TBDH	9A661001	Tablet Display Holder
USCB	9A656001	Adjustable Utility Shelf
VESA	9A660001	VESA Mounting Plate
VIPL	9A649001	VOIP Shelf-LH Mount
VIPR	9A649002	VOIP Shelf-RH Mount



Midmark 6281 Workstation

- 22" height adjustment
- 19.5" W x 14" D worksurface
- Slide-out mouse surface
- Up to 45" reach off wall
- Weight capacity: 0-8 lb

- Folds to 9" flat to wall
- Unit weight: 33 lb
- Shipping weight: 40 lb

Model

6281-001

Description

Laptop Wall Mounted Workstation

Accessories

Plant Installed	Field Installed	
n/a	9A588001	20" Extension Arm
n/a	9A591001	Track Mount Bar Code Holder
n/a	9A590001	Tamper Proof PC Holder
n/a	9A589001	Adjustable PC Holder



Midmark 6282 Workstation

- 22" worksurface height adjustment
- 25" W x 10" D worksurface
- Max monitor capacity is 20 lb
- 9.5" independent monitor height adjustment
- Independent 180° monitor rotation and 0-45° tilt

- VESA mount 75 mm and 100 mm
- Folds to 9.47" flat to wall
- Weight capacity: 5 - 20 lb
- Product weight: 43 lb
- Shipping weight: 50 lb

Model

6282-001

Description

Standard Duty Wall Mounted Workstation

Accessories

Plant Installed	Field Installed	
n/a	9A593001	Dual Monitor Arm
n/a	9A588001	20" Extension Arm
n/a	9A592001	VESA Mount Bar Code Holder
n/a	9A591001	Track Mount Bar Code Holder
n/a	9A590001	Tamper Proof PC Holder
n/a	9A589001	Adjustable PC Holder



Midmark 6283 Workstation

- 22" worksurface height adjustment
- 25" W x 10" D worksurface
- Max monitor capacity is 30 lb
- 9.5" independent monitor height adjustment
- Independent 180° monitor rotation and 0-45° tilt
- VESA mount 75 mm and 100 mm
- Folds to 9.47" flat to wall
- Weight capacity: 20-35 lb
- Product weight: 43 lb
- Shipping weight: 50 lb

Model

6283-001

Description

Heavy Duty Wall Mounted Workstation

Accessories

Plant Installed	Field Installed	
n/a	9A593001	Dual Monitor Arm
n/a	9A588001	20" Extension Arm
n/a	9A592001	VESA Mount Bar Code Holder
n/a	9A591001	Track Mount Bar Code Holder
n/a	9A590001	Tamper Proof PC Holder
n/a	9A589001	Adjustable PC Holder



*shown with optional large
Utility Basket

Care Exchange® Midmark 6231

Workstation with Equipment Pole and Rear Basket
(Discontinued 12/31/2018)

- 29" - 47" Height adjustment
- 20" W x 20" D Worksurface
- 20" W x 18" D Base Dimension
- Weight Capacity: 30 lb
- Unit weight: 37 lb
- Shipping weight: 62 lb
- Equipment pole included
- Rear basket included
- Cuff holder included
- Left and Right Spiro holsters included
- Left and Right Cord Disc's included

Model

6231-001-xxx¹

Description

Procedure Workstation

¹ - See color code chart on page 85
- See warranty on page 5

Accessories

Plant Installed	Field Installed	
62XX-501-38	9A489001	Scale Holder
62XX-501-05	9A440001	Multi-Purpose Shelf 20" W
62XX-501-04	9A439001	Utility Basket
62XX-501-12	9A450001	Sharps Container & Bracket
62XX-501-21	9A466001	Large Utility Basket
62XX-501-49	9A596001	Multi-Purpose Shelf Narrow 18" W
n/a	3-004-2008	Digital Vital Signs Equipment Pole Mount



6201-002-xxx shown with optional
coiled power cord w/ (3) outlets

Care Exchange® Midmark 6201 Workstation (Discontinued 12/31/2018)

- 29" to 47" height adjustment
- Keyed laptop security system
- 20" D x 20" W worksurface
- 9" D x 20" W secondary shelf
- 9" D x 8" W mouse surface
- 20" D x 18" W base dimension
- Weight capacity: 30 lb
- Unit weight: 44 lb
- Shipping weight: 69 lb

Model

6201-001-xxx¹
6201-002-xxx¹

¹ - See color code chart on page 85

- See warranty on page 5

Description

Traditional Style Laptop Workstation
Traditional Laptop Workstation w/Cord Tray

Accessories

Plant Installed	Field Installed	
**62XX-501-02	9A437001	Cable Management Holder
62XX-501-04	9A439001	Utility Basket
62XX-501-05	9A440001	Multi-Purpose Shelf 20" W
62XX-501-06	9A441001	Bar Code Holder
62XX-501-07	9A468001	Left Hand Swing-a-Mouse
62XX-501-08	9A443001	Slim Line CPU Holder
62XX-501-09	9A444001	Wide CPU Holder
62XX-501-46	9A595001	Equipment Pole / ECG Holder (Rake Style)
62XX-501-12	9A450001	Sharps Container & Bracket
**62XX-501-13	9A451001	Lockable Drawer
n/a	9A460001	Single Flat Panel Monitor Arm
62XX-501-21	9A466001	Large Utility Basket
**62XX-501-22	9A467001	Rear Basket
62XX-501-23	n/a	Med Drawer w/ Lock
62XX-501-24	n/a	Med Drawer No Lock
62XX-501-25	n/a	Med Drawer Chart Pocket
62XX-501-29	9A480001	Cordless Bar Code Holder
*62XX-501-31	9A482001	Cord Tray Rear Basket
62XX-501-32	9A483001	Thin Client CPU Holder
62XX-501-34	9A485001	Upper Work Surface Extension Brackets
62XX-501-39	9A490001	Cart Handle
*62XX-501-40	9A493001	Coiled Power Cord w/(3) Outlets
62XX-501-41	9A494001	Cart Counterweight
62XX-503-03	n/a	Heavy Duty Casters
n/a	3-004-6030	Spiro Holster - for Equipment Pole
n/a	3-004-2008	Digital Vital Signs Equipment Pole Mount
**62XX-501-45	9A594001	Coiled Power Cord
62XX-501-38	9A489001	Scale Holder
62XX-501-49	9A596001	Multi-Purpose Shelf Narrow 18" W

*Only available on 6201-002-XXX models

**Only available on 6201-001-XXX models



Care Exchange® Midmark 6202 Workstation (Discontinued 12/31/2018)

- 29" to 47" height adjustment
- Slim-line PC holder*
- Monitor arm capacity:
10 lb max -001; 20 lb max -002
- 20" D x 20" W worksurface
- 9" D x 8" W mouse surface
- 20" D x 18" W base dimension
- Weight capacity: 30 lb
- Unit weight: 51 lb
- Shipping weight: 76 lb
- VESA mount 75 mm and 100 mm

Model

6202-001-xxx¹
6202-002-xxx¹

¹ - See color code chart on page 85

- See warranty on page 5

Description

Flat Panel PC Workstation
HD Flat Panel Workstation

Accessories

Plant Installed	Field Installed	
62XX-501-04	9A439001	Utility Basket
62XX-501-05	9A440001	Multi-Purpose Shelf 20" W
62XX-501-06	9A441001	Bar Code Holder
62XX-501-07	9A468001	Left Hand Swing-a-Mouse
62XX-501-08*	9A443001	Slim Line CPU Holder
62XX-501-09	9A444001	Wide CPU Holder
62XX-501-46	9A595001	Equipment Pole / ECG Holder (Rake Style)
62XX-501-12	9A450001	Sharps Container & Bracket
62XX-501-13	9A451001	Lockable Drawer
62XX-501-19	9A452001	Cable Management Holder
62XX-501-22	9A467001	Rear Basket
62XX-501-26	n/a	Med Drawer w/ Lock
62XX-501-27	n/a	Med Drawer No Lock
62XX-501-28	n/a	Med Drawer Chart Pocket
62XX-501-29	9A480001	Cordless Bar Code Holder
62XX-501-30	9A481001	Monitor Arm TC CPU Holder
62XX-501-32	9A483001	Thin Client CPU Holder
62XX-501-39	9A490001	Cart Handle
62XX-501-41	9A494001	Cart Counterweight
62XX-503-03	n/a	Heavy Duty Casters
n/a	3-004-6030	Spiro Holster - for Equipment Pole
n/a	3-004-2008	Digital Vital Signs Equipment Pole Mount
62XX-501-45	9A594001	Coiled Power Cord
62XX-501-38	9A489001	Scale Holder
62XX-501-49	9A596001	Multi-Purpose Shelf Narrow 18" W

*Standard on 6202-001-xxx only



6203-002-xxx shown with optional
coiled power cord w/ (3) outlets

Care Exchange® Midmark 6203 Workstation (Discontinued 12/31/2018)

- 29" to 47" height adjustment
- Keyed laptop security system
- 23" D x 20" W worksurface
- 12" D x 20" W secondary shelf
- 9" D x 8" W mouse surface
- Monitor arm capacity: 10 lb max
- 20" D x 18" W base dimension
- Weight capacity: 30 lb
- Unit weight: 51 lb
- Shipping weight: 76 lb
- VESA mount 75mm and 100mm

Model

6203-001-xxx¹
6203-002-xxx¹

¹ - See color code chart on page 85

- See warranty on page 5

Description

Flat Panel EMR Transition Workstation
Flat Panel EMR Workstation w/Cord Tray

Accessories

Plant Installed	Field Installed	
**62XX-501-02	9A437001	Cable Management Holder
62XX-501-04	9A439001	Utility Basket
62XX-501-05	9A440001	Multi-Purpose Shelf 20" W
62XX-501-06	9A441001	Bar Code Holder
62XX-501-07	9A468001	Left Hand Swing-a-Mouse
62XX-501-08	9A443001	Slim Line CPU Holder
62XX-501-09	9A444001	Wide CPU Holder
62XX-501-46	9A595001	Equipment Pole / ECG Holder (Rake Style)
62XX-501-12	9A450001	Sharps Container & Bracket
**62XX-501-13	9A451001	Lockable Drawer
62XX-501-21	9A466001	Large Utility Basket
62XX-501-23	n/a	Med Drawer w/ Lock
62XX-501-24	n/a	Med Drawer No Lock
62XX-501-25	n/a	Med Drawer Chart Pocket
62XX-501-29	9A480001	Cordless Bar Code Holder
62XX-501-30	9A481001	Monitor Arm TC CPU Holder
*62XX-501-31	9A482001	Cord Tray Rear Basket
62XX-501-32	9A483001	Thin Client CPU Holder
62XX-501-34	9A485001	Upper Work Surface Extension Brackets
62XX-501-39	9A490001	Cart Handle
*62XX-501-40	9A493001	Coiled Power Cord w/(3) Outlets
62XX-501-41	9A494001	Cart Counterweight
62XX-503-03	n/a	Heavy Duty Casters
n/a	3-004-6030	Spiro Holster - for Equipment Pole
n/a	3-004-2008	Digital Vital Signs Equipment Pole Mount
**62XX-501-45	9A594001	Coiled Power Cord
62XX-501-38	9A489001	Scale Holder
62XX-501-49	9A596001	Multi-Purpose Shelf Narrow 18" W

*Only available on 6203-002-XXX models

**Only available on 6203-001-XXX models



Care Exchange® Midmark 6204 Workstation (Discontinued 12/31/2018)

- 29" to 47" height adjustment
- 15" D x 26" W worksurface
- 20" D x 18" W base dimension
- Weight capacity: 10 lb
- Unit weight: 31 lb
- Shipping weight: 56 lb

Model

6204-001-xxx¹

¹ - See color code chart on page 85

- See warranty on page 5

Description

Compact Rectangular Workstation

Accessories

Plant Installed	Field Installed	
62XX-501-04	9A439001	Utility Basket
62XX-501-05	9A440001	Multi-Purpose Shelf 20" W
62XX-501-06	9A441001	Bar Code Holder
62XX-501-08	9A443001	Slim Line CPU Holder
62XX-501-09	9A444001	Wide CPU Holder
62XX-501-46	9A595001	Equipment Pole / ECG Holder (Rake Style)
62XX-501-12	9A450001	Sharps Container & Bracket
62XX-501-19	9A452001	Cable Management Holder
62XX-501-21	9A466001	Large Utility Basket
62XX-501-32	9A483001	Thin Client CPU Holder
62XX-501-41	9A494001	Cart Counterweight
62XX-503-03	n/a	Heavy Duty Casters
n/a	3-004-6030	Spiro Holster - for Equipment Pole
n/a	3-004-2008	Digital Vital Signs Equipment Pole Mount
62XX-501-38	9A489001	Scale Holder
62XX-501-49	9A596001	Multi-Purpose Shelf Narrow 18" W



6205-001-xxx



6205-002-xxx

Care Exchange® Midmark 6205 Workstation (Discontinued 12/31/2018)

- 29" to 47" height adjustment
- 18" D x 28" W worksurface
- 20" D x 18" W base dimension
- Monitor arm capacity:
10lb max on -002
- Weight capacity: 10 lb
- Unit weight: 34 lb
- Shipping weight: 59 lb

Model

6205-001-xxx¹
6205-002-xxx¹

¹ - See color code chart on page 85

- See warranty on page 5

Description

Compact Kidney Workstation
Flat Panel Compact Kidney Workstation

Accessories

Plant Installed	Field Installed	
62XX-501-04	9A439001	Utility Basket
62XX-501-05	9A440001	Multi-Purpose Shelf 20" W
62XX-501-06	9A441001	Bar Code Holder
62XX-501-08	9A443001	Slim Line CPU Holder
62XX-501-09	9A444001	Wide CPU Holder
62XX-501-46	9A595001	Equipment Pole / ECG Holder (Rake Style)
62XX-501-12	9A450001	Sharps Container & Bracket
62XX-501-19	9A452001	Cable Management Holder
62XX-501-21	9A466001	Large Utility Basket
*62XX-501-30	9A481001	Monitor Arm CPU Holder
62XX-501-32	9A483001	Thin Client CPU Holder
62XX-501-41	9A494001	Cart Counterweight
62XX-503-03	n/a	Heavy Duty Casters
n/a	3-004-6030	Spiro Holster - for Equipment Pole
n/a	3-004-2008	Digital Vital Signs Equipment Pole Mount
62XX-501-38	9A489001	Scale Holder
62XX-501-49	9A596001	Multi-Purpose Shelf Narrow 18" W

*Only available on 6205-002-XXX models



Care Exchange® Midmark 6206 Workstation with Security System (Discontinued 12/31/2018)

- 29" to 47" height adjustment
- Keyed laptop security cable system
- 22" D x 32" W worksurface
- 20" D x 18" W base dimension
- Weight capacity: 30 lb
- Unit weight: 44 lb
- Shipping weight: 69 lb

Model

6206-001-xxx¹

¹ - See color code chart on page 85

- See warranty on page 5

Description

Kidney Style Laptop Workstation with Security System

Accessories

Plant Installed	Field Installed	
62XX-501-04	9A439001	Utility Basket
62XX-501-05	9A440001	Multi-Purpose Shelf 20" W
62XX-501-06	9A441001	Bar Code Holder
62XX-501-08	9A443001	Slim Line CPU Holder
62XX-501-09	9A444001	Wide CPU Holder
62XX-501-46	9A595001	Equipment Pole / ECG Holder (Rake Style)
62XX-501-12	9A450001	Sharps Container & Bracket
62XX-501-20	9A453001	Cable Management Holder
62XX-501-21	9A466001	Large Utility Basket
62XX-501-29	9A480001	Cordless Bar Code Holder
62XX-501-32	9A483001	Thin Client CPU Holder
62XX-501-41	9A494001	Cart Counterweight
62XX-503-03	n/a	Heavy Duty Casters
n/a	3-004-6030	Spiro Holster - for Equipment Pole
n/a	3-004-2008	Digital Vital Signs Equipment Pole Mount
62XX-501-45	9A594001	Coiled Power Cord
62XX-501-38	9A489001	Scale Holder
62XX-501-49	9A596001	Multi-Purpose Shelf Narrow 18" W



Care Exchange® Midmark 6207 Workstation (Discontinued 12/31/2018)

- 29" to 47" height adjustment
- 22" D x 32" W worksurface
- 20" D x 18" W base dimension
- Weight capacity: 10 lb
- Unit weight: 37 lb
- Shipping weight: 62 lb

Model

6207-001-xxx¹
¹ - See color code chart on page 85
- See warranty on page 5

Description

Kidney Style Laptop Workstation

Accessories

Plant Installed	Field Installed	
62XX-501-04	9A439001	Utility Basket
62XX-501-05	9A440001	Multi-Purpose Shelf 20" W
62XX-501-08	9A443001	Slim Line CPU Holder
62XX-501-09	9A444001	Wide CPU Holder
62XX-501-12	9A450001	Sharps Container & Bracket
62XX-501-21	9A466001	Large Utility Basket
62XX-501-32	9A483001	Thin Client CPU Holder
62XX-501-41	9A494001	Cart Counterweight
62XX-503-03	n/a	Heavy Duty Casters
62XX-501-38	9A489001	Scale Holder
62XX-501-49	9A596001	Multi-Purpose Shelf Narrow 18" W



Care Exchange® Midmark 6208 Workstation (Discontinued 12/31/2018)

- 29" to 47" height adjustment
- 22" D x 32" W worksurface
- 20" D x 18" W base dimension
- Monitor arm capacity:
10 lb max -001; 20 lb max -002
- Weight capacity: 30 lb
- Unit weight: 51 lb
- Shipping weight: 76 lb

Model

6208-001-xxx¹
6208-002-xxx¹

¹- See color code chart on page 85

- See warranty on page 5

Description

Flat Panel Kidney Style Workstation
HD Flat Panel Kidney Style Workstation

Accessories

Plant Installed	Field Installed	
62XX-501-04	9A439001	Utility Basket
62XX-501-05	9A440001	Multi-Purpose Shelf 20" W
62XX-501-06	9A441001	Bar Code Holder
62XX-501-08	9A443001	Slim Line CPU Holder
62XX-501-09	9A444001	Wide CPU Holder
62XX-501-46	9A595001	Equipment Pole / ECG Holder (Rake Style)
62XX-501-12	9A450001	Sharps Container & Bracket
62XX-501-19	9A452001	Cable Management Holder
62XX-501-21	9A466001	Large Utility Basket
62XX-501-29	9A480001	Cordless Bar Code Holder
62XX-501-30	9A481001	Monitor Arm CPU Holder
62XX-501-32	9A483001	Thin Client CPU Holder
62XX-501-41	9A494001	Cart Counterweight
62XX-503-03	n/a	Heavy Duty Casters
n/a	3-004-6030	Spiro Holster - for Equipment Pole
n/a	3-004-2008	Digital Vital Signs Equipment Pole Mount
62XX-501-45	9A594001	Coiled Power Cord
62XX-501-38	9A489001	Scale Holder
62XX-501-49	9A596001	Multi-Purpose Shelf Narrow 18" W



Care Exchange® Midmark 6261 Workstation (Discontinued 12/31/2018)

- Up to 12 hours of run time between charges
- 29" to 47" height adjustment
- Keyed laptop security system
- 20" D x 20" W worksurface
- 9" D x 20" W secondary shelf
- Coiled power cord
- 9" D x 8" W mouse surface
- 20" D x 18" W base dimension
- Weight capacity: 30 lb
- Unit weight: 90 lb for LiFe option and 127 lb for SLA option
- Shipping weight: 115 lb for LiFe option and 143 for SLA option

Model

6261-001-XXX¹

6261-002-XXX¹

¹ - See color code chart on page 85

- See warranty on page 5

Description

AC Traditional Laptop Cart w/LiFe Battery

AC Traditional Laptop Cart w/SLA Battery

Accessories

Plant Installed	Field Installed	
62XX-501-06	9A441001	Bar Code Holder
62XX-501-07	9A468001	Left Hand Swing-a-Mouse
62XX-501-46	9A595001	Equipment Pole / ECG Holder (Rake Style)
62XX-501-23	n/a	Med Drawer w/ Lock
62XX-501-24	n/a	Med Drawer No Lock
62XX-501-25	n/a	Med Drawer Chart Pocket
62XX-501-29	9A480001	Cordless Bar Code Holder
62XX-501-31	9A482001	Cord Tray Rear Basket
62XX-501-34	9A485001	Upper Work Surface Extension Brackets
n/a	3-004-6030	Spiro Holster - for Equipment Pole
n/a	3-004-2008	Digital Vital Signs Equipment Pole Mount
62XX-501-38	9A489001	Scale Holder



Care Exchange® Midmark 6263 Workstation (Discontinued 12/31/2018)

- Up to 12 hours of run time between charges
- 29" to 47" height adjustment
- Monitor arm: 10 lb max -001 & -002;
20 lb max -003 & -004
- 20" D x 20" W worksurface
- Coiled power cord
- 9" D x 8" W mouse surface
- 20" D x 18" W base dimension
- Weight capacity: 30 lb
- Unit weight: 93 lb for LiFe option
and 130 lb for SLA option
- Shipping weight: 118 lb for
LiFe option and 146 for SLA option
- VESA mount 75 mm and 100 mm

Model

6263-001-XXX¹
6263-002-XXX¹
6263-003-XXX¹
6263-004-XXX¹

Description

AC Flat Panel PC Cart w/LiFe Battery
AC Flat Panel PC Cart w/SLA Battery
AC HD Flat Panel PC Cart w/LiFe Battery
AC HD Flat Panel PC Cart w/SLA Battery

¹ - See color code chart on page 85

- See warranty on page 5

Accessories

Plant Installed	Field Installed	
62XX-501-06	9A441001	Bar Code Holder
62XX-501-07	9A468001	Left Hand Swing-a-Mouse
62XX-501-46	9A595001	Equipment Pole / ECG Holder (Rake Style)
62XX-501-23	n/a	Med Drawer w/ Lock
62XX-501-24	n/a	Med Drawer No Lock
62XX-501-25	n/a	Med Drawer Chart Pocket
62XX-501-29	9A480001	Cordless Bar Code Holder
62XX-501-30	9A481001	Monitor Arm TC CPU Holder
62XX-501-31	9A482001	Cord Tray Rear Basket
62XX-501-34	9A485001	Upper Work Surface Extension Brackets
n/a	3-004-6030	Spiro Holster - for Equipment Pole
n/a	3-004-2008	Digital Vital Signs Equipment Pole Mount
62XX-501-38	9A489001	Scale Holder



Care Exchange® Midmark 6271

Workstation (Discontinued 12/31/2018)

- 29" to 47" height adjustment
- 9" D x 25" W worksurface
- Max monitor weight: 20 lb
- 14" monitor height adjustment
- Worksurface depth adjustable from 14"-18"
- Folds to 5.5" flat to wall
- VESA mount 75 mm and 100 mm
- Product weight: 30 lb
- Shipping weight: 45 lb

Model

6271-001-XXX¹

¹ - See color code chart on page 85
² - Shown with optional shelf
- See warranty on page 5

Description

Wall Mount Adj. Flip-Up Keyboard/
Adj. Flat Screen Monitor

Accessories

Plant Installed	Field Installed	
n/a	9A439001	Utility Basket
62XX-501-06	9A441001	Bar Code Holder
n/a	9A443001	Slim Line CPU Holder
n/a	9A444001	Wide CPU Holder
n/a	9A483001	Thin Client CPU Holder
62XX-501-33	9A484001 ¹	Multi-Purpose Shelf
62XX-501-36	9A487001	CPU Holder



Care Exchange® Midmark 6272 Workstation (Discontinued 12/31/2018)

- Fixed height
- 9" D x 25" W worksurface
- Max monitor weight is 20 lb
- 14" Monitor height adjustment
- Worksurface depth adjustable from 14"-18"
- Folds to 5.5" flat to wall
- VESA mount 75 mm and 100 mm
- Product weight: 15 lb
- Shipping weight: 23 lb

Model

6272-001-XXX¹

¹- See color code chart on page 85

- See warranty on page 5

Description

Wall Mount Adj. Screen, Flip-Up Keyboard
Monitor Mount

Accessories

Plant Installed	Field Installed	
n/a	9A439001	Utility Basket
n/a	9A443001	Slim Line CPU Holder
n/a	9A483001	Thin Client CPU Holder
62XX-501-33	9A484001 ¹	Multi-Purpose Shelf
62XX-501-36	9A487001	CPU Holder



Care Exchange® Midmark 6273 Wall Mount Monitor Bracket (Discontinued 12/31/2018)

- VESA mount 75 mm x 100 mm
- Product weight: 1 lb
- Shipping weight: 2 lb
- Monitor weight capacity is 20 lb

Model

6273-001

- See warranty on page 5

Description

Wall Mount Bracket



Care Exchange® Midmark 6274

Wall Mount Flip-up Keyboard (Discontinued 12/31/2018)

- Worksurface depth adjustable from 14"-18" depth
- Folds to 5.5" from wall
- Product weight: 8 lb
- Shipping weight: 11 lb

Model	Description
6274-001-XXX ¹	Wall Mount Flip-Up Keyboard
¹ - See color code chart on page 85	
- See warranty on page 5	

Accessories

Plant Installed	Field Installed	
n/a	9A439001	Utility Basket
n/a	9A443001	Slim Line CPU Holder
n/a	9A483001	Thin Client CPU Holder
n/a	9A487001	Wall Column CPU Holder



Care Exchange® Midmark 6275

Wall Mount Adj. Flat Screen Monitor Mount (Discontinued 12/31/2018)

- 14" Adj. height
- Rotate and tilt
- VESA mount 75 mm and 100 mm
- Monitor weight capacity: 20 lb
- Product weight: 7 lb
- Shipping weight: 10 lb

Model	Description
6275-001	Wall Mount Adj. Flat Screen Monitor Monitor Mount
- See warranty on page 5	

Accessories

Plant Installed	Field Installed	
n/a	9A439001	Utility Basket
n/a	9A443001	Slim Line CPU Holder
n/a	9A483001	Thin Client CPU Holder
n/a	9A487001	Wall Column CPU Holder



Midmark Care Exchange® 6276-001 Workstation (Discontinued 12/31/2018)

- 29" - 47" height adjustment
- 13" x 25" work surface
- Weight capacity: 20 lb
- Folds to 5.5" to wall
- Product weight: 25 lb
- Shipping weight: 33 lb

Model

6276-001-XXX¹

¹ - See color code chart on page 85

- See warranty on page 5

Description

Wall Mount Adj. Flip-Up Keyboard

6276 Accessories

Plant Installed	Field Installed	
n/a	9A439001	Utility Basket
62XX-501-06	9A441001	Bar Code Holder
n/a	9A443001	Slim Line CPU Holder
n/a	9A444001	Wide CPU Holder
n/a	9A483001	Thin Client CPU Holder

Synthesis® Cabinetry Collection by Midmark: Limited Warranty

SCOPE OF WARRANTY

Midmark Corporation warrants to the original purchaser its Midmark Modular Casework products (Synthesis® Casework Collection) manufactured by Midmark to be free from defects in material and workmanship under normal interior use and service. Midmark's obligation under this warranty is limited to the repair or replacement, at Midmark's option, of the parts or the products the defects of which are reported to Midmark within the applicable warranty period and which upon examination by Midmark prove to be defective.

APPLICABLE WARRANTY PERIOD

The applicable warranty period, measured from the date of delivery to the original user, shall be as follows:

- A. Five years for all Midmark brand casework products and components (except for the shorter periods in B and C), including door and drawer fronts and casters
- B. Three years for task lights
- C. One year for countertop surfaces and accessories

EXCLUSIONS

This warranty does not cover sinks, faucets and plumbing accessories. Any warranties on these items are extended directly by the manufacturer of these items to the original purchaser. Information on these manufacturers' warranties will be enclosed with the applicable products. Also, Midmark will furnish copies of any warranties extended by any such manufacturers upon request.

This warranty does not cover and Midmark shall not be liable for: (1) repairs and replacements required because of misuse, abuse, negligence, alteration, accident, freight damage or tampering; (2) matching of color, grain or texture except to commercially acceptable standards; (3) changes in color caused by natural or artificial light; (4) products which are not installed, used or properly cleaned as required in Midmark's written Installation Manual; (5) specially manufactured products; (6) products considered to be of a consumable nature such as light bulbs and surge suppression products; (7) accessories or parts not manufactured by Midmark; (8) charges by anyone (including Midmark's authorized dealers) for adjustments, repairs, replacement parts, installation or other work performed upon or in connection with such products which is not expressly authorized in writing in advance by Midmark.

EXCLUSIVE REMEDY

Midmark's only obligation under this warranty is the repair or replacement of defective parts. Midmark shall not be liable for any direct, special, indirect, incidental, exemplary or consequential damages or delay.

NO AUTHORIZATION


No person or firm is authorized to create for Midmark any other obligation or liability in connection with the products.

THIS WARRANTY IS MIDMARK'S ONLY WARRANTY AND IS IN LIEU OF ALL OTHER WARRANTIES, EXPRESSED OR IMPLIED. MIDMARK MAKES NO IMPLIED WARRANTIES OF ANY KIND INCLUDING ANY WARRANTIES OF MERCHANTABILITY OR FITNESS FOR ANY PARTICULAR PURPOSE. THIS WARRANTY IS LIMITED TO THE REPAIR OR REPLACEMENT OF DEFECTIVE PARTS.

*Warranty information included in product manuals


Midmark Synthesis® Cabinetry Color-Specific Part Number Explanation and Standard Finish Options

E1RX . 803 . 362



Example: E1 36" H Exam Room Package, Right Hand Sink, No Options, Henna, Concrete Stone

N2A . 803



Example: Dressing Nook, LH Door, Keyed Locks, Henna



Flare Handle
Serenity Panel

Frame / Panel Color:

Solids:

216 Pebble
490 Fawn
649 Whisper
791 Frost

Metallics:

804 Venus Silver (M)
805 Pearl Essence (M)
806 Earthen Bronze (M)
807 Flat Iron (M)

Others:

793 Pebble PVC-Free

Woodgrains:

459 Acorn
487 Earth
794 Dune
795 Flax
796 Nest
797 Hewn
798 Path
799 Sediment
800 Radiance
801 Timber
802 Storm
803 Henna

Laminate Countertop Color:

362 Concrete Stone
443 Kalahari Topaz
364 Graphite Nebula
360 Evening Tigris
453 Milano Quartz
290 Green Tigris
354 Bronze Legacy
274 Pebble Grey
711 Colorado Slate



Exam Room Cabinetry Package E1

Exam Base: 36" H (39" to top of backsplash) x 48" W x 18" D
(2) Doors w/Waste Drop, (2) 4" Drawers and (2) Doors

Exam Overhead: 48" W x 24" H x 13" D w/3 position Glove,
Towel & Cup Dispenser

- 49" W x 19" D Postform laminate top
- Stainless steel drop-in sink
- Gooseneck facet w/ 4" wrist-blade handles
- Waste can
- Serenity panel style w/flare handles
- Shipping Weight: 327 lb

Product Information

Step 1: Configuration

E1

Step 2: Sink Location

L Left Hand

R Right Hand

Step 3: Option Package

X None

A Keyed Lock

Step 4: Cabinetry Finish*

XXX Solid Color Cabinet Finish
(Excludes PVC-Free)

XXX Woodgrain, Metallic & PVC-Free
Cabinet Finish

Step 5: Countertop Color*

XXX Countertop Color

Configuration Example: E1LX.803.362

*See page 85 for Cabinetry Finish and Countertop Color



Exam Room Cabinetry Package E2

Exam Base: 36" H (39" to top of backsplash) x 48" W x 18" D
(2) Doors w/Waste Drop, (3) 4" Drawers and (1) 8" Drawer

Exam Overhead: 48" W x 24" H x 13" D w/3 position Glove,
Towel & Cup Dispenser

- 49" W x 19" D Postform laminate top
- Stainless steel drop-in sink
- Gooseneck facet w/ 4" wrist-blade handles
- Waste can
- Serenity panel style w/flare handles
- Shipping Weight: 331 lb

Product Information

Step 1: Configuration

E2

Step 2: Sink Location

L Left Hand

R Right Hand

Step 3: Option Package

X None

A Keyed Lock

Step 4: Cabinetry Finish*

XXX Solid Color Cabinet Finish
(Excludes PVC-Free)

XXX Woodgrain, Metallic & PVC-Free
Cabinet Finish

Step 5: Countertop Color*

XXX Countertop Color

Configuration Example: E2LX.803.362

*See page 85 for Cabinetry Finish and Countertop Color



Exam Room Cabinetry Package E3

Overall Exam Base Dimensions: 36" H x 54" W x 18" D w/Support Leg

Exam Base Cabinet: 36" H (39" to top of backsplash) x 24" W x 18" D
(2) Doors w/Waste Drop

Exam Cart: 33" H x 24" W x 18" D (3) 4" Drawers, (1) 8" Drawer
w/Kydex Contour top and 60mm casters

Exam Overhead: 24" H x 48" W x 13" D w/(3) position Glove,
Towel & Cup Dispenser

- 55" W x 19" D Postform laminate top
- Stainless steel drop-in sink
- Gooseneck facet w/ 4" wrist-blade handles
- Waste can
- Serenity panel style w/flare handles
- Shipping Weight: 437 lb

Product Information

Step 1: Configuration

E3

Step 2: Sink Location

L Left Hand

R Right Hand

Step 3: Option Package

X None

A Keyed Lock

Step 4: Cabinetry Finish*

XXX Solid Color Cabinet Finish
(Excludes PVC-Free)

XXX Woodgrain, Metallic & PVC-Free
Cabinet Finish

Step 5: Countertop Color*

XXX Countertop Color

Configuration Example: E3LX.803.362

*See page 85 for Cabinetry Finish and Countertop Color



Exam Room Cabinetry Package E4

Overall Exam Base Dimensions: 36" H (39" to top of backsplash)
x 48" W x 18" D

Exam Base Cabinet: 36" H x 24" W x 18" D (2) 4" Drawers & (2) Doors
Wall Mounted Sink: 24" W X 18" D

Exam Overhead: 24" H & 48" W x 13" D w/(3) position Glove,
Towel & Cup Dispenser

- 49" W x 19" D Postform laminate top
- Stainless steel drop-in wall mounted sink
- Gooseneck facet w/ 4" wrist-blade handles
- Serenity panel style w/flare handles
- Shipping Weight: 307 lb

Product Information

Step 1: Configuration

E4

Step 2: Sink Location

L Left Hand

R Right Hand

Step 3: Option Package

X None

A Keyed Lock

Step 4: Cabinetry Finish*

XXX Solid Color Cabinet Finish
(Excludes PVC-Free)

XXX Woodgrain, Metallic & PVC-Free
Cabinet Finish

Step 5: Countertop Color*

XXX Countertop Color

Configuration Example: E4LX.791.362

*See page 85 for Cabinetry Finish and Countertop Color



Exam Room Cabinetry Package E5

Overall Exam Base Dimensions: 36" H (39" to top of backsplash)
x 54" W x 18" D w/Support Leg

Exam Cart: 33" H x 24" W x 18" D Cart (3) 4" Drawers, (1) 8" Drawer
w/Kydex Contour Top and 60mm Casters
Wall Mounted Sink: 24" W X 18" D

Exam Overhead: 24" H & 48" W x 13" D w/(3) position Glove,
Towel & Cup Dispenser

- 55" W x 19" D Postform laminate top
- Stainless steel drop-in wall mounted sink
- Gooseneck facet w/ 4" wrist-blade handles
- Serenity panel style w/flare handles
- Shipping Weight: 401 lb

Product Information

Step 1: Configuration

E5

Step 2: Sink Location

L Left Hand

R Right Hand

Step 3: Option Package

X None

A Keyed Lock

Step 4: Cabinetry Finish*

XXX Solid Color Cabinet Finish
(Excludes PVC-Free)

XXX Woodgrain, Metallic & PVC-Free
Cabinet Finish

Step 5: Countertop Color*

XXX Countertop Color

Configuration Example: E5LX.803.362

*See page 85 for Cabinetry Finish and Countertop Color



Exam Room ADA Cabinetry Package E6

Overall Exam Base Dimensions: 34" H ADA (37" to top of backsplash)
x 48" W x 18" D

Exam Base Cabinet: 34" H ADA x 18" W x 18" D (2) 4" Drawers
& (1) 8" Drawer

ADA Wall Mounted Sink: 30" W X 18" D

Exam Overhead: 24" H x 48" W x 13" D w/(3) position Glove,
Towel & Cup Dispenser

- 49" W x 19" D Postform laminate top
- Stainless steel drop-in wall mounted sink
- Gooseneck facet w/ 4" wrist-blade handles
- Serenity panel style w/flare handles
- Shipping Weight: 294 lb

*Customer is responsible for ensuring product meets Federal, state and local code.

Product Information

Step 1: Configuration

E6

Step 2: Sink Location

L Left Hand

R Right Hand

Step 3: Option Package

X None

A Keyed Lock

Step 4: Cabinetry Finish*

XXX Solid Color Cabinet Finish
(Excludes PVC-Free)

XXX Woodgrain, Metallic & PVC-Free
Cabinet Finish

Step 5: Countertop Color*

XXX Countertop Color

Configuration Example: E6LX.803.362

*See page 85 for Cabinetry Finish and Countertop Color



Exam Room ADA Cabinetry Package E7

Overall Exam Base Dimensions: 34" H ADA (37" to top of backsplash)
x 60" W x 18" D w/Support Leg

Exam Cart: 31" H x 24" W x 18" D Cart (2) 4" Drawers, (1) 8" Drawer
w/Kydex® Contour Top and 60 mm Casters
Wall Mounted Sink: 30" W X 18" D

Exam Overhead: 24" H x 48" W x 13" D w/(3) position Glove,
Towel & Cup Dispenser

- 61" W x 19" D Postform laminate top
- Stainless steel drop-in wall mounted sink
- Gooseneck facet w/ 4" wrist-blade handles
- Serenity panel style w/flare handles
- Shipping Weight: 294 lb

*Customer is responsible for ensuring product meets Federal, state and local code.

Product Information

Step 1: Configuration

E7

Step 2: Sink Location

L Left Hand

R Right Hand

Step 3: Option Package

X None

A Keyed Lock

Step 4: Cabinetry Finish*

XXX Solid Color Cabinet Finish
(Excludes PVC-Free)

XXX Woodgrain, Metallic & PVC-Free
Cabinet Finish

Step 5: Countertop Color*

XXX Countertop Color

Configuration Example: E7LX.803.362

*See page 85 for Cabinetry Finish and Countertop Color



* Shown with Kydex®
contour top

Mobile Treatment Cabinet M2

33" H x 27" W x 18" D (2) 4" Drawers & (2) Doors

- Serenity panel style w/flare handles
- Shipping Weight: 155 lb

Product Information

Step 1: Configuration

M2

Step 2: Depth

1 18" Deep

Step 3: Option Package

A 60mm Casters

B 4" Quiet Casters

C 60mm Casters & Keyed Locks

D 4" Quiet Casters & Keyed Locks

Step 4: Cabinetry Finish*

XXX Solid Color Cabinet Finish (Excludes PVC-Free)

XXX Woodgrain, Metallic & PVC-Free Cabinet Finish

Step 5: Countertop Type

KCN Kydex Contour Top

SEB Soft Edge Bumper Top

SEH Soft Edge Handle Top (29"W)

Configuration Example: M21A.791.KCN

*See page 85 for Cabinetry Finish



* Shown with soft edge bumper top

Mobile Treatment Cabinet M4

- 33" H x 27" W x 18" D (3) 4" Drawers & (1) 8" Drawer
- Serenity panel style w/flare handles
 - Shipping Weight: 152 lb

Product Information
Step 1: Configuration
M4
Step 2. Depth
1 18" Deep
Step 3. Option Package
A 60mm Casters
B 4" Quiet Casters
C 60mm Casters & Keyed Locks
D 4" Quiet Casters & Keyed Locks
Step 4. Cabinetry Finish*
XXX Solid Color Cabinet Finish (Excludes PVC-Free)
XXX Woodgrain, Metallic & PVC-Free Cabinet Finish
Step 5. Countertop Type
KCN Kydex® Contour Top
SEB Soft Edge Bumper Top
SEH Soft Edge Handle Top (29"W)
Configuration Example: M41A.791.SEB

*See page 85 for Cabinetry Finish



Mobile Treatment Cabinet M5

33" H x 27" W x 18" D (5) 4" Drawers

- Serenity panel style w/flare handles
- Shipping Weight: 159 lb

* Shown with soft edge handle top

Product Information

Step 1: Configuration

M5

Step 2. Depth

1 18" Deep

Step 3. Option Package

A 60mm Casters

B 4" Quiet Casters

C 60mm Casters & Keyed Locks

D 4" Quiet Casters & Keyed Locks

Step 4. Cabinetry Finish*

XXX Solid Color Cabinet Finish (Excludes PVC-Free)

XXX Woodgrain, Metallic & PVC-Free Cabinet Finish

Step 5. Countertop Type

KCN Kydex® Contour Top

SEB Soft Edge Bumper Top

SEH Soft Edge Handle Top (29"W)

Configuration Example: M51A.791.SEH

*See page 85 for Cabinetry Finish

Glove, Towel & Cup Dispenser



* D2 shown

- 18" W x 8" D x 24" H
- Shipping weight: 37 lb
- Serenity panel style w/flare handle



* D4 shown

Product Information

Step 1. Configuration

- D1** Glove/Towel/Cup Dispenser - (2) Opening RH Door
- D2** Glove/Towel/Cup Dispenser - (2) Opening LH Door
- D3** Glove/Towel/Cup Dispenser - (3) Opening RH Door
- D4** Glove/Towel/Cup Dispenser - (3) Opening LH Door

Step 2. Option Package

- X** No Lock

Step 3. Cabinetry Finish*

- XXX** Solid Color Cabinet Finish (Excludes PVC-Free)
- XXX** Woodgrain, Metallic & PVC-Free Cabinet Finish

Configuration Example: D2X.803

*See page 85 for Cabinetry Finish



Dressing Nook

- Left-hand door pictured
- 24" W x 7" D x 63" H
- Interior mirror, garment hooks and shoe storage included
- Serenity panel style w/flare handle
- Shipping weight: 138 lb



Product Information

Step 1. Configuration

N1 Dressing Nook - RH Door

N2 Dressing Nook - LH Door

Step 2. Option Package

X No Lock

A Keyed Locks

Step 3. CabinetryFinish*

XXX Solid Color Cabinet Finish (Excludes PVC-Free)

XXX Woodgrain, Metallic & PVC-Free Cabinet Finish

Configuration Example: N2X.803

*See page 85 for Cabinetry Finish



Waste Receptacle

- 12" W x 12" D x 24" H
- Capacity: 33 liquid qts
- Serenity panel style
- Shipping weight: 35 lb

Product Information

Step 1. Configuration

W1 Waste Receptacle

Step 2. Option Package

X No Lock

Step 3. Cabinetry Finish*

XXX Solid Color Cabinet Finish (Excludes PVC-Free)

XXX Woodgrain, Metallic & PVC-Free Cabinet Finish

Configuration Example: W1X.803

*See page 85 for Cabinetry Finish



Instrument Panel

- Instruments not included
- Shipping weight: 17 lb

Product Information

Step 1. Configuration

I1 Instrument Panel - 30" x 20"

I2 Instrument Panel - 34" x 20"

Step 2. Option Package

X No Lock

Step 3. Cabinetry Finish*

XXX Solid Color Cabinet Finish (Excludes PVC-Free)

XXX Woodgrain, Metallic & PVC-Free Cabinet Finish

Configuration Example: I2X.803

*See page 85 for Cabinetry Finish



Ritter Waste Receptacle

- With bag retaining ring
- 14" W x 16" D x 25" H
- Shipping weight: 17 lb
- Capacity: 43 liquid qts.
- Color: Pebble

Model	Description
261-001	Waste receptacle



Ritter M11

Automatic Sterilizer

- Automatic Door
- Chamber size: 11" x 18"
- Meets UL 61010-1, 2nd Edition, IEC 61010-2-040, 1st Edition
- Meets CAN/CSA C22.2 No. 61010-1, 2nd Edition
- Meets CSA C22.2 No. 61010-2-040-07 Part 2-040, 1st Edition
- ASME Certified
- Sterilizer weight: 99 lb (44.9 kg)
- Shipping weight: 131 lb (59.4 kg)

Model

M11-022

Description

Ritter M11 - large table-top automatic sterilizer with automatic door



Ritter M9

Automatic Sterilizer

- Automatic Door
- Chamber size: 9" x 15"
- Meets UL 61010-1, 2nd Edition, IEC 61010-2-040, 1st Edition
- Meets CAN/CSA C22.2 No. 61010-1, 2nd Edition
- Meets CSA C22.2 No. 61010-2-040-07 Part 2-040, 1st Edition
- ASME Certified
- Sterilizer weight: 73 lb (33.1 kg)
- Shipping weight: 81 lb (36.7 kg)

Model

M9-022

Description

Ritter M9 - fully automatic sterilizer with automatic door



Ritter M9D

Manual Door

- Chamber size: 9" x 15"
- Meets UL 61010-1, 2nd Edition, IEC 61010-2-040, 1st Edition
- Meets CAN/CSA C22.2 No. 61010-1, 2nd Edition
- Meets CSA C22.2 No. 61010-2-040-07 Part 2-040 1st Edition
- ASME Certified
- Sterilizer weight: 73 lb (33.1 kg)
- Shipping weight: 81 lb (36.7 kg)

Model

M9D-022

Description

Ritter M9D - automatic sterilizer
with manual door



Ritter M3

Automatic Sterilizer

- Automatic Door
- Chamber size: 12.1" L x 7.6" W x 1.6" D (0.49 gal.)
- Meets UL 61010-1 2nd edition, IEC 61010-2-040
- Meets CAN/CSA C22.2 No. 60110-1 2nd Edition
- Meets CSA C22.2 No. 61010-2-040
- FCC Part 15, Sub-part B
- ASME Certified
- Sterilizer weight: 71 lb (32.2 kg)
- Shipping weight: 80 lb (36.3 kg)

Model

M3-003

Description

Ritter M3 - compact, fast, automatic sterilizer

Sterilizer Accessories

Model	Description
002-0396-00	Speed-Clean (16 oz. single bottle) (C)
002-0396-05	Speed-Clean (12 - 16 oz. bottles)
9A226001	M9/M11 Pouch Rack (D)
9A259001	M9/M11 Impact Dot Matrix Printer (B)
9A307001	M9/M11 Cool Hand Tool (Tray and Cassette Removal Tool) (A)
9A401001	M3 Printer (E)
9A402001	M3 Door Tray (F)
9A404001	M3 Top Cover Protector (G)

(A) Cool Handtool

8.5" L



(B) Impact Dot Matrix Printer

7.5" L x 4.4" W x 2.25" H




(C) SpeedClean

8" H



(D) Pouch Rack



E) M3 Printer



(F) M3 Door Tray



(G) M3 Top Cover Protector





Single VistaCool™

Direct-to-Drain System

- Designed to be compatible with the Ritter M3 sterilizer
- Tank dimensions w/fittings: 7.6" W x 7.6" L x 13.75" H
(19.3 cm x 19.3 cm x 34.9 cm)
- Tank volume: 1.50 gal (5.70 liters)
- Tank weight w/water and fittings: 15.0 lb (6.8 kg)
- Approximate shipping weight: 7.0 lb (3.2 kg)
- Max coolant water pressure: 100 psi (6.9 bar)

Model

9A586001

Description

VistaCool single unit for connecting to one M3 sterilizer



Double VistaCool™

Direct-to-Drain System

- Designed to be compatible with Ritter M11 or M9 sterilizer
- Tank dimensions w/fittings: 4.0" W x 13.0" L x 17.0" H
(10.2 cm x 33.1 cm x 43.2 cm)
- Tank volume: 3.00 gal (11.40 liters)
- Tank weight w/water and fittings: 28.0 lb (12.7 kg)
- Approximate shipping weight: 10.0 lb (4.5 kg)
- Max coolant water pressure: 100 psi (6.9 bar)

Model

9A586002*

Description

VistaCool double unit for connecting to any two current Ritter sterilizers

*Required adapter kits included. Ritter M9 and M11 only compatible with the double unit.

Accessories

Model

Description

- | | |
|----------|---|
| 9A586004 | Air Gap Kit (needed only if required by local code) |
| 9A586005 | Single VistaCool adjustable tank stand |
| 9A586006 | Double VistaCool adjustable tank stand |

VistaCool™ is a trademark of Vista Research Group, Ashland, OH



QuickClean™ Ultrasonic Cleaner QC6/QC6R

- Tank Capacity: 6.6 gal/25L
- Tank Size: 19.9" L x 11.8" D x 7.9" H (50.5 x 30 x 20 cm)
- Overall Size: 21.6" L x 16.5" D x 14.3" H (55 x 42 x 36.5 cm)
- Table Top Unit Weight: 29.6 lb (13.44 kg)
- Table Top Unit Shipping Weight: 32.6 lb (14.78 kg)
- Recessed Unit Weight: 32.9 lb (14.92 kg)
- Recessed Unit Shipping Weight: 35.8 lb (16.26 kg)
- Safety basket INCLUDED
- UL 61010-1:2012
- CAN/CSA C22.2 No. 61010-1:2012

Model	Description
QC6-01	QuickClean Ultrasonic Cleaner – 6.6 gal Tabletop Unit
QC6R-01	QuickClean Ultrasonic Cleaner – 6.6 gal Recessed Unit



QuickClean™ Ultrasonic Cleaner QC3/QC3R

- Tank Capacity: 3.3 gal/12.5L
- Tank Size: 13" L x 11.8" W x 5.9" H (33 x 30 x 15 cm)
- Overall Size: 14.7" L x 16.5" W x 10.4" H (37.5 x 42 x 26.5 cm)
- Table Top Unit Weight: 18.6 lb (8.42 kg)
- Table Top Unit Shipping Weight: 20.9 lb (9.50 kg)
- Recessed Unit Weight: 22.3 lb (10.10 kg)
- Recessed Unit Shipping Weight: 24.6 lb (11.18 kg)
- Safety basket INCLUDED
- UL 61010-1:2012
- CAN/CSA C22.2 No. 61010-1:2012

Model	Description
QC3-01	QuickClean Ultrasonic Cleaner – 3.3 gal Tabletop Unit
QC3R-01	QuickClean Ultrasonic Cleaner – 3.3 gal Recessed Unit



QuickClean™

Ultrasonic Cleaner QC1

- Tank Capacity: 1.2 gal/4.5L
- Tank Size: 11.8" L x 5.9" W x 5.9" H (30 x 15 x 15 cm)
- Overall Size: 13.5" L x 10.2" W x 10.4" H (30 x 26 x 26.5 cm)
- Table Top Unit Weight: 12.7 lb (5.76 kg)
- Table Top Unit Shipping Weight: 14.6 lb (6.64 kg)
- Safety basket INCLUDED
- UL 61010-1:2012
- CAN/CSA C22.2 No. 61010-1:2012

Model	Description
QC1-01	QuickClean Ultrasonic Cleaner – 1.2 gal Tabletop Unit

Ultrasonic Parts and Accessories

9A612001	QC1 QuickClean Beaker Accessory Kit
9A613001	QC3/QC3R QuickClean Beaker Accessory Kit
9A614001	QC6/QC6R QuickClean Beaker Accessory Kit
9A298001	Midmark Enzymatic Cleaner (32 oz)
9A296001	Midmark General Purpose Cleaner (32 oz)
9A297001	Midmark Tartar and Stain Remover (32 oz)
002-10007-00	QC1 Safety Basket (1 basket comes standard with a QC1 unit)
002-10008-00	QC3/QC3R Safety Basket (1 basket comes standard with QC3/QC3R unit)
002-10009-00	QC6/QC6R Safety Basket (1 basket comes standard with QC6/QC6R unit)



Designing better care.™

CARB 93120.2 Phase 2 Compliant

Midmark is an ISO 13485 and ISO 9001 Certified Company.
For more information or a demonstration, contact your Midmark
dealer or call: 1-800-MIDMARK Fax: 1-800-365-8631
or visit our website at midmark.com

© 2008 Midmark Corporation
Midmark Corporation, Dayton, OH
Products subject to improvement changes without notice
Litho in U.S.A. MKT00176

