

Physician Stools

COMPARISON CHART

Features	Ritter® Value Series Stools				Ritter® Classic Series Stools						Midmark® Stools		
	270	271	272	273	274	275	276	277	278	279	425	427	
Seat height	18" - 24.5"	18" - 24.5"	18" - 24.5"	18" - 24.5"	18" - 24.5"	18" - 24.5"	18" - 24.5"	18" - 24.5"	18" - 24.5"	18" - 24.5"	19" - 25.5"	17.5" - 24"	17" - 23"
Seat diameter	16.25"	16.25"	16.25"	16.25"	16.25"	16.25"	16.25"	16.25"	16.25"	16.25"	16.25"	16.875"	16.875"
Seat cushion thickness	3.25"	3.25"	3.25"	3.25"	3.25"	3.25"	3.25"	3.25"	3.25"	3.25"	3.25"	4.16"	4.16"
Weight capacity	300 lb	225 lb	300 lb	225 lb	350 lb	225 lb	350 lb	225 lb	225 lb	225 lb	225 lb	255 lb	255 lb
Composite black 5-leg caster base and hood	•	•	•	•								Composite Grey Base	Composite Grey Base
Polished aluminum/black 5-leg caster base and hood					•	•	•	•	•	•			
Caster base diameter	23"	23"	23"	23"	21.5"	21.5"	21.5"	21.5"	21.5"	21.5"	21.5"	23"	23"
Backrest		•		•		•		•				Optional	Optional
Backrest dimension		14" W x 7.5" H		14" W x 7.5" H		14" W x 7.5" H		14" W x 7.5" H				14" W x 7.5" H	14" W x 7.5" H
Armrests									Surgical	Surgical	Optional	Optional	Optional
Height adjustment mechanism	Screw adjust	Screw adjust	Air lift	Air lift	Screw adjust	Screw adjust	Air lift	Air lift	Air lift	Air lift	Air lift	Air lift hand-ring	Air lift foot-lever
Height control	Manual	Manual	Single lever	Single lever	Manual	Manual	8" Hand ring	8" Hand ring	Foot lever	8" Hand ring	8" Hand ring	10.5" Hand Ring	Foot Lever
Upholstery color options	30	30	30	30	30	30	30	30	30	30	30	30	30
Optional accessories available	Ritter Value Series Stools Locking casters, Auto-locking casters, Soft rubber casters, Glides, 6" lab height extension, Foot ring				Ritter Classic Series Stools Locking casters, Auto-locking casters, Soft rubber casters, Glides (not available on 278), 6" lab height extension, Foot ring						Midmark Series Stools Backrest, Arms, 6" lab height extension, Foot ring		

For product inquiries, please contact the Midmark Medical Customer Experience Department at 1.844.856.1233, option 1. For more information, call 1.800.MIDMARK or visit: midmark.com

MKT-00266 © 2019 Midmark Corporation, Miamisburg, Ohio USA

