

# 1000 Mobile Dental Delivery System



Who doesn't want the best tools for their clinic?

Since oral care is so important for your patients and your practice, we've designed the right tools to help you provide the best care.

The Midmark 1000 Mobile Dental Delivery System uses oil-free technology that eliminates the oil that can contaminate your handpieces. You get a clean, dry aseptic work environment that's good for you, your team, and your patients. Best of all, you don't have to do any more oil changes!

## DESIGNED FOR VETERINARY CLINICS

Staff turnover is a huge issue affecting many veterinary clinics. It's important to get new team members up to speed fast, and intuitive tools make the job easier. Users can switch from handpiece to handpiece without touching a button. The Midmark 1000 automatically adjusts as each handpiece is lifted from its cradle. The Midmark 1000 also features one of the quietest compressors on the market, keeping your work environment more peaceful during procedures.

## ERGONOMICALLY CORRECT

Strain on hands and wrists is a major health concern for veterinarians and techs performing dentistry. When your wrist is already aching after one procedure, you know you're in for a long day.

That's why using the correct handpieces is so important—they can increase your day-to-day stamina and career longevity. The swivel handpiece is designed to work with your hand and wrist movements and sit comfortably in your hand. This design, coupled with the mobility of the Midmark 1000 system, maximizes your work space and movements to give you a wider range of motion while you work on patients.

## CONSTRUCTED FOR DURABILITY

Hair and debris can interfere with some units when you're doing procedure after procedure. You need a unit that can withstand the rigors of a clinical setting. The Midmark 1000 is made from sturdy ABS plastic, protects its internal mechanics with a completely enclosed system, and protects its wheels from debris with medical-grade casters.



## YOU'RE SUPPORTED

Midmark is 100% committed to offering you and your team the training you'll need to become totally proficient on your new equipment. With the free resources at [midmark.com](http://midmark.com), you'll be able to take advantage of some of the clearest, most comprehensive training in the industry. The results will mean better care for patients and a better-looking bottom line.

# MIDMARK 1000 SPECIFICATIONS

<b>Overall Depth</b>	18.5" (46.99 cm)
<b>Overall Width</b>	23" (58.42 cm)
<b>Overall Height</b>	36.5" (92.71 cm)
<b>Casters</b>	3" (7.62 cm) medical grade
<b>Compressor</b>	Oil-free 1/3 hp motor Working pressure: 70-100 psi Flow rate: 2.5 CFM at 40 psi Operational sound level: ≤50 decibels Air reservoir: 1 gallon Manual air tank purge
<b>Water Module</b>	Two 1-liter reservoirs (toggled) Pressurized system: 30-40 psi Aseptic, non-retracting HPC coolant Manual handpiece flush system External quick connects for water and air
<b>Electrical</b>	Voltage: 115 VAC, 6 Amps Frequency: 50/60 Hz Powercord: 10' hospital grade



## MIDMARK 1000 DENTAL DELIVERY SYSTEMS

- 8000-003** Midmark 1000 with non-lighted integrated scaler and fiber optics
- 8000-007S** Midmark 1000 with integrated LED scaler and 2-position fiber optic swivel high-speed handpieces
- 8000-008** Midmark 1000 with integrated LED scaler and fiber optics
- 8000-008S** Midmark 1000 with integrated LED scaler and fiber optic swivel high-speed handpiece

## POPULAR COMPANION PRODUCTS

- 00** Venturi Suction Upgrade Kit
- 01** High- and Low-Speed Swivel Handpieces
- 02** Meisinger Bur Starter Kit
- 03** Optional CROSSTEX Twist Prophy Angle
- 04** Rodent Dental Bur Kit
- 05** Newtron® LED Scaler  
Antimicrobial Cover (not shown)



For a complete listing of available dental delivery supplies and accessories, visit [midmarkserviceparts.com](http://midmarkserviceparts.com).

Midmark is an ISO 13485 and ISO 9001 Certified Company.\*

\*Certain products are not included. See the complete list at: [midmark.com/ISO](http://midmark.com/ISO)

Newtron and Meisinger are trademarks or registered trademarks of their respective companies.

For more information, contact your Midmark dealer or call: 1.800.MIDMARK Outside the USA call: 1.937.526.3662 or visit our website: [midmark.com](http://midmark.com)

© 2019 Midmark Corporation, Miamisburg, Ohio USA  
Products subject to improvement changes without notice. Litho in USA.

007-10160-00 Rev. A1 (7/19)

