

QuickClean™ Ultrasonic Cleaner



Recessed units not shown

POWERFUL CLEANING. CONSISTENT RESULTS.

Providing powerful, effective cleaning and consistent results, the QuickClean Ultrasonic Cleaner is one piece of equipment your team shouldn't be without. The QuickClean's cutting-edge technology can help you create a safer, more efficient work environment by decreasing worker exposure to contaminants and potential sharps injuries while reducing the time and effort needed for cleaning.

- 150% higher efficacy than Soniclean
- Customized cycle times from 1 to 99 minutes
- Full-size safety basket included
- 3-year limited warranty

LEAP

Powerful Technology

The QuickClean unit cleans thoroughly and consistently throughout the entire bath using Frequency-Leap technology. This advancement in cleaning efficiency helps ensure all of your instruments are fully cleaned the first time, every time, no matter where they are placed in the basket.



Simple Operation

The QuickClean is easy to use right out of the box, so you can get your staff up and running with minimal training time. Features like options for preset or customized settings, intelligent software that remembers the last cycle settings, a digital LCD control panel that displays prompts, and a fluid sensor that ensures the bath always has enough water all help simplify your processes and increase productivity.



Range of Sizes and Design Options

The QuickClean comes in three table top sizes (1.2, 3.3 and 6.6 gal) and two recessed options (3.3 and 6.6 gal), so you can choose the option that best fits your configuration and workflow. The stainless steel finish and attractive design coordinate well with any décor.

SPECIFICATIONS

QuickClean Model Numbers:	Table Top Unit: QC1-01	Table Top Unit: QC3-01 Recessed Unit: QC3R-01	Table Top Unit: QC6-01 Recessed Unit: QC6R-01
Capacity:	1.2 gal/4.5 L	3.3 gal/12.5 L	6.6 gal/25 L
Table Top Units Overall Length:	13.5" (34.5 cm)	14.7" (37.5 cm)	21.6" (55 cm)
Overall Width:	10.2" (26 cm)	16.5" (42 cm)	16.5" (42 cm)
Overall Height:	10.4" (26.5 cm)	10.4" (26.5 cm)	14.3" (36.5 cm)
Tank Size:	11.8" L x 5.9" W x 5.9" H (30 x 15 x 15 cm)	13" L x 11.8" W x 5.9" H (33 x 30 x 15 cm)	19.9" L x 11.8" W x 7.9" H (50.5 x 30 x 20 cm)
Recessed unit cut out size:	N/A	15.55" x 14.37" (39.5 x 36.5 cm)	22.48" x 14.37" (57.1 x 36.5 cm)
Table Top Unit Weight:	12.7 lb (5.76 kg)	18.6 lb (8.42 kg)	29.6 lb (13.44 kg)
Table Top Shipping Weight:	14.6 lb (6.64 kg)	20.9 lb (9.50 kg)	32.6 lb (14.78 kg)
Recessed Unit Weight:	N/A	22.3 lb (10.10 kg)	32.9 lb (14.92 kg)
Recessed Unit Shipping Weight:	N/A	24.6 lb (11.18 kg)	35.8 lb (16.26 kg)
Operating Frequency:	32-38 kHz	32-38 kHz	32-38 kHz
Temperature up to:	140°F / 60°C	140°F / 60°C	140°F / 60°C
Stainless Steel Safety Basket:	Included	Included	Included
Electrical Requirements	Input: 115 V +/- 10%, 60 Hz, 110 W, 1.0 A Output: 100 W	Input: 115 V +/- 10%, 60 Hz, 210 W, 1.8 A Output: 200 W	Input: 115 V +/- 10%, 60 Hz, 410 W, 3.6 A Output: 400 W
Standards	UL 61010-1:2012; CAN/CSA C22.2 No. 61010-1:2012; EN 60601-1-2:2007; EN 61326:2013; FCC CFR 47 Part 15 part B	UL 61010-1:2012; CAN/CSA C22.2 No. 61010-1:2012; EN 60601-1-2:2007; EN 61326:2013; FCC CFR 47 Part 15 part B	UL 61010-1:2012; CAN/CSA C22.2 No. 61010-1:2012; EN 60601-1-2:2007; EN 61326:2013; FCC CFR 47 Part 15 part B
Limited Warranty	3 Years	3 Years	3 Years

ACCESSORIES AND SERVICE PARTS



Beaker Accessories

- 9A612001 - QC1
Two-Beaker Accessory
- 9A613001 - QC3/QC3R
Four-Beaker Accessory
- 9A614001 - QC6/QC6R
Six-Beaker Accessory



Extra Safety Baskets

- 002-10007-00 - QC1 Safety Basket
- 002-10008-00 - QC3/QC3R Safety Basket
- 002-10009-00 - QC6/QC6R Safety Basket



Cleaning Solutions

- 9A296001 Midmark General Purpose
Cleaner (32 oz)
- 9A297001 Midmark Tartar and Stain
Remover (32 oz)
- 9A298001 Midmark Enzymatic
Cleaner (32 oz)