



Greater Versatility Proper Ergonomics Increased Efficiency

Consider how often you reach for tools and technologies at the point of care.

Designed with caregivers in mind, Midmark workstations support proper ergonomics and integrate with many veterinary care technologies. You'll have what you need within easy reach.



VERSATILITY SUPPORTS PATIENT-CLIENT INTERACTIONS

Better face-to-face interactions between caregivers, patients and clients occur when equipment barriers are removed. Adjustable workstations allow you to position your worksurface to maintain eye contact and share content, without impacting your focus on the patient.

PROPER ERGONOMICS REDUCE CAREGIVER STRAIN

The repetitive motions of lifting, reaching and twisting during patient care take a toll on many veterinarians and their care teams. Midmark workstations are designed to support ergonomics for greater staff safety and comfort.

POINT-OF-CARE ACCESS ENHANCES EFFICIENCY

If your care areas are less than spacious, mobile and wall-mounted workstations afford the same clinical efficiencies and maneuverability, yet easily stow when not in use. From exam rooms and treatment areas to surgical suites, you'll have what you need, where you need it.

1. HEIGHT ADJUSTMENT

26" of height adjustment allow you to work from an ergonomically-correct position whether seated, standing, stationary or mobile.

2. STURDY CASTER BASE

Effortlessly move the workstation from location to location with the sturdy caster base design.

3. EASY POSITIONING

One-touch lever and intuitive guidance label support easy height adjustment for each user.

MIDMARK WORKSTATION SPECIFICATIONS



6212	6215	6282
<p>Larger workstation supports flat-panel PC use. Perfect for recovery areas, charting in treatment areas and easily moving between care zones.</p> <p>Worksurface Dimensions: 23" D x 26" W</p>	<p>Smaller, compact workstation is ideal for users with wireless devices and laptops. Three worksurface dimensions available to fit your unique application.</p> <p>Worksurface Dimensions: 6215-001— 21" D x 26" W (Rectangular) 6215-002— 17.5" D x 23.5" W (Compact Rectangular) 6215-003— 17" D x 25" W (Compact Kidney)</p>	<p>Wall-mounted workstation maximizes floor space in exam rooms and recovery areas. Station folds easily for extra space.</p> <p>Worksurface Dimensions: 10" D x 25" W</p>
<p>Base Dimensions: 22.5" D x 20.5" W</p>	<p>Base Dimensions: 19.5" D x 20.5" W</p>	<p>N/A</p>
<ul style="list-style-type: none"> • 3.94" oversized casters • 2 total-locking casters front; 2 non-locking casters rear • Polished aluminum column and caster base 	<ul style="list-style-type: none"> • 180° of left-to-right monitor rotation • 22 worksurface colors available 	<ul style="list-style-type: none"> • AC power supply • UL/IEC 60601-1 (3rd Ed.) UL/IEC 60950-1
<p>Weight Limit: 30 lb</p>	<p>Weight Limit: 10 lb</p>	<p>Weight Limit: 20 lb</p>

MOBILE WORKSTATION ACCESSORIES:

AC Power System



Column Basket



Coiled Power Cord



Drawer



Adjustable Utility Shelf



Equipment Pole



Rear Basket



VOIP Shelf



WALL UNIT ACCESSORIES:

Dual Monitor Arm



Extension Arm



For a complete list of available mobile workstation accessories, scan the code to the left.