



# Midmark Big Top Grooming Table

## Built for Big Jobs and **Bold Ideas**



With an oversized tabletop and thoughtful ergonomic features, the Big Top™ Grooming Table gives you the space and support to work the way you want. It lifts pets securely and adjusts smoothly to match your grooming style. The pivoting grooming arm offers plenty of flexibility to suit different techniques and preferences. Built to last, it's a durable, reliable tool for your busiest days.

- Accommodates larger breeds and complex tasks with an oversized 48" L x 26" W tabletop.
- Provides adjustable support with a grooming arm that pivots 180° around the table end, locks in five horizontal positions and adjusts vertically from 22" to 39".
- Lifts smoothly and quietly with dual-side foot pedal access—electric or hydraulic lift options match workflow and space needs.



Reduces the risk of injury with rounded tabletop corners and a silicone ball at the arm's end.

# MIDMARK BIG TOP GROOMING TABLE SPECIFICATIONS

| Detail              | Description  |
|---------------------|--|
| Material            | Powder-coated heavy-gauge steel frame; polyurethane-coated table surface; stainless steel grooming arm   |
| Dimensions          | Tabletop: 48" L x 26" W  |
| Lift Range (Height) | 20" to 41" (Electric and Hydraulic models)   |
| Arm Adjustments     | <ul style="list-style-type: none"><li>• Rotates 180°</li><li>• Locks in 5 horizontal positions</li><li>• Adjusts vertically from 22" to 39"</li><li>• Loop hook with 6 position choices</li><li>• Optional rear arm available separately</li></ul> |
| Maximum Capacity    | 200 lb   |
| Lifting Mechanism   | Electric or Hydraulic options  |
| Table Weight        | 190 lb (Electric)<br>145 lb (Hydraulic)  |
| Certifications      | MET Certified/UL Listed - Electric lift models only  |

## ORDERING GUIDE

| Part Number | Description                 |
|-------------|-----------------------------|
| 903-3220-53 | Electric Lift, Plum         |
| 903-3220-55 | Electric Lift, Green        |
| 903-3220-57 | Electric Lift, Orange       |
| 903-3220-59 | Electric Lift, Pink         |
| 903-3220-61 | Electric Lift, Silver Vein  |
| 903-3220-63 | Electric Lift, Blue         |
| 903-3220-77 | Electric Lift, Grey         |
| 903-3220-65 | Hydraulic Lift, Plum        |
| 903-3220-66 | Hydraulic Lift, Green       |
| 903-3220-67 | Hydraulic Lift, Orange      |
| 903-3220-68 | Hydraulic Lift, Pink        |
| 903-3220-69 | Hydraulic Lift, Silver Vein |
| 903-3220-70 | Hydraulic Lift, Blue        |
| 903-3220-74 | Hydraulic Lift, Grey        |



## AVAILABLE COLORS



Scan for additional grooming solutions

In the USA call: 1.800.MIDMARK. Outside the USA call: 1.937.526.3662 or visit our website: [midmark.com](http://midmark.com)

Midmark is an ISO 13485 and ISO 9001 Certified Company. Certain products are not included. See the complete list at: [midmark.com/ISO](http://midmark.com/ISO)

© 2025 Midmark Corporation, Versailles, Ohio USA  
Products subject to improvement changes without notice.

007-10534-00 Rev A1 (7/25)

