



Midmark Elite Grooming Table

Designed
for Comfort.
Built to Last.



The Elite Grooming Table brings together everyday comfort and long-lasting performance. Its coated surface is stable and secure, while the adjustable grooming arm helps you work efficiently and comfortably. Available in electric or hydraulic models, this table is built with durable materials to perform reliably in high-volume grooming environments.

- Accommodates small to mid-sized breeds with a large 42.4" L x 26" W tabletop.
- Provides adjustable support with a grooming arm that pivots 180° around the table end, locks in five horizontal positions, and adjusts vertically from 22" to 39".
- Lifts smoothly and quietly with dual-side foot pedal access. Electric or hydraulic lift options match workflow and space needs.



Reduces the risk of injury with rounded tabletop corners and a silicone ball at the arm's end.

MIDMARK ELITE GROOMING TABLE SPECIFICATIONS

Detail	Description
Material	Powder-coated heavy-gauge steel frame; polyurethane-coated table surface; stainless steel grooming arm
Dimensions	Tabletop: 42.4" L x 26" W
Lift Range (Height)	20" to 41" (Electric and Hydraulic models)
Arm Adjustments	<ul style="list-style-type: none">• Rotates 180°• Locks in 5 horizontal positions• Adjusts vertically from 22" to 39"• Loop hook with 6 position choices• Optional rear arm available separately
Maximum Capacity	200 lb
Lifting Mechanism	Electric or Hydraulic options
Table Weight	187 lb (Electric) 137 lb (Hydraulic)
Certifications	MET Certified/UL Listed - Electric lift models only

ORDERING GUIDE

Part Number	Description
903-3220-19	Electric Lift, Plum
903-3220-20	Electric Lift, Green
903-3220-21	Electric Lift, Orange
903-3220-22	Electric Lift, Pink
903-3220-24	Electric Lift, Silver Vein
903-3220-42	Electric Lift, Blue
903-3220-75	Electric Lift, Grey
903-3220-37	Hydraulic Lift, Plum
903-3220-38	Hydraulic Lift, Green
903-3220-39	Hydraulic Lift, Orange
903-3220-40	Hydraulic Lift, Pink
903-3220-41	Hydraulic Lift, Silver Vein
903-3220-44	Hydraulic Lift, Blue
903-3220-73	Hydraulic Lift, Grey



AVAILABLE COLORS



Scan for additional grooming solutions

In the USA call: 1.800.MIDMARK. Outside the USA call: 1.937.526.3662 or visit our website: midmark.com

Midmark is an ISO 13485 and ISO 9001 Certified Company. Certain products are not included. See the complete list at: midmark.com/ISO

© 2025 Midmark Corporation, Versailles, Ohio USA
Products subject to improvement changes without notice.

007-10533-00 Rev A1 (7/25)

