

Fold-Down Wall-Mount Exam Table

A Better Way to Use Your Space



Combines smooth operation, an integrated work surface and a clean, finished appearance.

The Midmark Fold-Down Wall-Mount Veterinary Exam Table—our best-selling wall-mounted veterinary exam table—helps maximize usable floor space. The integrated 16-inch stationary shelf keeps charts, instruments and supplies within reach to support workflow and organization.

- **Support Flexible Room Use:**
Fold-away design helps create an open, adaptable exam room.
- **Raise and Lower with Ease:**
Gas-assist lifting mechanism provides controlled raising and lowering.
- **Maintain a Clean Appearance:**
Finished stainless steel tabletop remains visible when folded.
- **Enhance Patient Comfort:**
Compatible traction mats are available to support patients during exams.



MIDMARK FOLD-DOWN WALL-MOUNT EXAM TABLE SPECIFICATIONS



Detail	Description
Overall Table in Raised Position	42" L x 24" W x 1½" D
Exam Surface	26" L x 24" W
Stored Depth (Folded)	16" from wall
Table Weight	113 lb
Patient Capacity	125 lb

ORDERING GUIDE

Part Number	Description
8012-00-001	Fold Down Exam Table
9A549001	Traction Mat Accessory



Scan for additional table solutions

Note: Optional mounting bracket (002-1692-00, sold separately) or additional backing structure may be required for steel stud applications.