

Fold-Up Wall-Mount Exam Table

Big
Function.
Small
Footprint.



Designed to maximize usable space with a slim wall profile.

The Midmark Fold-Up Wall-Mount Veterinary Exam Table helps maximize usable space without sacrificing functionality. It stores just 4.5 inches from the wall when not in use, yet provides a durable stainless steel exam surface when needed. Its streamlined design is ideal for compact exam rooms.

- **Maximize Exam Room:**
Slim wall profile helps keep valuable floor space clear.
- **Raise and Lower Smoothly:**
Counterbalance mechanism provides controlled, quiet operation.
- **360° Patient Access:**
Supports efficient exams and patient handling.
- **Simplify Cleanup:**
Raised marine edge helps contain fluids.
- **Built for Daily Use:**
Durable stainless steel construction stands up to the demands of veterinary practice.




MIDMARK FOLD-UP WALL-MOUNT EXAM TABLE SPECIFICATIONS



Detail	Description
Worksurface Lengths	38" and 44" models
Worksurface Width	22"
From Wall	4.5" folded
Materials	– Stainless steel tabletop with raised anti-drip edges – Wall-mount plate with corrosion-resistant baked-on finish
Table Weight	130 lb
Patient Capacity	150 lb

ORDERING GUIDE

Part Number	Description
101-2640-00	44" Length, Stainless Steel Top
171-2940-00	38" Length, Stainless Steel Top
803-0201-00	Stat-Mat, Blue, 20"W x 44"L x ½"H



Scan for additional table solutions