

Synthesis® Fixed-Height Dry Treatment Tables

Built for the Demands of Daily Care



The Midmark Synthesis Fixed-Height Dry Treatment Table is built to support daily patient care and preparation. Designed for efficiency and cleanability, it integrates seamlessly with Midmark cabinetry to create a professional, cohesive environment. With multiple configurations and storage options, your team can help stay organized and focused on delivering excellent care.

- Choose a peninsula-style wall or chase mount, or freestanding island configuration to fit your space and workflow.
- Accommodate a variety of procedures with a durable stainless steel or laminate work surface.
- Maintain a quiet, low-stress environment with soft-close hinges and slides.
- Simplify upkeep with thermofoil cabinet panels that resist moisture and chemical damage.
- Support infection control protocols with optional antimicrobial pulls and surfaces.*




Peninsula-style Configuration



Island Configuration

SYNTHESIS® FIXED-HEIGHT DRY TREATMENT TABLE SPECIFICATIONS



TECHNICAL SPECIFICATIONS

Height	36"
Length	48" or 60" (includes 18" overhang)
Width	26"
Weight	250 lb (48" table) 275 lb (60" table)
Weight Capacity	Capacity: 300 lb
Frame	18-gauge stainless steel construction. Thermofoil over electrostatic, powder-coated painted steel (woodgrain and metallic colors); electrostatic, powder-coated painted steel (solid colors)
Doors and Drawers	Sealed ¾" medium-density fiberboard with PVC thermofoil and electrostatic, powder-coated paint
Configurations	<ul style="list-style-type: none"> – Stainless steel or laminate surface – Left- hand or right-hand layout – Multiple cabinet door and drawer options – Multiple cabinet colors and surface options – Peninsula-style or island configurations



- 01** Durable 18-gauge steel frame protects against rust, while sealed cabinet panels resist moisture and chemical damage.
- 02** Seamless polystyrene drawers with rounded corners and one-piece construction allow drawers to contain spills.
- 03** Modular drawer and door layouts support supply organization.

ORDERING GUIDE

For quick ordering, explore the many pre-designed table options in the Midmark Equipment Catalog, or contact your Midmark sales representative to custom-configure the right solution for your space and care needs.



Scan QR code to view the Midmark Equipment Catalog.

*EPA registration numbers 85353-2 for pulls and 84542-7 for surfaces.

© 2026 Midmark Corporation, Versailles, Ohio USA. All rights reserved.
007-10550-00 Rev A1 (1/26)

