



- Energy efficient 22-watt power consumption
- 21.5" monitor available
- 5-year limited warranty

ADVANCED LIGHTING DESIGN

Midmark lights are engineered to deliver the true, reduced-shadow white light you need to properly match shades, identify details and diagnose tissue. The automatic on/off function and up-front touchpad control system make the lights easy to use. Midmark lights are available in various mounting configurations, including the Track Light Monitor, to meet your unique needs.



LED OPERATORY LIGHT

Offers energy efficiency and delivers true, reduced-shadow white light that allows you to properly identify shades, see details and evaluate tissues.



TRACK LIGHT MONITOR

Combines excellent dental lighting with highly flexible monitor positioning, enhancing the patient experience and easing patient acceptance of procedures.



MIDMARK 255 LED PROCEDURE LIGHT

Delivers the superior white, shadow-free illumination you need for procedures whether you practice in an office environment or specialized procedure area.

LIGHTING SPECIFICATIONS

Dental Lights	LED Operatory Light	Midmark® 255 LED Procedure Light
Operate Mode Color Temperature	4,500 Kelvin 4,750 Kelvin 5,000 Kelvin	4,400 Kelvin
Operate Mode Intensity, Foot Candles (lx)	1,400 fc (15,000 lx) 1,950 fc (21,000 lx) 2,500 fc (27,000 lx)	7,500 fc (80,700 lx)
Composite Mode Color Temperature	5,700 Kelvin	Not Applicable
Composite Mode Intensity, Foot Candles (lx)	550 fc (5,900 lx)	Not Applicable
Auto On/Off	✓	Not Applicable
Operate (Composite-Safe) Mode	✓	Not Applicable
Composite Mode	✓	Not Applicable
Touchpad Controls	✓	✓
Low-Wattage Power Consumption	✓	✓
LED Technology	✓	✓

MOUNTING OPTIONS



Left/Right



Console



Cabinet



Wall



Ceiling



Track



Track Light Monitor



Universal



Single Ceiling Mount



Dual Ceiling Mount



Wall Mount



Mobile

For more information, contact your Midmark dealer or call 1.800.MIDMARK.
Outside the USA, call 1.937.526.3662 or visit our website at midmark.com

© 2023 Midmark Corporation, Versailles, Ohio USA.
Products subject to improvement changes without notice. Litho in USA.

007-10046-00 Rev E1 (4/23)