



DESIGNED FROM THE INSIDE OUT TO SUPPORT THE WAY YOU WORK

Whether you need a compact delivery system like the Asepsis 21®, the flexibility of the Procenter® or the enhanced ergonomics of the Elevance®, you can rely on our proven technology. At the heart of each Midmark dental delivery system is our patented Kink Valve® technology, which offers an independent control system for each handpiece for an unprecedented level of control with minimal maintenance. Midmark delivery systems are available in multiple configurations with accessories and options to meet your unique needs.



ELEVANCE CONTINENTAL INSTRUMENT DELIVERY SYSTEM

For over-the-patient instrument delivery designed with ergonomic balance, look no further than the Elevance Continental Delivery System. Its Integration Control Module helps you effortlessly adjust or return settings to a desired range, and 10 fully customizable presets let you tailor your air and water adjustments for up to three separate users. Integrate accessories in up to five positions in addition to the syringe.



ASEPSIS 21 INSTRUMENT DELIVERY SYSTEM

The Asepsis 21 provides reliable delivery in a compact, value-priced package. This unit's surfaces are smooth and easy to clean, and its air and water controls are located beneath a removable top to help protect them from contamination. Integrate accessories in up to five positions in addition to the syringe.



PROCENTER INSTRUMENT DELIVERY SYSTEM

The Procenter is a flexible instrument delivery unit that provides line-of-sight access to air and water controls and fine needle-valve adjustments through an exclusive pop-up lid. Keep your instruments within reach by integrating your scaler, curing light, intraoral camera, air-driven handpieces and up to two electric motors. Integrate accessories in up to five positions in addition to the syringe.

DELIVERY UNIT SPECIFICATIONS

Delivery	Elevance®	Procenter®	Asepsis 21®
Left-Right Delivery Option	✓	✓	✓
Console Delivery Option	✓	✓	✓
12:00 Rear Delivery Option		✓	✓
Continental Unit Available	✓		✓
Assistant's Unit Available	✓	✓	✓
5-Year Limited Warranty	✓	✓	✓
Kink-Valve® Technology (with 10-year limited warranty)	✓	✓	✓
Chair Control Option	✓	✓	✓
Handpiece Integration	✓	✓	
3-Way Syringe Included	✓	✓	✓
Air-Lock Flex Arm		✓	
Programmable Presets	✓		
Touch-Screen Controls (Available in English, Spanish and French languages)	✓		
Needle-Valve Air/Water Adjustments		✓	✓

AVAILABLE ACCESSORIES:

Air Driven Options

- Handpiece Illumination System
- Tubing and Lamp Fiber Optic Handpiece (5-hole)
- Tubing for ISO-C Fiber Optic Handpiece (6-pin)

Intraoral Cameras

(Not available on Asepsis 21)

- SOPRO® Camera USB2 Docking Station
- SOPRO® 617 Camera
- SOPRO® 717 First Camera
- SoproLIFE® Camera
- SoproCARE® Camera

Bien-Air Electric Handpiece Solutions

(Not available on Asepsis 21)

- Bien-Air® MX2 Plus LED, Single Motor System
- Bien-Air® MX2 Plus LED, Dual Motor System
- Bien-Air® CA 1:1 NOVA Microseries Low-Speed Handpiece
- Bien-Air® CA 1:5 NOVA Microseries High-Speed Handpiece
- Bien-Air® HP 1:1 Straight Handpiece

Curing Light

(Not available on Asepsis 21)

- SATELEC® Mini LED Curing Light

Ultrasonic Scalers

- Dentsply Cavitron® Steri-mate® 360 Ultrasonic Scaler (G139)
- SATELEC® Pure Newtron® SLIM Ultrasonic Scaler
- SATELEC® Pure Newtron® SLIM B.LED Ultrasonic Scaler

-
- 10-year limited warranty on kink valves
 - 5-year limited product warranty

Cavitron, Steri-mate, Bien-Air, SATELEC, Newtron, SOPRO, SoproLIFE and SoproCARE are not trademarks of Midmark Corporation and are the sole property of their respective owners.

For more information, contact your Midmark dealer or call 1.800.MIDMARK. Outside the USA, call 1.937.526.3662 or visit our website at midmark.com

© 2022 Midmark Corporation, Miamisburg, Ohio USA
Products subject to improvement changes without notice. Litho in USA.
007-10047-00 Rev D1 (2/22)

