


DENTAL

Midmark PowerVac Dry Vacuum



SPECIFICATIONS

Model/Order No.	P3 (SINGLE)	P5 (SINGLE)	P6 (TWIN)	P7 (SINGLE)	P10 (TWIN)	P14 (TWIN)
No. Users	3-5	5-7	6-10	7-10	10-14	14-20
No. Vacuums	1	1	2	1	2	2
Total Horsepower	2	3	6	3	6	6
Vacuum (inches Hg)	up to 18	up to 18	up to 18	up to 18	up to 18	up to 18
Voltage	208-230 V 60 Hz	208-230 V 60 Hz	208-230 V 60 Hz	208-230 V 60 Hz	208-230 V 60 Hz	208-230 V 60 Hz
Total Amps	11.4	11.4	11.4 x 2	11.4	11.4 x 2	11.4 x 2
Circuit Breaker Amps	30	30	30 x 2	30	30 x 2	30 x 2
Dimensions Overall (H x W x D)	58" x 22" x 22"	58" x 22" x 22"	58" x 56" x 22"	58" x 22" x 22"	58" x 56" x 22"	58" x 56" x 22"
Base (H x W x D) Tank (H x W x D)	25" x 22" x 20" 34" x 22" x 22"	25" x 22" x 20" 34" x 22" x 22"	25" x 22" x 20" 34" x 22" x 22"	25" x 22" x 20" 34" x 22" x 22"	25" x 22" x 20" 34" x 22" x 22"	25" x 22" x 20" 34" x 22" x 22"
Shipping Weight	Base - 402 lb Tank - 90 lb	Base - 402 lb Tank - 90 lb	Base - 402 lb x 2 Tank - 90 lb x 2	Base - 402 lb Tank - 90 lb	Base - 402 lb x 2 Tank - 90 lb x 2	Base - 402 lb x 2 Tank - 90 lb x 2