



Sterilization Compliance Simplified

Adding a data logger to your Midmark M9 or M11 Steam Sterilizer streamlines and simplifies documentation of the physical and mechanical performance (sterilization time and temperature) of every sterilization cycle.

Now you can effortlessly record important data to help maintain compliance. Save time during audits and easily trouble-shoot sterilizer chemical or biological indicator failures.

Simple USB enabled electronic records can now be collected, transferred and organized as the user desires. Important information regarding the temperatures and sterilization times achieved for every sterilization cycle is at your fingertips.

The Midmark Sterilizer Data Logger offers many important features, including:

- High capacity, detailed sterilizer data storage
- Compatible with all M9 and M11 steam sterilizer models
- Reliable USB based data transfer

FEATURES



High-capacity, USB-transferable data storage

The data storage capacity of the Sterilizer Data Logger handles years of operation. Data is automatically loaded and backed up for safe transfer and easy access.



Compatible with all M9 and M11 Sterilizers

Designed to fit all models of M9 and M11 sterilizers, the Sterilizer Data Logger expands the documentation capabilities of your M9 or M11 sterilizer.

SPECIFICATIONS

Data Logger Compatibility	All M9 and M11 Steam Sterilizer models
Added Height to Sterilizer Unit	1.44 inches (3.66 cm)
Data Storage Capacity	>10,000 cycles; 8 gB
USB Storage Capacity	>10,000 cycles; 4 gB
USB Compatibility	USB 2.0/3.0
Computer Software Requirements	Windows ME, XP, 2000 or more recent versions

9A682001 Midmark Sterilizer Data Logger for M9 and M11 Sterilizers; 002-10502-00 Data Logger USB Device