



BETTER VERSATILITY RIGHT AT YOUR FINGERTIPS

With changing technology comes different needs. Midmark Workstations offer a variety of options for the integration of technology in the dental space. Configure your workstation the way you need it, where you need it.


Midmark Workstations are designed to improve flexibility and efficiency, allowing you to work from a variety of ergonomically correct positions while maintaining eye contact with patients.



01



02



03

- 01** Height adjustment feature brings technology to the point of care.
- 02** Two total-locking casters in the front lock both the wheel and swivel simultaneously.
- 03** Ample work surface to transport laptop and supplies room to room.

Laptop, monitor, keyboard and mouse not included.

MOBILE WORKSTATION SPECIFICATIONS

Model/Order No.	6212	6215	6215-002	6215-003
Worksurface Dimensions	23" D x 26" W	21" D x 26" W	17.5" D x 23.5" W	17" D x 25" W
Base Dimensions	22.5" D x 20.5" W	19.5" D x 20.5" W	19.5" D x 20.5" W	19.5" D x 20.5" W
Weight Limit	30 lb.	20 lb.	20 lb.	20 lb.
Monitor Arm (Model 6212)	75 mm and 100 mm VESA mounting patterns 9.5" independent vertical adjustment + up to 40 degrees positive tilt			
Casters	3.94" (100 mm) oversized casters			
Height Adjustment	23" to 49"			
Standards	ANSI/BIFMA Compliant			



6212



6215



6215-002



6215-003