

Explore the next generation of Midmark® sterilizers

Designed for better:

- DURABILITY
- WORKFLOW
- COMPLIANCE



Increase Longevity

25,000 CYCLES
Gain 2.5-times longer cycle life than previous models.

Reduce Lifetime Cost

UP TO US\$11,800 IN SAVINGS
Save thousands compared to competitive models with a sterilizer designed to require less frequent routine technician-performed service.

Enhance Workflow



3 LANGUAGES
Select English, French or Spanish for your preferred display language.



LIGHTER, LARGER TRAYS
Increase usable space with a larger but lighter tray than previous models.



5-INCH TOUCHSCREEN
Navigate setup with a clear, graphical display you can utilize while wearing gloves.

Improve Efficiency

AUDIO + VISUAL INDICATORS
Determine cycle status—cycle in progress, cycle complete and cycle error—from a distance with a progressing color-coded LED light bar, large cycle-countdown clock and audible cues.



CYCLE IN PROGRESS

CYCLE COMPLETE

CYCLE ERROR



TIME-SAVING DOCUMENTATION
Help save time, simplify teammate training and build instrument processing confidence with automated mechanical monitoring and chemical indicator reporting.



OPTIONAL AUTOFILL + AUTO-DRAIN
Get fresh water for every cycle without burdening your team.

Reduce Stress



ONSCREEN GUIDANCE
Improve safety with touchscreen prompts that lead your team through a standardized process.



ON-DEVICE TUTORIALS
Reduce the burden of training with step-by-step tutorials.



STREAMLINED COMPLIANCE RECORDKEEPING
Stay audit-ready with device reminders and notifications, user authentication, unlimited storage of routine care events for the life of the sterilizer, and automated cycle recordkeeping.



MORE CYCLE OPTIONS
Simplify clinical workflow with dynamic-air-removal steam sterilizers that provide unique preset cycles to meet both the CDC and AAMI sterilization cycle recommendations for wrapped instruments.



WRAPPED
270°F, 4 MIN



DELICATE
250°F, 30 MIN

NEW



WRAPPED
275°F, 3 MIN



UNWRAPPED
270°F, 3 MIN



5 CUSTOM CYCLES*
250°F–275°F, 3 MIN–45 MIN

*These cycles are not FDA-cleared, and validation of the sterility of items processed using a custom cycle is the responsibility of the user.

Simplify sterilization and get your team back to the work they love. Choose technology designed with them in mind.

[LEARN MORE >](#)



© 2025 Midmark Corporation, Versailles, Ohio USA
MKF00573

midmark®