

# Bien-Air® Lubricare 2®

## Dental Handpiece Care + Maintenance System

### LUBRICARE 2

Quickly perform dental handpiece maintenance for up to four handpieces at a time with the Lubricare 2 system, an essential addition to your instrument processing workflow. The clean air-blow system ensures thorough cleaning and lubrication inside the handpieces and allows for consistent, repeatable high standards without the variability that may occur during manual cleaning and lubrication.

### SPRAYNET® + LUBRIFLUID®

Use Spraynet and Lubrifluid for cleaning and lubrication to maintain impeccably hygienic instruments with an unrivaled service life.



## FAST, EFFECTIVE COST-SAVING MAINTENANCE



### 01 FAST

Clean a turbine in approximately 20 seconds and a contra-angle handpiece in 40 seconds with preprogrammed operating modes that enable automatic maintenance in the Lubricare 2 system.

### 02 EFFECTIVE

Achieve best cleaning and lubrication outcomes with separate dispensing of cleaner and lubricant for more precise application than systems that dispense the solutions together. Reduce 99.99% of the microbial load with Spraynet cleaning solution.



### 03 COST-SAVING

Reduce maintenance oil consumption with a more consistent process compared to manual maintenance procedures.

## ENHANCE YOUR INSTRUMENT PROCESSING WORKFLOW



### LUBRICARE 2® SPECIFICATIONS

<b>Compatibility</b>	E-type attachments (Additional couplings may be required)
<b>Cycle Duration</b>	20 s turbine, 40 s contra angle
<b>Capacity</b>	4 handpiece lanes (2 cleaning, 2 lubricating) and 1 chuck lane
<b>Rating</b>	AC 100 V–240 V 50–60 Hz
<b>Power Consumption</b>	10 VA–25 VA
<b>Air Flow Rate</b>	60 NL/min–100 NL/min (0.3 MPa–0.5 MPa)
<b>Air Pressure</b>	43.5 psi–72.5 psi (0.3 MPa–0.5 MPa) 58 psi (0.4 MPa) recommended
<b>Weight</b>	Approx. 17 lb (8 kg)
<b>Outer Dimensions</b>	11.81" L x 11.81" W x 14.37" H (300 mm x 300 mm x 365 mm)

Bien-Air, Lubricare 2, Spraynet and Lubrifluid are not trademarks of Midmark Corporation and are the sole property of their owner.

For more information, contact your Midmark dealer or call 1.800.MIDMARK. Outside the USA, call 1.937.526.3662 or visit our website at [midmark.com](http://midmark.com)

© 2023 Midmark Corporation, Versailles, Ohio USA.

Products subject to improvement changes without notice. Litho in USA. 007-10414-00 Rev B1 (4/23)

