



DENTAL

Operator Cabinetry



Designing Better Care



We believe design is about not just the space, but the purpose of the space. Specialized care requires specialized solutions. Choose cabinetry designed for dentistry.



Choose Cabinetry Designed for Dentistry

Midmark offers two collections of clinical-grade cabinetry—Synthesis® and Artizan® Expressions—specifically designed for the rigors of the healthcare environment and crafted to fit your space, style and workflow.

06
BETTER INSIGHT
FOR BETTER CARE

12
CHOOSE CLINICAL
GRADE CABINETRY FOR
YOUR PRACTICE

18
MIDMARK PRODUCTS +
SOLUTIONS

28
TREATMENT STATIONS +
SIDE CABINETS

36
MIDMARK FINISHES +
DETAILS

42
ACCENT UPHOLSTERY
COLORS

08
OPTIMIZE
WORKFLOWS

14
COMPARE QUALITY

20
SYNTHESIS
CABINETRY

30
WALL-HUNG
CABINETS

38
CABINET FINISHES

43
RESIN PANELS

10
SUPPORT
MUSCULOSKELETAL
HEALTH

16
THINGS TO
CONSIDER

24
ARTIZAN EXPRESSIONS
CABINETRY

32
SPECIALTY
CABINETS

40
HANDLE + PANEL
STYLE OPTIONS

44
SEE WHAT OUR
CUSTOMERS ARE SAYING

26
CENTRAL STATIONS

34
MOBILE WORKSTATIONS

41
COUNTERTOP COLORS

46
DESIGN WITH MIDMARK



SCAN OR CLICK TO VIEW
OPERATORY CABINETRY

Better Insight for Better Care

Know how your equipment is performing so you can focus on the people who matter most.

Midmark Smart View™ is a cloud-based service* designed to maximize the value of Midmark products by providing

- real-time performance insights,
- proactive maintenance tracking,
- cycle records and reporting,
- and software updates

through the myMidmark portal.

With continuous enhancements and feature upgrades, Smart View can help eliminate the risks of outdated technology, while enhancing performance to support long-term efficiency and compliance.


** Some features may require an active Smart View plan*

EXPERIENCE THE NEW MIDMARK DENTAL

With new connected technology and enhanced operator, instrument processing and mechanical room solutions, the new Midmark Dental is helping set a new standard and advance what's possible in dental care. **Discover how we can design better care together.**



**SCAN OR CLICK TO LEARN
HOW TO CONNECT**



SEE THE WHOLE PICTURE

Access usage trends, compliance data and equipment status for your connected Midmark solutions from one secure dashboard, whether you manage one location or many.



ENHANCE PERFORMANCE

Connect your dental chairs, delivery systems and LED lights for ongoing software updates designed to enhance equipment performance and workflow.



STREAMLINE EQUIPMENT MANAGEMENT

Manage connected sterilizers across your organization to help improve workflow efficiency, simplify compliance and protect the lifespan of your equipment.



KEEP CARE MOVING

Receive routine maintenance notifications and critical alerts for your connected air compressors and dry vacuums to help minimize downtime and keep care moving.

Optimize Workflows



DESIGNED TO STREAMLINE YOUR PRACTICE


To support optimal workflow and ergonomic efficiency, clinical design best practices recommend allocating approximately 120 square feet per operatory—measuring 10 feet, 4 inches wide by 12 feet deep.

These dimensions support efficient cabinetry and equipment placement, provide maneuvering space around the chair and ensure ADA clearance compliance. They also help keep worksurfaces, instruments and the patient’s oral cavity within easy reach during procedures.

In contrast, smaller operatories can compromise efficiency by limiting worksurfaces, reducing storage capacity and positioning essential components out of reach—ultimately disrupting your clinical workflow.

PLAN FOR EFFICIENCY

Start with the dental chair as your anchor, then layer in essential clearances, worksurfaces and equipment. This dual-entry design supports ergonomic access, streamlined workflows and optimal use of space.



*Adjustable backrest: 77 in

Support Musculoskeletal Health





DESIGNED FOR BETTER ERGONOMICS

Reducing musculoskeletal pain and disorders starts with equipment and environments designed to work with the body, not against it. For better ergonomics, choose cabinetry that helps position supplies, work surfaces and technology where your team needs them most.

Side cabinets and treatment stations can help create convenient work zones near the point of care, while wall-hung cabinets with lower upper-cabinet heights and shallower depths can improve access to stored supplies. Optional angled shelving and organized bins help keep items visible and easier to retrieve.

Midmark clinical-grade cabinetry offers these design choices to help reduce unnecessary reaching, leaning, twisting and wasted motion, supporting a more ergonomic workflow.



01 A Midmark left/right doctor's delivery unit and rear chair-hub-mounted assistant's unit combine with a mobile cabinet and pivoting worksurface on the treatment station to help improve ergonomic reach and workflow.

Choose Clinical-Grade Cabinetry for Your Dental Practice

Midmark Synthesis® and Artizan® Expressions cabinetry collections are engineered specifically for the unique demands of dentistry, bringing together clinical performance, design flexibility and long-term reliability.

Durable materials and construction help resist the moisture, stringent cleaning and heavy use common in dental operatories, while asepsis-friendly surfaces support cleanability and infection prevention. With flexible configurations, colors and finishes, Midmark cabinetry can be tailored to support your workflow, reflect your brand and help create a comfortable, welcoming patient experience.

- 01 Integrated Electrical
- 02 Solid Surface Countertops
- 03 Thermally Fused Melamine Finish
- 04 Resin Finish
- 05 Upholstery Accent
- 06 Stainless Steel Base



SYNTHESIS® CABINETRY












ARTIZAN® EXPRESSIONS CABINETRY

Compare Quality

MIDMARK® CLINICAL-GRADE CABINETRY VS. CONSUMER-GRADE CABINETRY

Most generic or local cabinetry is not designed for the demanding dental environment. Consumer-grade cabinetry can delaminate over time, develop unsealed corners susceptible to bacteria and moisture, and warp from the humidity inherent in a dental practice.

	CABINET FRAME	PANEL SUBSTRATES	BASE MATERIAL		FINISHES	TYPES OF ASSEMBLY	DRAWERS	DRAWER SLIDES	HANDLES
MIDMARK SYNTHESIS® CABINETRY A steel foundation provides strength and durability.	 18-gauge cold-rolled steel - Modular design	 - LEED® compliant, Medium-density fiberboard 0.75 in panels	 - Thermofoil over electrostatic, powder-coat painted steel (woodgrain and metallic colors) - Electrostatic, powder-coat painted steel (solid colors) - Integrated adjustable levelers		 - PVC thermofoil (woodgrain and metallic colors) - PVC thermofoil and electrostatic, powder-coat paint (solid colors)	 Mechanical fasteners (screws, pop rivets) Tog-L-Loc® sheet metal joining system Mechanical fasteners and Tog-L-Loc provide reliable strength.	 Seamless, one-piece polystyrene Seamless polystyrene drawers with rounded corners contain spills and simplify cleaning.	 - Full-extension ball bearing drawer system - Steel mounting foundation - Soft-close feature	 - Integrated and insert options - Brushed nickel options - Antimicrobial option Synthesis pull options support an aseptic environment.
MIDMARK ARTIZAN® EXPRESSIONS CABINETRY Industrial-grade strength prolongs durability.	 0.75 in industrial-grade particle board - Customizable design that can be configured to specific requirements	 - Industrial-grade particle board 0.75 in panels	 - Seamless stainless steel - Integrated adjustable levelers		 - Factory pre-applied, thermally fused melamine 1.2 mm PVC edgebanding	 100% glue-dowel assembly Glue-dowel assembly adds strength to the frame.	 - Wood bottom and backs using factory, pre-applied thermally fused melamine - Metal slides with front panel adjustment and removal feature	 - Full-extension ball bearing drawer system - Soft-close feature	 Brushed nickel Brushed nickel pulls are easy to clean for better asepsis.
CONSUMER-GRADE CABINETRY* Millwork options are not specifically designed to endure the clinical dental environment.	 0.5 in to 0.625 in low-density particle board or plywood (typically) - Basic or limited design	 - Low-density particle board - Plywood	 - Basic plywood or particle board covered with laminate - Wood shims for leveling		 - High-pressure laminate – Not always properly applied – May use low-grade substrate material - Less than 2 mm edgebanding - Varnish, paint, or left unfinished	 Staple fasteners (frequently) The staple fasteners typically used in consumer-grade millwork are unreliable for clinical-grade cabinetry.	 Insides often varnished, painted or unfinished Seams and sharp inner corners can result in difficulty cleaning and a build-up of contaminants.	 - Single, undermounted runner (monorail glide—typically used) - Slides common in kitchens (less durable)	 - C-style shape - Plastic or aluminum Consumer-grade millwork pulls are not made to support an aseptic clinical environment.



Things to Consider



ERGONOMIC WORKFLOW

Is your storage designed for the clinicians who use it most?

Women represent about 96%¹ of dental hygienists and more than half of recent dental school graduates,² yet many cabinets place supplies out of comfortable reach. Overly tall or deep storage can make routine tasks harder on the body.

Consider wall-hung cabinetry designed to improve reach and visibility. Midmark Artizan® Expressions and Synthesis® wall-hung cabinets feature a lower upper-cabinet height and thinner depth to help reduce stretching and leaning.



DESIGN FLEXIBILITY

Does your cabinetry need to fit a unique room size, configuration or workflow?

Standard cabinetry layouts may not support the way every practice works. When space, processes or design goals vary, flexibility becomes essential.

Consider clinical-grade cabinetry with customizable options. Midmark Artizan Expressions Cabinetry offers layouts, configurations and finishes that can be tailored to your workflow, brand and overall practice design.



ROOM SIZE + LAYOUT

Does your layout support efficient movement and ADA clearances?

The right cabinetry plan can help reduce wasted motion, improve efficiency and support equipment integration, ergonomic movement and ADA considerations. But every design decision can affect how the room functions.

Consider guidance from an experienced dental practice designer. Midmark Design Experts can help create beautiful, functional operatories that support workflow and enhance patient comfort.



ASEPSIS + DURABILITY

Is your cabinetry designed for the demands of the dental environment?

Moisture, stringent cleaning, equipment integration and daily use can take a toll on cabinetry. In clinical settings, consumer-grade materials may be more prone to warping, delamination and bacterial growth.

Consider clinical-grade steel cabinetry. Midmark Synthesis Cabinetry combines steel-on-steel construction, nonporous seamless surfaces, modular flexibility and asepsis-friendly design options for better cleanability and long-term performance.



MOBILE TECHNOLOGY

Do you share technology, supplies or accessories between operatories?

When tools need to move from room to room, the right mobile solution can help reduce setup time and keep essential items close to the point of care.

Consider sturdy carts and mobile workstations. Midmark mobile cabinets and workstations help you transport technology and supplies where they are needed, supporting flexible room use and efficient workflows.



THE PATIENT EXPERIENCE

Have you considered how cabinetry finishes can affect the patient experience?

Research suggests color, texture and tone can impact patient mood and shape perceptions of the overall healthcare experience.³ With careful selection, these design options may help reduce stress, enhance comfort and improve patient satisfaction.

Consider cabinetry with colors and finishes selected to support a calm, comfortable environment. Midmark Synthesis and Artizan Expressions cabinetry collections offer soft, cool tones, balanced neutrals and natural wood finishes to help shape operatories that feel welcoming, soothing and beautifully composed.

1, 2, 3 See back page for details and sources

Midmark® Products + Solutions




SCAN OR CLICK TO
VIEW PRODUCTS

Midmark® Synthesis® Cabinetry

DESIGNED FOR DECADES OF CLINICAL PERFORMANCE

Steel-on-steel construction, seamless surfaces and modular configurations help Synthesis cabinetry deliver lasting durability, cleanability and flexibility that supports efficient dental workflows and long-term performance. Synthesis cabinetry can also contribute to earning credits in the LEED® Green Building Rating System.



01 BINS + DIVIDERS

Bins and dividers help improve organization so you can easily find the supplies you need.

02 FLOW SHELVING

Gravity-fed, angled flow-shelving improves visibility and access to supplies, even for items stored at the back of the cabinet.



03 THERMOFOIL PANELS

Seamless thermofoil panels are designed to resist spills, chemical degradation and delamination, unlike common consumer-grade products.

Supporting a Healthier Environment

SYNTHESIS® CABINETRY



The unique product attributes of Midmark Synthesis cabinetry have been designed to assist you in obtaining credits within the LEED® family of rating systems. We choose to be proactive because we know that LEED contributes to healthy environments through improved air quality, and we want to be your partner in the facilitation of a safer working environment.

DESIGN A HEALTHY ENVIRONMENT TOGETHER

Midmark Synthesis cabinetry can contribute to specific strategies and LEED credits that can impact building certification including water efficiency, construction waste management, material reuse, recycled content, local/regional materials and low-emitting materials. See the table for a comprehensive list of credits.

LEED® v5 Credit Mapping

Credit	Contribution	Possible Points
MRc1 Interior Materials Reuse	ID+C: Commercial Interior Modular, reconfigurable cabinetry supports on-site reuse within the same project, supporting adaptable healthcare environments and reducing the need for new materials.	1-4
MRc1 Interior Materials Reuse	ID+C: Commercial Interior Modular, reconfigurable cabinetry supports off-site reuse in new project locations, supporting material circularity by integrating salvaged cabinetry into evolving clinical environments, extending product life and reducing demolition waste and new manufacturing demand.	1-3
MRc3 Low-Emitting Materials	BD+C: New Build + ID+C: Commercial Interiors Low-emitting cabinetry components support healthier, human-centered clinical environments by reducing VOC exposure and improving indoor air quality for patients, clinicians and staff.	1-4
MRc4 Building Product Selection and Procurement	BD+C: New Build + ID+C: Commercial Interiors Select cabinetry components enable informed, sustainability-driven material selection through transparent life cycle impact data, supporting project-level environmental performance goals.	1-5
MRc5 Waste Diversion	BD+C: New Build + ID+C: Commercial Interiors Recyclable cabinetry packaging supports waste diversion by facilitating efficient sorting and reducing construction waste sent to a landfill.	1-3
WEc2 Water Efficiency	BD+C: New Build + ID+C: Commercial Interiors Integrated sensor-operated, low-flow fixtures enhance water efficiency by reducing potable water use and supporting project-level water reduction strategies.	1-6
EQc3 Accessibility + Inclusion	BD+C: New Build + ID+C: Commercial Interiors Cabinetry and sink configurations support inclusive, human-centered clinical environments by enhancing usability, accessibility and workflow efficiency for diverse patient populations and clinical staff, fostering comfort, dignity and a sense of belonging beyond minimum code requirements.	1-2

Midmark® Artizan® Expressions Cabinetry

DESIGNED FOR STYLE, STRENGTH + FLEXIBILITY

Blending modern aesthetics with clinical functionality, Artizan Expressions cabinetry helps create an organized, inviting operatory that supports patient comfort and efficient care delivery. Customizable configurations and finishes help tailor the space to your workflow, brand and design vision.



01 LIFT-UP DOORS

Lift-up doors provide greater clearance in tight spaces and improve visibility to all upper contents when open.

02 SOFT-CLOSE HINGES

Concealed Blum® soft-close steel hinges adjust easily, open 110° and prevent damaging door slamming.

03 MOISTURE RESISTANCE


Non-corroding brushed nickel handles and 1.2 mm PVC edgebanding, which seals panel edges and corners from moisture, support durability and asepsis.

Central Stations

SYNTHESIS® CABINETRY



01 PASSTHROUGH CABINET
Make efficient use of your space with flexible passthrough options for your intraoral X-ray solutions.




02 DUAL-SIDE STORAGE
Save valuable floor space with Synthesis Central Stations ideal for multi-room plans, providing storage on both sides for the dentist and assistant.




03 FULL-EXTENSION DRAWERS
Improve durability with full-extension ball-bearing drawer slides and slide attachments designed for better withdrawal strength.


ARTIZAN® EXPRESSIONS CABINETRY



01 CONCEALED STORAGE
Bring a modern look to the operator by keeping clinical supplies tucked out of view of patients yet easily accessible by the dental team.




02 PASSTHROUGH CABINET
Share one intraoral X-ray between two operator entries with an Artizan Central Station.



03 SHARED SINK
Optimize space with one sink shared between two operatories and add sophistication with sleek upholstered panels.


Treatment Stations + Side Cabinets

SYNTHESIS® CABINETRY



- 01 SEAMLESS DRAWER INTERIORS**
Seamless, one-piece polystyrene drawers with rounded corners contain spills and simplify cleaning.
- 02 STRONG DRAWERS**
Synthesis drawers are engineered for enhanced strength and durability.
- 03 SLIDE-OUT STORAGE**
Slide-out storage can be prepped with commonly needed equipment, keeping clinicians in the operatory and making your equipment easily accessible during procedures.

ARTIZAN® EXPRESSIONS CABINETRY



- 01 FULL-EXTENSION DRAWERS**
Full-extension soft-close drawers can withstand the demands of the dental practice and help promote a calmer environment by preventing slamming drawers.
- 02 TOE KICK**
A 4.25 in tall toe kick allows clinician seating to roll underneath cabinet doors and drawers while preventing damage to the cabinetry.
- 03 CUSTOMIZABLE OPTIONS**
Create your own style with customizable options like Ultraleather® upholstered panels for your treatment station.

Wall-Hung Cabinets

SYNTHESIS® CABINETRY



Designed around the clinicians who use cabinetry most, Synthesis wall-hung cabinets help bring supplies into clearer view and easier reach. With women representing 96% of dental hygienists¹ and more than half of recent dental school graduates², Synthesis supports the needs of today's dental teams with a lower upper-cabinet height and shallower depth. Angled flow shelving and stackable bins help simplify supply retrieval and support more efficient workflows.



- 01 ANGLED FLOW SHELVING**
Gravity-fed angled flow shelving helps keep supplies visible and accessible, even at the back of the cabinet.
- 02 WALL-TO-WALL ACCESS**
Wall-to-wall floor access supports cleaning and disinfection to help maintain a safe, aseptic environment.

ARTIZAN® EXPRESSIONS CABINETRY



Create a streamlined operatory with wall-hung storage designed to complement your space and workflow. Artizan Expressions wall-hung cabinets help keep supplies close while supporting a clean, modern look and flexible design.

- 01 ACCESSIBLE STORAGE**
Wall-mounted side cabinets provide accessible storage and a secondary work surface for added efficiency at the point of care.
- 02 WALL-TO-WALL ACCESS**
A wall-hung design leaves the floor clear for wall-to-wall cleaning and better asepsis.

Specialty Cabinets

SYNTHESIS® CABINETRY

01 MOBILE PROCEDURE CART
A no-tip base design enables sturdy room-to-room maneuverability, and flexible options allow you to create the cart you need.

- A Transit Handle
- B Disinfectant Wipe Bracket
- C Single Glove Holder
- D 6 x 6 Storage Bins
- E Soft Edge Contour Top with Raised Corners
- F Sharps Bracket
- G Waste Basket Bracket



ARTIZAN® EXPRESSIONS CABINETRY

- 01 MOBILE TREATMENT CABINET**
Bring worksurfaces and supplies within ergonomic reach for teammates of all sizes.
- 02 ORTHO DELIVERY CABINET**
Select a cabinet designed for the unique storage and workflow requirements of your orthodontic practice.



02



02 MOBILE TREATMENT CABINET
Make the most of the space you have. Bring the supplies you need to the point of care with mobile cabinets that optimize space.

03



03 WASTE CART
Available in 31 in and 33 in heights, waste carts have a no-touch front design with front lock for easy waste removal.

04



04 LINEN CART
Available in 31 in and 33 in heights, linen carts can assist with the removal of dirty linens, especially from your oral surgery rooms.

02



Mobile Workstations

DESIGNED TO MOVE WITH YOUR WORKFLOW

Enhance patient care through better ergonomics and versatility.
Ensure your technology is where you need it when you need it.



01

01 6282 WALL MOUNTED WORKSTATION
This wall-mounted workstation is designed with independent height adjustment for standard PC based technology components.

02 6215 COMPACT WORKSTATION
This workstation is ideal for smaller work areas and for wireless users who carry their laptop and supplies room to room, yet need a work area in the operatory space.

03 6212 FLAT PANEL PC
This flat panel PC workstation brings technology to the point of care with a convenient mount for a flat panel monitor.



02



03

Midmark[®] Finishes + Details

"Strong yet stylish, the steel Midmark cabinetry line [Synthesis] was the right choice for my new office. With well-thought-out functionality in a contemporary design, it shows every patient that our office invests in the absolute best for their care."

Dr. Brian S. Nylaan

Grand Rapids, MI















SCAN OR CLICK TO
VIEW ARTIZAN EXPRESSIONS
STYLE GUIDE



SCAN OR CLICK TO
VIEW SYNTHESIS STYLE GUIDE












Synthesis[®] Cabinet Finishes

Synthesis[®] door and drawer panels are sealed with membrane-pressed thermofoils in woodgrain and metallic colors. Thermofoil is a dimensionally flexible material that provides a continuous seamless surface, eliminating the need for panel edgebanding and offering a high level of scratch and stain resistance.

				
Mist PVC Free 980	Venus Silver (metallic) 804	Flat Iron (metallic) 807	Storm 802	Riverstone 967
				
Pearl Essence (metallic) 805	Ember 968	Tundra 966	Frost 791, Frost PVC Free 931	Fog PVC Free 969
				
Pebble Grey 216, Pebble Grey PVC Free 793	Sandbar 963	Prairie 622	Dune 794	Flax 795
				
Oat 611	Russet 964	Timber 801	Sable 965	

Artizan[®] Expressions Cabinet Finishes

Available in a variety of Arauco[®] colors, Artizan Expressions uses thermally fused melamine panels. These panels are created at high temperature and pressure to form a durable, peel-resistant finish. Matching edgebanding is applied to completely seal the corners from moisture.

				
White 337	Silver Frost 757	Pewter Frost 825	Charcoal 750	Black 312
				
Bergen Ash 970	Presten Ash 978	Sahalie Pine 973	Shadow Frost 979	Pelee Island Pear 975
				
Karuna Ash 971	Talas Cherry 974	Licorice Groovz 977	Koosha Pine 972	Sable Glow 976


Handle + Panel Style Options

SYNTHESIS® CABINETRY

Handle Styles




Panel Styles



*EPA Registration Number 85353-2
**Available in Pebble Grey or Pebble Grey PVC Free only

ARTIZAN® EXPRESSIONS CABINETRY

Handle Style



Panel Style








Artizan Standard

Countertop Colors

Solid Surface (SS) is a unique material made from an advanced blend of natural minerals and acrylic polymer. It has proven itself to be a remarkably durable, versatile material that’s an ideal surface for dental environments.

Quartz (Q) is virtually maintenance free and, because it is one of nature’s strongest materials, it is also incredibly durable, exceptionally tough and scratch resistant.

				
SS - Sparkling White	SS - Antarctica	SS - Arrowroot	SS - Silver Birch	SS - Lava Rock
				
SS - Neutral Aggregate	SS - Dove	SS - Ash Aggregate	SS - Juniper	SS - Weathered Concrete
				
SS - Linen	SS - Clam Shell	SS - Sandstone	SS - Hazelnut	SS - Willow
				
SS - Canyon	SS - Deep Night Sky	Q - Aria	Q - Everest	Q - Minuet
				
Q - Carrara Delphi	Q - Chateau	Q - Moon Gaze	Q - Moorland Fog	Q - Umber






SS = Solid Surface
Q = Quartz

Accent Upholstery Colors

Complete your distinctive design with the same Ultraleather® upholstery you've chosen for your treatment chairs, stools and side chairs. Match accents on the slim profile midsection of your treatment station and vanity front of your central station. Upholstery is available in your choice of many standard colors and additional colors by special order.




				
Mica (pearlized) 398	Pewter (pearlized) 618	Lake (pearlized) 749	Midnight (pearlized) 326	Milkweed 395
				
Papyrus 227	Cashmere 368	Fudge 751	Arctic 244	Stone 871
				
Eucalyptus 371	Riviera 400	Raven Wing 228	Charcoal 750	Dove Grey 372

				
Silver Pearl 732	Sassy 731	Apricot 262	Red 730	Nile 369
				
Admiral 401	Bayou 722	Cyan 723	Sonora 436	Parrot 729

Resin Panels

Select an artful resin panel to overlay the front of your treatment station's upper cabinet for a distinctive style all your own. Lumicor® EcoShield panels are made of Lumiform, a PETG resin that is 40 times stronger than glass, impact resistant and comprised of a minimum of 40% recycled content. This resin fully encapsulates hand-laid botanicals, organics, woven textiles or graphics.

				
Copper Spun 809	Silver Spun 808	Ivory Taipei 771	Natural Leaves 781	Coastal Light 786
				
Kauri 788	Nido 770			

See What Our Customers Are Saying

"It has almost been ten years since we had our cabinets from Midmark installed, and they still look and perform just as well as when they were first installed! We loved being able to go to the Midmark showroom and see all of their products in work settings. It clearly demonstrated to us our options and allowed us to see and choose exactly which products would best suit our needs. We wouldn't go anywhere else."

Dr. Mark Armstrong
Armstrong Dentistry, Troy, OH

"I really like how modern the Artizan Expressions look can be. It allows a lot of versatility to customize a space to your own style."

Dr. Jenna Triggiani
Two Rivers Family Dentistry,
Lisbon, ME

Design with Midmark®

DESIGNED FOR YOU

Our design experts are ready to help you personalize and equip your entire clinical environment—from the layout of every room to each sterilizer, chair, door and drawer, down to the smallest detail. Let us help you design the best environment for your practice.



CONNECT VIRTUALLY FOR LIVE DESIGN

Connect with Midmark Live Design services to bring your vision to life. We help you visualize your space, understand how it flows and spot potential challenges so you can make confident decisions to support better care.



**SCAN OR CLICK
TO SCHEDULE AN
APPOINTMENT.**



JOIN US IN PERSON FOR A MIDMARK EXPERIENCE TRIP

Explore equipment, finishes and floor plan options during a hands-on Midmark Experience trip. With expert guidance and a clear picture of your options, you'll be able to make more confident, informed decisions for your practice.



**SCAN OR CLICK TO
REQUEST A VISIT**



STREAMLINE YOUR DELIVERY + INSTALLATION

Ensure a smooth, on-time setup with Midmark Delivery Services. Our team coordinates the shipping, delivery, staging and setup of your Midmark equipment to help keep your project on track.



**SCAN OR CLICK TO LEARN
MORE ABOUT MIDMARK
DELIVERY SERVICES**

SOURCES

- 1 CareerExplorer, "Dental Hygienist Demographics in the United States." <https://www.careerexplorer.com/careers/dental-hygienist/demographics/>
- 2 ADA Health Policy Institute, "Dental School Grads and Advanced Program Enrollment by Gender, 2021-22." https://www.ada.org/-/media/project/ada-organization/ada/ada-org/files/resources/research/hpi/dental_school_grads_gender_2021-22.pdf?rev=6654eee768e74b8593d4b729dc1c6c95&hash=7B9AA598FDCEACDB9A29B623D6A23557
- 3 Tabanejad, Z. "The Impact of Color in Healthcare Environments: A Systematic Review and Research Agenda." Shiraz E-Medical Journal, vol. 25, no. 12, 2024. <https://brieflands.com/journals/semj/articles/146676.pdf>



Designing better care.®



CARB 93120.2 Phase 2 Compliant and TSCA Title VI Compliant

The color examples shown are the best representation of the original material. Actual color may vary slightly. We strongly recommend that you contact Midmark Customer Experience at 1.800.MIDMARK to request a sample before placing your order.

No product can guarantee points. Points are cumulative for all effected materials for a LEED® project. For further information on LEED, please visit the USGBC website: usgbc.org
LEED is a registered trademark of the US Green Building Council.

USGBC® and the related logo are trademarks owned by the U.S. Green Building Council® and are used with permission.

Ultraleather, Tog-L-Loc, Arauco, Formica and Lumicor are not trademarks of Midmark Corporation and are the sole property of their respective owners.

Midmark is an ISO 13485 and ISO 9001 Certified Company. Certain products are not included. See the complete list at: midmark.com/ISO

For more information, contact your authorized Midmark distributor or call: 1.800.MIDMARK Outside the USA call: 1.937.526.3662 or visit our website: midmark.com

© 2026 Midmark Corporation, Versailles, Ohio USA.
Products subject to improvement changes without notice.

007-10010-00 Rev H1 (7/26)

