



DENTAL

# Imaging Solutions



# Designing Better Care



We believe who's using our equipment and technology is every bit as important as how and why it's being used. Enhance your clinical outcomes with imaging equipment designed to be easy, comfortable and efficient to use to meet your needs and the needs of your patients.



# Diagnose with Confidence

## You can treat only what you can see.

Midmark® imaging solutions are designed to improve patient comfort while helping you see more, do more and provide better care. With the help of digital imaging, you can more clearly identify problems and provide accurate treatment early before those issues become serious.

06  
DESIGN FOR SAFETY

08  
DESIGN FOR  
ERGONOMICS

10  
THINGS TO CONSIDER

12  
MIDMARK PRODUCTS +  
SOLUTIONS

14  
MIDMARK PREVA  
INTRAORAL X-RAY

17  
ACCESSORIES

18  
DESIGN WITH MIDMARK



SCAN OR CLICK TO VIEW  
IMAGING SOLUTIONS

# Design for Safety

While radiation from dental imaging may be fairly low compared to everyday environmental exposures, the cumulative effects of radiation are well documented.

Help ensure patients receive radiation “as low as reasonably achievable” (ALARA) by following ADA and American Academy of Oral and Maxillofacial Radiology (AAOMR) recommendations<sup>1</sup> for patient-specific imaging. These guidelines emphasize evaluating each patient’s condition and diagnostic needs to determine appropriate imaging while minimizing radiation exposure.

## USE THE RIGHT DOSE.

Support safe, effective imaging by following a patient-specific, risk-based approach.

- Define the diagnostic need
- Evaluate risk, history and clinical findings
- Consider age and dentition stage
- Review prior imaging before exposure
- Limit exposure to the area of interest

**Apply the ALARA principle—select the lowest dose that achieves diagnostic quality.**




Proudly designed,  
assembled + supported  
in the USA.



# Design for Ergonomics

Improve ergonomics with an intraoral X-ray you can position effortlessly with one hand and rely on to keep that position without drifting.

The **Midmark Preva Intraoral X-Ray** system is designed with multiple points of adjustment and a distinctive ergonomic handle for effortless movement and fast one-handed positioning. The factory-balanced arm mechanics and proprietary dual-arm braking system help the tubehead stay exactly where you position it without sagging or drifting.



## LET US HELP YOU DESIGN YOUR IMAGING WORKFLOW

Midmark design experts are ready to help you create a beautiful and functional space optimized for the way you work. From proper equipment sizing to layout, we'll work with you to design the best environment for your practice.

# Consistent. Durable. Unobtrusive.

The best imaging equipment doesn't demand attention—it simply works the way it should without unnecessary complications or interruptions. Midmark Preva is engineered for dependable performance and long-term reliability to help keep your operatory running smoothly. Easy to install, easy to maintain and easy to operate, it's built for years of consistent use.

## Things to Consider

### CONSISTENCY

Can you rely on your imaging equipment to deliver clear results every time?

**Clear, consistent imaging is essential for accurate diagnosis and treatment planning.** The **Midmark Preva Intraoral X-Ray** is engineered to perform reliably day after day, patient after patient, supporting confident clinical decisions and helping reduce the need for retakes.

### ACCURACY

Does your X-ray system support precise positioning without drifting or sagging?

**Precise positioning is critical for capturing the right images the first time.** With smooth arm movement and drift-free stability, the **Midmark Preva** supports accurate alignment to help minimize repositioning, improve efficiency and reduce variability.

### RELIABILITY

Can your equipment keep up with daily demands without interruption?

**Equipment should support your workflow—not interrupt it.** **Midmark Preva** is designed to deliver dependable performance throughout the day and long-term durability over years of use.



# Midmark® Products + Solutions



SCAN OR CLICK TO  
VIEW PRODUCTS

# Midmark® Preva™ Intraoral X-Ray

## DESIGNED TO CAPTURE SHARP IMAGES THE FIRST TIME

Detect dental caries, visualize tooth roots and evaluate surrounding bone with clear, consistent imaging. Designed for ease of use, Midmark Preva supports efficient workflows with smooth, one-handed positioning and reliable, drift-free performance.

## Features + Options





**Midmark Preva Intraoral X-Ray is the leading wall-mounted intraoral X-ray unit used by dental practitioners\* to:**

- Identify cavities and visualize tooth roots
- Examine the status of developing teeth
- Monitor the overall health of the bony areas surrounding teeth

**01** Select the configuration that best fits your operatory with wall-mounted (1A), cabinet-mounted (1B) and mobile-mounted (1C) options, plus multiple arm lengths.

**02** Support consistent, repeatable imaging with reliable exposure and stable positioning.



# What Are the Differences?

The **Midmark Preva™ Intraoral X-Ray** is available in two models—Preva DC and Preva 2.0. Both deliver the reliable performance you expect from Preva, but a few key differences can help you choose the right model.



01




02

**01 Preva DC is the base Preva unit.**

- Use Preva DC to capture images and transmit them to a sensor plugged into your workstation.
- Preva DC is available in either white or grey.

**02 Preva 2.0 features cabling in the arm and an integrated high-speed USB connector at the elbow so you can integrate your sensor directly into your X-ray unit.**

- Keep workstations and wires out of pathways with sensors plugged directly into the unit.
- Preva 2.0 is available in white only.



**For more information about all imaging products available from Midmark, please visit [midmark.com/dental](https://midmark.com/dental) or contact your Midmark representative at 1.800.MIDMARK.**

# Design with Midmark®

## DESIGNED FOR YOU

Our design experts are ready to help you personalize and equip your entire clinical environment—from the layout of every room to each sterilizer, chair, door and drawer, down to the smallest detail. Let us help you design the best environment for your practice.



## CONNECT VIRTUALLY FOR LIVE DESIGN

Connect with Midmark Live Design services to bring your vision to life. We help you visualize your space, understand how it flows and spot potential challenges so you can make confident decisions to support better care.



**SCAN OR CLICK  
TO SCHEDULE AN  
APPOINTMENT.**



## JOIN US IN PERSON FOR A MIDMARK EXPERIENCE TRIP

Explore equipment, finishes and floor plan options during a hands-on Midmark Experience trip. With expert guidance and a clear picture of your options, you'll be able to make more confident, informed decisions for your practice.



**SCAN OR CLICK TO  
REQUEST A VISIT**



## STREAMLINE YOUR DELIVERY + INSTALLATION

Ensure a smooth, on-time setup with Midmark Delivery Services. Our team coordinates the shipping, delivery, staging and setup of your Midmark equipment to help keep your project on track.



**SCAN OR CLICK TO LEARN  
MORE ABOUT MIDMARK  
DELIVERY SERVICES**

## SOURCES

<sup>1</sup> The Journal of the American Dental Association (JADA): [https://jada.ada.org/article/S0002-8177\(25\)00631-2/fulltext](https://jada.ada.org/article/S0002-8177(25)00631-2/fulltext). Accessed August 2026.



Designing better care.®



CARB 93120.2 Phase 2 Compliant and TSCA Title VI Compliant

Midmark intraoral X-ray systems are proudly designed, assembled and supported in the USA.

Midmark is an ISO 13485 and ISO 9001 Certified Company. Certain products are not included. See the complete list at: [midmark.com/ISO](http://midmark.com/ISO)

For more information, contact your authorized Midmark distributor or call: 1.800.MIDMARK Outside the USA call: 1.937.526.3662 or visit our website: [midmark.com](http://midmark.com)

© 2026 Midmark Corporation, Versailles, Ohio USA.  
Products subject to improvement changes without notice.

007-10053-00 Rev F1 (8/26)

