


DENTAL

Midmark Smart G-Vacuums



SPECIFICATIONS

Model	Rated users	Base unit HxWxD (inch)	Separator HxWxD (inch)	Actual Weight (lb)	Total HP	Voltage	Amps	Hertz	Recommended breaker size (Amps)	Inlet connection size (inch)	Fuses	Drain connection size (inch)
G3-050	3-5	28 x 22 x 20	36 x 22 x 22	Base unit 375 Separator 56	3	230	20	60	30	1 1/2	2.5 A, 250V Slow-Blow	1 1/2
G5-050	5-7	28 x 22 x 20	36 x 22 x 22	Base unit 375 Separator 56	3	230	20	60	30	1 1/2	2.5 A, 250V Slow-Blow	1 1/2
G7-050	7-10	28 x 22 x 20	36 x 22 x 22	Base unit 375 Separator 56	3	230	20	60	30	1 1/2	2.5 A, 250V Slow-Blow	1 1/2
G6-050	6-10	28 x 22 x 20 Each	36 x 22 x 22 Each	Base unit 375 Separator 56 Each	6	230	40	60	2 x 30	2	2.5 A, 250V Slow-Blow	1 1/2
G10-050	10-14	28 x 22 x 20 Each	36 x 22 x 22 Each	Base unit 375 Separator 56 Each	6	230	40	60	2 x 30	2	2.5 A, 250V Slow-Blow	1 1/2
G14-050	14-20	28 x 22 x 20 Each	36 x 22 x 22 Each	Base unit 375 Separator 56 Each	6	230	40	60	2 x 30	2	2.5 A, 250V Slow-Blow	1 1/2

CONNECTIVITY

Network Protocols			Link Layer: IEEE 802.3, IEEE 802.11 Internet Layer: IP Transport: TCP, UDP App Layer: TLS, MQTT, NTP, HTTP, etc.
Network Data	LAN Physical Connection	RJ-45	
	LAN Cabling	Ethernet (min. CAT 5)	
	LAN Max Speeds	10/100 Mbps	
	LAN Protocols	IEEE 802.3u	
	Wi-Fi Connection	The device supports connectivity to Wi-Fi using the following standards: 802.11a/b/g/n/ac	
	Operational Frequency	2.4 GHz or 5.0 GHz	
	Wi-Fi Security Standards	When connecting to Wi-Fi, the device uses the following security standards: WPA2-PSK, WPA2 PEAP/MSCHAPv2, WPA2 EAP-TLS	
	Network Interface Card Parameters	Not Applicable	